
JUN 2 9 2005



PREMARKET NOTIFICATION 510(K) SUMMARY

STALIF TT™ SYSTEM

Company:

Surgicraft, LTD

8 Padgets Lane South Moons Moat

Redditch, United Kingdom

B98 ORA

Contact:

Steve Trotman

Proposed Proprietary

Trade Name:

STALIF[™] TT

Classification Name:

888.3060 Orthopedics

FDA Product Code

Classification:

MQP, Class II

Device Description:

The Surgicraft STALIFTM TT is a radiolucent vertebral body replacement device used in conjunction with cancellous bone screws and is designed to restore biomechanical integrity from thoracic (T9 to T12) and lumbar spine (L1 to L5) following vertebrectomy or corpectomy for patients with spine tumors or trauma/ fracture. The system provides anterior column support both immediately after surgery and for prolonged periods in the absence of bone fusion.

Material:

PEEK-OPTIMA LT1

Intended Use:

The STALIF™ TT device is indicated for the replacement, partial vertebrectomy and / or augmentation of a vertebral body due to destruction or partial destruction by tumor or fracture/trauma leading to the restoration of planar alignment, restoration of the height of the collapsed vertebral body and indirectly facilitating neural decompression. The device is designed to restore biomechanical integrity of the anterior, middle and posterior spinal column even in the absence of fusion for a prolonged period. The device is intended for use as a vertebral body replacement or partial replacement in the thoracolumbar spine (from T9 to L5). The STALIFTM TT device may be used with bone graft.

Predicate Device:

DePuy AcroMed Stackable Cage[™] System Interpore Cross International GEO Structure

Substantial Equivalence:

Performance data were provided.



JUN 2 9 2005

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Surgicraft L.T. D. C/o Ms. Janet M. Webb MEDVantage Incorporated 950 N. Michigan Avenue #2202 Chicago, Illinois 60611

Re: K051027

Trade/Device Name: STALIF TT Regulation Number: 21 CFR 888.3060

Regulation Name: Spinal intervertebral body fixation orthosis

Regulatory Class: II Product Code: MQP Dated: June 13, 2005 Received: June 14, 2005

Dear Ms. Webb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Miriam C. Provost, Ph.D.

Acting Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

INDICATIONS FOR USE

Device Name: STALIF TT TM Indications for Use:
The STALIF TM TT device is indicated for the replacement, partial vertebrectomy and / or augmentation of a vertebral body due to destruction or partial destruction by tumor or fracture/trauma leading to the restoration of planar alignment, restoration of the height of the collapsed vertebral body and indirectly facilitating neural decompression. The device is designed to restore biomechanical integrity of the anterior, middle and posterior spinal column even in the absence of fusion for a prolonged period. The device is intended for use as a vertebral body replacement or partial replacement in the thoracolumbar spine (from T9 to L5). The STALIF TM TT device may be used with bone graft.
Prescription UseX or Over-The-Counter Use
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
(Division Sign-Off) Division of General, Restorative, and Neurological Devices
510(k) Number K05/077