

JUL 22 2005

Premarket Notification 510 (k)  
**Zirox**

**WIELAND**  
Dental + Technik

Wieland Dental + Technik  
GmbH & Co. KG  
Schwenninger Straße 13  
D-75179 Pforzheim  
Postfach 10 20 40  
D-75120 Pforzheim  
Telefon 0 72 31 / 37 05 -0  
Telefax 0 72 31 / 35 79 59

K051249

**5. 510 (k) Summary**

**Submitter of 510(k):** Wieland Dental + Technik GmbH & Co. KG  
Schwenninger Str. 13  
D-75179 Pforzheim  
Germany  
Phone: +49-7231-3705-0

**Contact person:** Dr. Gerhard Polzer  
**Phone:** +49-7231-3705-219  
**Fax:** +49-7231-357959  
**e-mail:** gerhard.polzer@wieland-dental.de

**Date of Summary:** 2005-04-27

**Trade name:** Zirox®

**Classification name:** Porcelain powder for clinical use  
**Product code:** EIH  
**C.D.R section:** 872.6660  
**Classification:** Class II

**Legally marketed  
equivalent device:** Cerabien ZR Dental Porcelain System

**510(k) number:** K 031968

## 510 (k) Summary

### Device description

**Zirox®** is a dental porcelain system that consists of 127 different ceramic powders.

It is intended to be used by professional dental technicians to manufacture all-ceramic dental appliances for the sole use of particular patients.

**Zirox®** is recommended for veneering zirconium oxide ( $ZrO_2$ ) frameworks, in which it provides an easy to use dental restorative material to fabricate dental restorations with the best possible aesthetic results.

The coefficient of thermal expansion [CTE  $_{(25-500^{\circ}C)}$ ] of these zirconium oxide ( $ZrO_2$ ) frameworks had to be approximately  $10 \times 10^{-6} K^{-1}$ .

**Zirox®** meets all applicable requirements of the standard ISO 6872: 1995 "Dental ceramic".

Type of Powder	Shades
Liner	A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; C3; C4; D2; D3; D4; White, Yellow, Violet, Orange, Brown, Gum
Dentine	A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; C3; C4; D2; D3; D4; Crystal 1 light; Crystal 2 pearl; Crystal 3 creme; Gum 1; Gum 2; Gum 3; Gum 4, Gum 5
Chroma-Dentine Chromatix	A1; A2; A3; A3.5; A4; B1; B2; B3; B4; C1; C2; C3; C4; D2; D3; D4
Dentine Modifier	Mocca; Corn; Mango; Caramel; Khaki; Brown; Yellow; Ivory;
Flu-Dentines	Flamingo; Straw; Bright; Sunny; Crystal;
Incisal / Enamel	Incisal 1, Incisal 2, Incisal 3; Incisal 4; Topas; Anthrazit; Amethyst; Aquamarin; Citrin; Rubin; Lemon; Melon; Transpa Neutral; Transpa Clear
Opale incisal	Opale Incisal 1, Opale Incisal 2, Opale Incisal 3, Opale Incisal 4, Crystal Blue, Frosty, Milky, Snow, Ice
Shoulder porcelains margin	High Flu, High 1, High 2, High 3, High 4, High Red, High Crystal; Low Flu, Low 1, Low 2, Low 3, Low 4, Low 5, Low Bleach
Stain	White; Black; Grey; Caramel; Orange; Ocker; Peach; Melon; Blue; Steel; Violet; Gum; Marone; Olive; Ivory; Yellow; Bodystain A; Bodystain B; Bodystain C; Bodystain D
Glaze	Glaze
Correction	Correction



JUL 22 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Gerhard Polzer  
Wieland Dental & Technik GMBH & Company KG  
Schwenninger Strabe 13  
D-75179 Pforzheim  
GERMANY

Re: K051249  
Trade/Device Name: ZIROX  
Regulation Number: 21 CFR 872.6660  
Regulation Name: Porcelain powder for clinical use  
Regulatory Class: II  
Product Code: EIH  
Dated: May 4, 2005  
Received: May 16, 2005

Dear Mr. Polzer:


We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", with a stylized flourish at the end.

Chiu Lin Ph.D.

Director

Division of Anesthesiology, General Hospital

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known):

K051249

Device Name: **Zirox**

Indications For Use:

Zirox is a dental ceramic that can be used by dental technicians to fabricate all-ceramic restorations by veneering zirconium oxide ( $\text{ZrO}_2$ ) frameworks.

The coefficient of thermal expansion [CTE  $_{(25-500^\circ\text{C})}$ ] of these zirconium oxide ( $\text{ZrO}_2$ ) frameworks has to equal approximately  $10 \times 10^{-6} \text{K}^{-1}$ .

Prescription Use ☒  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Robert DNS for Dr. Susan Runner

(Division Sign-Off)  
Division of Anesthesiology, General Hospital,  
Infection Control, Dental Devices

510(k) Number: K051249

Page 1 of 1