510 (k) Summary

As Required by 21 section 807.92 (c)

1. Submitter Name: Siam Sempermed Corp., Ltd

JUL 1 9 2005

2. Address:

110 Kanjanavanit Road. Pathong Hatyai

Songkhla. Thailand 90230

3. Phone:

(+66) 74 291 648

4. Fax:

(+66) 74 291 650

5. Contract Person: Mrs. Sureerat Choosri (Product Manager)

6. Date summary prepared: Mar 31, 2005

7. Official Correspondent: Sempermed USA Inc.

8. Address:

13900 49th Street North

Clearwater, USA, FL 33762

9. Phone:

727 787 7250

10. Fax:

727 787 7558

11. Contact person:

Mr. William E. Harris

- 12. Device Trade or Proprietary Name: Sterile, powder-free latex examination gloves with polymer coating and protein claim (50 micrograms or less).
- 13. Device Common or usual name: Examination glove
- 14. Device Classification Name: Glove, Patient Examination, Latex

15. Description of the Device:

Sterile, powder-free latex examination gloves with polymer coating and protein claim (50 micrograms or less)

16. Intended use of the device:

This device is a disposable device intended for medical purpose that is worn on the examiner 's hand to prevent contamination between patient and examiner.

17. Summary of The Technological Characteristics of The devices :

Sterile, powder-free latex examination gloves with polymer coating and protein claim (50 micrograms or less) are summarized with the following technological characteristics compared to ASTM or equivalent standards.

CHARACTERISTICS	STANDARDS	DEVICE PERFORMANCE	
Dimensions	ASTM D 3578-2001-ae ²	Meets	
Physical Properties	ASTM D 3578-2001-ae ²	Meets	
Freedom from pinholes	ASTM D 3578-2001-ae ²	Meets	
Powder Free Residue	ASTM D 6124-2001	≤2 mg/glove	
Protein Level	ASTM D 5712-99	<50 µg/g	
Biocompatability	Primary Skin Irritation in Rabbits	Passes	
	Guinea Pig Sensitization	Passes	

18. Substantial Equivalents Based on Assessment of Non-Clinical Performance Data

The performance test data of the non-clinical tests that support a determination of substantial equivalence is the same as mentioned immediately above.


19. Conclusion

It can be concluded that the Sterile, powder-free latex examination gloves with polymer coating and protein claim (50 micrograms or less) will perform according to the glove performance standards referenced in section 17 above and meet ASTM standards, and FDA requirements. Consequently, this device is substantially equivalent to currently marketed devices. This device is safe and effective as the predicate device *Siam Sempermed Latex Examination Glove*, *Powder free*. Indeed, it is equivalent This is better expressed in the tabulated comparison as below.

Technical comparison of specific elements is attached in the main submission.

FDA file reference number	510k number: k970794	
Attachments inside notification	REFER TO APPENDIX 1	
submission file		
TECHNOLOGICAL		
CHARACTERISTICS	Comparison result	
	REFER TO ADDITIONAL TECHNICAL	
	COMPARATIVE TABLE WITHIN 510K	
	SUBMISSION	
Indications for use	Identical	
Target population	Identical	
Design	Similar	
Materials	Similar	
Performance	Identical	
Sterility	different	
Biocompatibility	Identical	
Mechanical safety	Identical	
Chemical safety	Identical	
Anatomical sites	Identical	
Human factors	Identical	
Energy used and/or delivered	Identical (Not applicable)	
Compatibility with environment	Identical	
and other devices		
Where used	Identical	
Standards met	Identical	
Electrical safety	Identical (not applicable)	
Thermal safety	Identical (not applicable)	
Radiation safety	Identical (not applicable)	

<u>نقيني کاري تاجو</u>





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JUL 1 9 2005

Siam Sempermed Corporation Limited... c/o Mr. William E. Harris Sempermed USA, Incorporated 13900 49th Street North Clearwater, Florida 33762

Re: K051329

Trade/Device Name: Sterile Powder-Free Latex Examination Glove with Polymer

Coating and Protein Claim (50 micrograms or less)

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient examination glove

Regulatory Class: I Product Code: LYY Dated: July 6, 2005 Received: July 11, 2005

Dear Mr. Harris:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): k051329

(50 microc	rams or less)	on Glove with Polymer Coating				
Indications For Use: A pa	ndications For Use: A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between					
patie	nt and examiner					
		•				
Prescription Use	AND/OR	Over-The-Counter Use _	X			
(Part 21 CFR 801 Subpart D)	_	(21 CFR 801 Subpart C)				
	RITE BELOW THIS LINE-	CONTINUE ON ANOTHER	PAGE IF			
NEEDED)						
			, , , , , , , , , , , , , , , , , , ,			
Concurrence of CDRH, Office of Device Evaluation (ODE)						
Shud IM my	lug 160 7/18/08					
(Division Sign-Off) Division of Anesthesiology, G Infection Control, Dental Dev	eneral Hospital, ices	Page 1 c	of			