

K051571

EXHIBIT #1

**510(K) SUMMARY**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.

The assigned 510(k) number is:\_\_\_\_\_.

**1. Submitter's Identification:**

Dmetec Co., Ltd  
402-803 Tchno-Park, 193 Yakdae-dong, Wonmi-ku,  
Bucheon City, 420-734 Korea  
Tel: +82-32-234-1441  
Fax: +82-32-234-1444  
Web site: [www.dmetec.com](http://www.dmetec.com)  
Contact: J.W, Park

Date Summary Prepared: June 13, 2005

**2. Name of the Device:**

Skylight LED Curing Light

**3. Common or Usual Name:**

Polymerization Light-Curing Device

**4. Predicate Device Information:**

Mini L.E.D by SATELEC (K032465)

**5. Device Description:**

Skylight by DMETEC Co., Ltd is classified as an Ultraviolet Activator for polymerization(21 C.F.R. 872.6070) because it is a device intended for the photo-polymerization light cured dental materials.

The Skylight is a universal photo-polymerization light curing source working in cordless conditions and producing visible blue light in the 440 to 490 nm waveband of spectrum with a power density comprised 1,000W/cm<sup>2</sup>, 13-11mm light guide.

These power densities are sufficient for photo-polymerization in the 440-490nm waveband of visible light cured (VLC) dental materials of restorative composite, orthodontic brackets, orthodontic bonding and sealing.

The exposure time can be set by one second to 99 seconds with 3 modes of operation: Standard, Pulse, and Soft Start. It has also memory storage functions for 5 exposure programs set by the practitioner.

**6. Intended Use:**

The Skylight LED Curing Light is intended to polymerize resinous dental materials, restorative composite materials, and orthodontic brackets, bonding, and sealing materials that are photo-polymerized in the 440-490 nm waveband of visible light.

**7. Comparison to Predicate Devices: Mini L.E.D by SATELEC:**

The DMETEC Skylight L.E.D curing light has similar design and performance characteristics and identical intended use of the 510(K) cleared light curing unit, the SATELEC Mini L.E.D (K032465) for the photo-polymerization of dental materials of restorative composite, bonding and sealing.

These devices are well established and determined to be safe and effective.

**8. Discussion of Non-Clinical Tests Performed for Determination of Substantial Equivalence are as follows:**

Testing was conducted in accordance with EN 60601-1: 1990 +A, +B; EN 60601-1-2:2001; EN 61000-3-2:2000 and EN 61000-3-3, 1998/A1;2001 and EN 55011: 1995/A1: 1999; Class B supports testing information demonstrating safety and effectiveness of the Skylight LED curing light in the intended environment of use.

None of the testing demonstrated any design characteristics that violated the requirements of the standards or resulted in any safety hazards.

**9. Discussion of Clinical Tests Performed:**

No clinical testing was conducted

**10. Conclusions:**

The Skylight LED curing light is substantially similar to the predicate intended use, operation, safety and function, and is effective for its' intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUL 19 2005

DMETEC Company, Ltd.  
C/O Ms. Carolann Kotula  
Official Correspondent  
MDI Consultants, Inc.  
55 Northern Boulevard, Suite 200  
Great Neck, New York 11021

Re: K051571

Trade/Device Name: Skylight LED Curing Light  
Regulation Number: 21 CFR 872.6070  
Regulation Name: Ultraviolet Activator for Polymerization  
Regulatory Class: II  
Product Code: EBZ  
Dated: June 13, 2005  
Received: June 14, 2005

Dear Ms. Kotula:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu S. Lin, PhD

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K051571

Device Name: Skylight LED curing light

Indications For Use:

The Skylight LED Curing Light is intended to polymerize resinous dental materials, restorative composite materials, and orthodontic brackets, bonding and sealing materials that are photo-polymerized in the 440-490nm waveband of visible light.

Prescription Use X  
Use \_\_\_\_\_

(Per 21 CFR 801 Subpart D)

OR

Over-The Counter

(21 CFR 807 Subpart C)

Page 1 of \_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

R. Betz, DDS for Dr. Susan Runner  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital,  
Infection Control, Dental Devices

510(k) Number: K051571