

JUL 7 - 2005

K051624

V/I

**510(k) Summary of Safety and Effectiveness:  
T2® Recon Nail System Line Extension**

**Submission Information**

Name and Address of the Sponsor  
of the 510(k) Submission:

Howmedica Osteonics Corp  
325 Corporate Drive  
Mahwah, NJ 07430

For Information contact:

Vivian Kelly, Regulatory Affairs Specialist  
Howmedica Osteonics Corp.  
325 Corporate Drive  
Mahwah, NJ 07430  
Phone: (201) 831-5581  
Fax: (201) 831-6038

Date Summary Prepared:

June 17, 2005

**Device Identification**

Proprietary Name:

T2® Recon Nail System

Common Name:

Intramedullary Nail

Classification Name and Reference:

Intramedullary Fixation Rod and Accessories,  
21 CFR §888.3020

Device Product Code:

87 HSB

**Description:**

The T2® Recon Nail System is a family of IM Nails for various types of femoral fractures. This Special 510(k) submission is a line extension to the T2® Recon Nail System to add an alternate design of Set Screw to the system. There is no change in intended use for the subject device, which is provided below.

**Intended Use:**

The subject T2® Recon Nail System is a fracture fixation device comprised of Femoral Nails and the related accessories such as Washers, Locking Screws, Set Screws, End Caps, and Lag Screws.

The subject and predicate devices are intended to provide strong and stable internal fracture fixation with minimal soft tissue irritation. This device is utilized as an aid to healing, not as a substitute for normal intact bone and tissue.

The T2® Recon Nail indications include fixation of subtrochanteric, intertrochanteric, ipsilateral neck/shaft, comminuted proximal femoral shaft fractures, femoral fixation required as a result of pathological disease, and temporary stabilization of fractures of the femoral shaft ranging from the femoral neck to the supracondylar regions of the femur.

**Statement of Technological Comparison:**

The subject device is fabricated from titanium alloy and PEEK while the predicate device is made from titanium alloy. Mechanical testing demonstrates the comparable mechanical properties of the subject T2® Recon Nail System to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUL 7 - 2005

Ms. Vivian Kelley, RAC  
Regulatory Affairs Specialist  
Howmedica Osteonics Corporation  
325 Corporate Drive  
Mahwah, New Jersey USA 07430

Re: K051624  
Trade/Device Name: T2® Recon Nail System  
Regulation Number: 21 CFR 888.3020  
Regulation Name: Intramedullary fixation rod  
Regulatory Class: II  
Product Code: HSB  
Dated: June 17, 2005  
Received: June 20, 2005

Dear Ms. Kelley:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Vivian Kelley, RAC

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120 . Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Miriam C. Provost", is written over the typed name.

Miriam C. Provost, Ph.D.  
Acting Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K

Device Name: T2® Recon Nail System

### Indications For Use:

The T2® Recon Nail indications include fixation of subtrochanteric, intertrochanteric, ipsilateral neck/shaft, comminuted proximal femoral shaft fractures, femoral fixation required as a result of pathological disease, and temporary stabilization of fractures of the femoral shaft ranging from the femoral neck to the supracondylar regions of the femur.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR


Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

  
(Division Sign-Off)  
**Division of General, Restorative,  
and Neurological Devices**

510(k) Number   K051624