OCT 2 1 2005

EXHIBIT 2

510(k) Summary

KaVo Dental Corporation 340 East Main Street Lake Zurich, Illinois 60047 Toll Free: 800 323 8029

Tel: 847 / 550 - 6800 Fax: 847 / 550 - 6825

e-mail: info@kavousa.com Contact: John Miller, Director of RA/QA July 11, 2005

1. Identification of the Device:

Proprietary-Trade Name: KaVo DIAGNOdent® 2190

Classification Name: Dental Hand Instrument, Laser Fluorescence Caries Detection Device,


Product Codes NTK

Common/Usual Name: Laser fluorescence caries detection device

2. Equivalent legally marketed device: DIAGNOdent® 2095, K983658.

3. Indications for Use (intended use): For use as an aid in the diagnosis of dental caries.

4. Description of the Device: This submission is for a modification of a device system cleared under K983658, the DIAGNOdent® 2095. The modification is in the form of a new method of packaging the device, in the form of a pen-like device instead of the former electronic unit with a fiber optic probe.



5. Safety and Effectiveness, comparison to predicate device:

r			
Comparison	DIAGNOdent® 2095,	DIAGNOdent® 2190	
Areas	K983658		
Indications	For use as an aid in the	SAME	
for use	Diagnosis of Dental Caries		
Probe	Fiber Optic with sapphire	Probe is integrated into	
technology	tip.	the body of the hand held	
		unit, with sapphire tip.	
Light Source	655 nm <1 mw Laser	650 nm, 190 nw	
		maximum average radiant	
		power	
Laser power	Class II	Class I	
class			
Returned	Fluorescence	Fluorescence	
light			
Sterilization	Probe tip only, autoclave	SAME	
User interface	Numeric and audible tones	SAME except	
	LED numbers	LCD numbers	
Power source	6- AA Alkaline battery	1- AA Alkaline battery	
Target	Dentists' offices	SAME	
population			

6. Conclusion: In all important respects, the DIAGNOdent® 2190 is substantially equivalent to the DIAGNOdent® K983658. This conclusion is based on indications for use, bench, invitro, and clinical studies, as well as EMC and electrical safety testing.



OCT 2 1 2005

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

KAVO AMERICA C/O Daniel Kamm Regulatory Engineer Kamm & Associates P.O. Box 7007 Deer Field, Illinois 60015

Re: K051909

Trade/Device Name: DIAGNODENT Regulation Number: 21 CFR 872.1745

Regulation Name: Laser Fluorescence caries detection device

Regulatory Class: II Product Code: NBL Dated: October 12, 2005 Received: October 13, 2005

Dear Mr. Kamm:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

sufte y Michanoms

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if k	cnown): <u>K05</u>	1909		
Device Name:	DIAGNOde	ent® 2190		
Indications For Use	:			
For use as an aid in	the diagnosis	of dental caries.		
Prescription Use		AND/OR	Over-The-Counter Us (21 CFR 807 S	eubpart C)
(PLEASE DO 1 NEEDED)	NOT WRITE E	BELOW THIS I	LINE-CONTINUE ON .	ANOTHER PAGE IF
	Concurrence o	of CDRH, Office	e of Device Evaluation ((ODE)
				•
				Page 1 of 1
\sim				<u>-</u>

(Svision Sign-Off)
Elivision of Aneshesiology, General Hospital, Infection Control, Dental Devices

510(k) Number:___