

JAN 3 2006



K 051933

SYBRON DENTAL SPECIALTIES

Section III - 510(k) Summary of Safety and Effectiveness

Submitter:

Sybron Dental Specialties, Inc.
1717 W. Collins Avenue
Orange, California 92867
(714) 516-7484 - Phone
(714) 516-7488 - Facsimile
Colleen Boswell - Contact Person

Date Summary Prepared: July 2005

Device Name:

- Trade Name - *Expasyl Power Applicator*
- Common Name - Dental Handpiece
- Classification Name - Dental Handpiece and Accessories, per 21 CFR § 872.4200
-

Devices for Which Substantial Equivalence is Claimed:

- Electro Medical Systems, *EMS Air-Flow handy 2*

Device Description:

The Expasyl Power Applicator is a dental handpiece intended to be used to dispense Expasyl, a temporary hemostatic/gingival retracting agent. The paste is dispensed through a tip attached to the cartridge of Expasyl through the use of the *Expasyl Power Applicator*. The Expasyl Power Applicator is connected to any micro-motor (air or electric) that has a coupling socket and a maximum speed of 40,000 rpm. When the capsule of Expasyl is empty, the plunger on the Expasyl Power Applicator rotates and the applicator makes a repetitive clicking noise which signals the micro-motor must be stopped. The capsule holder of the *Expasyl Power Applicator* is autoclaveable.

Intended Use of the Device:

The intended use of the *Expasyl Power Applicator* is to dispense Expasyl, a temporary hemostatic/gingival retracting agent.

Substantial Equivalence:

The *Expasyl Power Applicator* is substantially equivalent to other legally marketed devices in the United States. The *Expasyl Power Applicator* functions in a manner similar to and is intended for a similar use as the product Air-Flow handy 2 cleared for marketing for Electro Medical Systems.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 3 2006

Ms. Colleen Boswell
Director, Corporate Compliance
Sybron Dental Specialties, Incorporated
1717 W. Collins Avenue
Orange, California 92867

Re: K051933

Trade/Device Name: Expasyl Power Applicator
Regulation Number: 21 CFR 872.3660
Regulation Name: Impression Material
Regulatory Class: II
Product Code: ELW
Dated: December 14, 2005
Received: December 15, 2005

Dear Ms. Boswell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Chiu Lin', with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: *Expasyl Power Applicator*

Indications For Use:

The *Expasyl Power Applicator* is a dental handpiece intended to be used to dispense Expasyl, a hemostatic/gingival retracting agent. The paste is dispensed through a tip attached to the cartridge of Expasyl through the use of the *Expasyl Power Applicator*.

Prescription Use x
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Suzanne Pinner
(Signature)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K051933

Page 1 of 1