



SEP 30 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

SonoScape, Inc.  
% Mr. Bob Leiker  
Quality & Regulatory Services  
7263 Cronin Circle  
DUBLIN CA 94568

Re: K052042

Trade Name: SonoScape SSI-1000/SSI-5000™  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulation Number: 21 CFR 892.1560  
Regulation Name: Ultrasonic pulsed echo imaging system  
Regulation Number: 21 CFR 892.1570  
Regulation Name: Diagnostic ultrasonic transducer  
Regulatory Class: II  
Product Code: IYN, IYO, and ITX  
Dated: July 25, 2005  
Received: July 28, 2005

Dear Mr. Leiker:

This letter corrects our substantially equivalent letter of August 9, 2005 regarding a wrong transducer number.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the SonoScape SSI-1000/SSI-5000™, as described in your premarket notification:

Transducer Model Number

2P1, PA2.5 MHz Phased Array  
L741, LA7.5 MHz Linear Array

C344, CLA3.5 MHz Curved Linear Array  
MPTEE, 7-4 MHz/10mm/64 element Multi-Plane Transesophageal Phased Array  
5P1, 7-4 MHz/10mm/64 element Phased Array  
C611, 8-4 MHz/R11mm/128 element Micro-convex Array  
6V1, 8-4 MHz/R11mm/128 element Transvaginal Micro-convex Array  
10L1, 12-5 MHz/38mm/128 element Flat Linear Array

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,

  
for Nancy C. Brogdon

Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: SonoScape SSI-1000/SSI-5000

Diagnostic Ultrasound Pulsed Echo System

Diagnostic Ultrasound Pulsed Doppler Imaging System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal		P	P	P		P	P		Note 1	N
Ob/GYN		P	P	P		P	P		Note 1	N
IntraOperative										
Neurological										
Pediatric		N	N	N	N	N	N		Note 1	N
Small Organ (breast, thyroid, testes)		P	P	P		P	P		Note 1	N
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P	P		Note 1	N
Transesophageal		N	N	N	N	N	N		Note 1	N
Trans-Rectal		N	N	N		N	N		Note 1	N
Trans-Vaginal		N	N	N		N	N		Note 1	N
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular		P	P	P		P	P		Note 1	N
Laparoscopic										
Muscular-Skeletal Conventional		P	P	P		P	P		Note 1	N
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRL, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Section 4.3

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

Indications For Use

510(k) Number

Page 2 of 10

K052042

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: 2P1, PA2.5 MHz Phased Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Ob/GYN										
IntraOperative										
Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P	P		Note 1	P
Transesophageal										
Trans-Rectal										
Trans-Vaginal										
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular										
Laparoscopic										
Muscular-Skeletal Conventional										
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRE, Office of Device Evaluation (ODE)

(Division Sign-Off)

Prescription Use (Per 21 CFR 801.109) Division of Reproductive, Abdominal, and Radiological Devices

Section 4.3

510(k) Number

K052042

Page 3 of 10

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: L741, LA7.5 MHz Linear Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Ob/GYN										
IntraOperative										
Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)		P	P	P		P	P		Note 1	P
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Trans-Rectal										
Trans-Vaginal										
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular		P	P	P		P	P		Note 1	P
Laparoscopic										
Muscular-Skeletal Conventional		P	P	P		P	P		Note 1	P
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDREH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal  
and Radiological Devices510(k) Number K052042

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: C344, CLA3.5 MHz Curved Linear Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

## Mode of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal		P	P	P		P	P		Note 1	P
Ob/GYN		P	P	P		P	P		Note 1	P
IntraOperative										
Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Trans-Rectal										
Trans-Vaginal										
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular										
Laparoscopic										
Muscular-Skeletal Conventional										
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRR, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal  
and Pathological Devices

Prescription Use (Per 21 CFR 801.109)

Section 4.3

Indications For Use

K052042

Page 5 of 10

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: MPTEE, 7-4MHz/10mm/64 element Multi-Plane Transesophageal Phased Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Ob/GYN										
IntraOperative										
Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal		N	N	N	N	N	N		Note 1	N
Trans-Rectal										
Trans-Vaginal										
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular										
Laparoscopic										
Muscular-Skeletal Conventional										
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRE, Office of Device Evaluation (ODE)

*Nancy Brogdon*  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K052042

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: 5P1, 7-4MHz/10mm/64 element Phased Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

## Mode of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Opthalmic										
Fetal										
Abdominal										
Ob/GYN										
IntraOperative										
Neurological										
Pediatric		N	N	N	N	N	N		Note 1	N
Small Organ (breast, thyroid, testes)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N	N	N	N	N		Note 1	N
Transesophageal										
Trans-Rectal										
Trans-Vaginal										
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular										
Laparoscopic										
Muscular-Skeletal Conventional										
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRL, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Section 4.3

*Nancy Grodon*  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

Indications for Use Number K052042 Page 7 of 10



## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: C611, 8-4MHz/R11mm/128 element Micro-convex Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

## Mode of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Ophthalmic										
Fetal										
Abdominal										
Ob/GYN										
IntraOperative										
Neurological										
Pediatric		N	N	N		N	N		Note 1	N
Small Organ (breast, thyroid, testes)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N	N		N	N		Note 1	N
Transesophageal										
Trans-Rectal										
Trans-Vaginal										
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular										
Laparoscopic										
Muscular-Skeletal Conventional										
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Section 4.3

*Nancy C. Brogdon*  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 Indications for Use  
 510(k) Number K052042

Page 8 of 10

## Diagnostic Ultrasound System Indications for Use

510(k) Number: \_\_\_\_\_

Device Name: 6V1, 8-4MHz/R11mm/128 element Transvaginal Micro-convex Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

## Mode of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
Opthalmic										
Fetal										
Abdominal										
Ob/GYN		N	N	N		N	N		Note 1	N
IntraOperative										
Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Trans-Rectal		N	N	N		N	N		Note 1	N
Trans-Vaginal		N	N	N		N	N		Note 1	N
Trans-Urethral										
Intra-Vascular										
Peripheral Vascular										
Laparoscopic										
Muscular-Skeletal Conventional										
Muscular-Skeletal Superficial										
Others (Specify)										

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDRL, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Section 4.3

*Nancy C. Brogdon*  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
Indications for Use Number 5052042 Page 9 of 10

## Diagnostic Ultrasound System Indications for Use

510(k) Number: K052042

Device Name: 10L1, 12-5MHz/38mm/128 element Flat Linear Array  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation						Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined	Tissue Harmonic Imaging
	A	B	M	PWD	CWD						
Opthalmic											
Fetal											
Abdominal											
Ob/GYN											
IntraOperative											
Neurological											
Pediatric											
Small Organ (breast, thyroid, testes)		N	N	N			N	N		Note 1	N
Neonatal Cephalic											
Adult Cephalic											
Cardiac											
Transesophageal											
Trans-Rectal											
Trans-Vaginal											
Trans-Urethral											
Intra-Vascular											
Peripheral Vascular		N	N	N			N	N		Note 1	N
Laparoscopic											
Muscular-Skeletal Conventional		N	N	N			N	N		Note 1	N
Muscular-Skeletal Superficial											
Others (Specify)											

N = new indication P = previously cleared by FDA E = added under Appendix E

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Power Doppler; B/Color Doppler/PWD and B/Power Doppler/PWD

Additional Comments:

Concurrence of CDREK, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal, and  
Radiological Devices

Prescription Use (Per 21 CFR 801.109)

Section 4.3

510(k) Number K052042  
Indications For Use

Revised Page 10 of 10

AUG 9 - 2005

## PREMARKET NOTIFICATION [510(k)] Summary

K052042

**Trade Name:** SonoScape Ultrasound System, SSI-1000/SSI-5000™ with 6V1  
Micro-convex EndoCavity Linear Array, 10L1 Linear Array, MPTEE  
Multi-plane Trans-Esophageal Phased Array, 5P1 Phased Array, and  
C611 Micro-Convex Linear Array transducers.

**Common Name:** Diagnostic Ultrasound System and Transducers

**Classification Names:** Ultrasonic Pulsed Echo Imaging System, 90 IYO  
Ultrasonic Pulsed Doppler Imaging System, 90 IYN  
Diagnostic Ultrasound Transducer, 90 ITX

**Manufacturer's Name:** SonoScape Company Limited  
4/F., Yizhe Building, Yuquan Road,  
Nanshan, 518051, Shenzhen, China

**Contact:** Mr. Jinzhong Yao, President  
Telephone: (86) 755-26722890  
Fax: (86) 755-26722850

**U.S. Agent:** Bob Leiker  
Quality & Regulatory Services, Inc.  
Dublin, CA 94568

**Predicate Devices:** SonoScape SSI-1000/SSI-5000, K042369 and the GE LogiQ 500,  
K970901, K991611, and K010329

**Device Description:** The SonoScape SSI-1000/SSI-5000 ultrasound system is an integrated preprogrammed color ultrasound imaging system, capable of producing high detail resolution intended for clinical diagnostic imaging applications.

The all digital architecture with progressive dynamic receive focusing allows the system to maximize the utility of all imaging transducers to enhance the diagnostic utility and confidence provided by the system. The exam dependent default setting allows the user to have minimum adjustment for imaging the patient, while the in depth soft-menu control allows the advanced user to set the system for different situations.

The SonoScape System can be configured either as a portable (SSI-1000) model, or as a roll-around model on wheels (SSI-5000).

**Intended Use:** The device is a general-purpose ultrasonic imaging instrument intended for use by a qualified physician for evaluation of Abdomen; Pediatric; Small Organ (breast, tests, thyroid); heart soft tissue; Peripheral Vascular, Musculo-skeletal (conventional) and Urology.

### Technological Characteristics

<b>Display Modes</b>	Single and dual 2-D; Display of Duplex 2-D/M-mode; 2-D/Pulsed Doppler and Triplex 2-D/CD/Pulsed Doppler image formats; Dual B and Color in real time
<b>Measurements</b>	Distance; area; circumference; calipers; velocity, PI, RI. Cardiac. OB and Vascular package
<b>Principle of Operation</b>	Applying high voltage burst to the Piezoelectric material in the transducer and detect the reflected echo to construct the 2-D B-mode, Doppler color, and Doppler spectrum image for diagnostic purpose.
<b>Operating Controls</b>	<ul style="list-style-type: none"> <li>• TGC 8 slider, +/- 24dB</li> <li>• Depth Range: 3 to 24 cm</li> <li>• Image sector size: 32 lines to full B (256 lines)</li> <li>• Image Sector position: Steering within full maximum</li> <li>• B orientation flip: L/R key with marking on the screen</li> <li>• B Dynamic range control: preset 5 curves over 50-90 dB</li> <li>• Gray Scale Control: 8 Settings</li> <li>• Focal Number: 16 focal zone setting</li> <li>• B persistence: 30-90% recursive</li> <li>• Image Processing: Smoothing, edge enhancement</li> <li>• PW sweeping speed 2,4,8 sec over display.</li> <li>• PW Wall filter setting: 16 settings, 0.25 to 20% of PRF</li> <li>• PW sample volume: 0.5 to 10mm with 0.5mm step size.</li> <li>• PW/B update: with UPDATE key</li> <li>• PW cursor steering: Steer soft key</li> <li>• PW angle correction: 0 to 70 degree user control</li> <li>• PW trace: Peak, Mean</li> <li>• PW spectrum dynamic range: 5 preset curve over 15-48 dB</li> <li>• Spectrum baseline shift and invert</li> <li>• Color ROI setting: trackball and set key to control size and position</li> <li>• Color steering on flat probe: +, 0, -</li> <li>• Color Wall Filter: Color wall filter with 16 selection, 0.25-20% of PRF</li> <li>• Color &amp; B priority: C-B priority soft menu</li> <li>• Color Packet size: preset per Exam range from 8 to 12</li> <li>• Color spatial filter: preset per Exam, horizontal, vertical, off</li> <li>• Zoom factor: 1 to 10 continuously</li> <li>• Cine control: step, play backward, play continuously</li> </ul>
<b>Acoustic Output</b>	Track 3; MI, TIS, TIC, TIB Derated Ispta: 720mW/cm <sup>2</sup> maximum, TIS/TIB/TIC:0.1-4.0 Range, Mechanical Index: 1.9 Maximum, or Derated Isppa: 190 W/cm <sup>2</sup> max
<b>Display Annotations</b>	Time/date/year; transducer type; power output in %; frames per second; persistence; and compression settings; image depth; patient name and ID; institution name; focal position; TGC curve display; Doppler & M scale in sec; Doppler angle correction cursor; EKG trace; probe tip temperature / angle for MPTEE; free form annotation anywhere on image; trackball controlled; selective or global erase of the display annotations, body markers with transducer annotations
<b>Safety Compliance</b>	IEC601-1 International Electrotechnical Commission; Medical Electrical Equipment IEC60601-2 International Electrotechnical Commission; Electromagnetic Compatibility