

NOV - 8 2005



# CAPINTEC, INC.

September 7, 2005

K052595

RE: Summary of Safety and Effectiveness Information for the Capintec CRC 15Ultra.

The Capintec CRC 15Ultra is a multiple detector system, which may include up to eight sealed pressurized re-entrant well chambers, plus a high energy beta counter, and a well counter. The CRC 15Ultra combines the features of several well established Capintec models: CRC 35R, CRC 15R, Beta Counter, and Caprac R Counter into one unit with flexible configurations which are field upgradeable. The 15Ultra incorporates enhanced electronic circuitry and faster microprocessor capabilities, which extends its capability to network eight ionization chambers, as did the CRC 35R. However, the CRC 15Ultra can also add the detector assembly from the Beta C Counter, or the detector assembly from the Caprac R Well Counter, to the same universal readout. The basic design concepts, functionality, calculations, algorithms, and response are the same. The 15Ultra is a far more cost effective solution to a user who requires multiple chambers and detectors when compared to the cost of purchasing three separate units.

The predicate devices, upon which the CRC15Ultra is based, have a long history of safe and effective use in the field. The unit is intended for use by trained nuclear medicine technologists, nuclear medicine physicians, or medical physicists in clinical applications.

The CRC 15Ultra has also been tested and approved to the following safety standards for laboratory equipment:

- IEC 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use-Part 1 General Requirements
- IEC 61010-2-101 Safety requirements for electrical equipment for measurement, control, and laboratory use-Part 2 Particular Requirement for In Vitro Diagnostic Equipment
- IEC 61326 Electrical equipment for electrical equipment for measurement, control, and laboratory use-EMC requirements
- UL 61010A-1 Electrical Equipment for Laboratory Use; Part 1: General Requirements
- CAN/CSA-C22.2 No. 61010-1:2004 Standard for Safety Electrical Equipment for Measurement, Control, and Laboratory Use; Part 1: General Requirements
- CAN/CSA-C22.2 No. 61010-2-101-04 Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use; Part 2-101: Particular Requirements for In Vitro Diagnostics (IVD) Medical Equipment



NOV - 8 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Mary Ann Dell  
Vice President & General Manager  
CAPINTEC, INC.  
540 Alpha Drive  
PITTSBURGH PA 15238

Re.: K052595  
Trade/Device Name: CRC 15ULTRA Dose  
Calibrator  
Regulation Number: 21 CFR 892.1360  
Regulation Name: Radionuclide dose  
calibrator  
Regulatory Class: II  
Product Code: KPT  
  
Regulation Number: 21 CFR 892.5730  
Regulation Name: Radionuclide brachytherapy  
source  
Regulatory Class: II  
Product Code: JAQ  
Dated: September 6, 2005  
Received: September 21, 2005

Dear Ms. Dell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registrations, listing of devices, good manufacturing practice, labeling and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

21 CFR 876.xxxx	(Gastroenterology/Renal/Urology)	240-276-0115
21 CFR 884.xxxx	(Obstetrics/Gynecology)	240-276-0115
21 CFR 892.xxxx	(Radiology)	240-276-0120
Other		240-276-0100

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

Appendix # 2

Indications for Use Form-

Page \_\_\_\_ of \_\_\_\_

510(k) Number K052595

Device Name: CRC 15Ultra

Indications For Use: The CRC 15Ultra multiple chamber dose calibrator, which includes optional beta counter and well counter, is intended to be used by qualified nuclear medicine technologists to measure a wide range radiopharmaceuticals, including high energy beta and gamma emitters. It is also designed for use by trained medical physicists to measure the output of most radioactive brachytherapy sources, including HDR, LDR and IVBT sources. All brachytherapy sources must be measured in the appropriate source holder. The well counter is designed for measurement of low activity radioactive sources or solutions, including laboratory test applications. This device is used in numerous research applications for measurement of radioactive materials.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

(Optional Format 1-2-96)

Nancy C. Brogdon  
(Division Sign-Off)  
Division of Reproductive, Abdominal, and  
Radiological Devices  
510(k) Number K052595