

NOV 28 2005

K052900,1/4

Section 5 - 510(k) Summary

1. Applicant Contact:

Lois Smart
Director, Quality Assurance and Regulatory Affairs
Quill Medical, Inc.
2505 Meridian Drive, Suite 150
Research Triangle Park, NC 27713
Phone: 919-806-1961
Fax: 919-806-1953
Email: lsmart@quillmedical.com

Date Prepared:

- 2. Name of Device:** Quill® Nonabsorbable Nylon Barbed Suture
Common Name: Nonabsorbable Nylon Surgical Suture
Classification Name: Nonabsorbable Nylon Surgical Suture
Regulation 21 CFR 878.5020, Product Code GAR

3. Identification of device(s) to which the submitted claims equivalence:

The Quill® Nonabsorbable Nylon Barbed Suture is substantially equivalent to the following predicate devices:

- a. Predicate for Material:
- Nonabsorbable Nylon Surgical Sutures by Surgical Specialties Corp., 510(k) K930825
 - Grams Nylon Nonabsorbable Suture by GramsMed, LLC, 510(k) K003000
- b. Predicate for Indication for Use based on Technological Characteristics:
- Quill® Nonabsorbable Polypropylene Barbed Suture by Quill Medical, Inc., 510(k) K052373

4. Device Description:

The Quill® Nonabsorbable Nylon Barbed Suture is a monofilament, flexible thread prepared from long chain aliphatic polymers Nylon 6 and Nylon 6,6 (per 21 CFR 878.5020). It is available sterile, dyed black (logwood extract per 21 CFR 73.1410), dyed blue (FD&C Blue No. 2 per 21 CFR 74.3102) or undyed in various suture diameters, lengths and needle configurations in USP Sizes 2, 1, 0, 2-0, 3-0. Each suture has bi-directional barbs along the long axis of the suture monofilament.

Section 5 - 510(k) Summary (continued)

The Quill® Nonabsorbable Nylon Barbed Sutures approximate tissues by using the opposing barbs on the suture surface to imbed in the tissues after the surgeon precisely places the suture within the tissues. Each Quill® Nonabsorbable Nylon Barbed Suture pass provides the security of an interrupted suture strand without the added bulk of a knot. As with interrupted sutures, if the Quill® Nonabsorbable Nylon Barbed Suture breaks, the remaining suture passes will hold the wound edges in approximation.

5. Intended Use of the Device:

Quill® Nonabsorbable Nylon Barbed Sutures are indicated for use in soft tissue approximation excluding closure of the epidermis.

6. Characteristics of the device in comparison to those of the predicate device(s)**Indication for Use and Technology Comparison:**

The Quill® Nonabsorbable Nylon Barbed Suture is equivalent to the Quill® Nonabsorbable Polypropylene Barbed Suture in its intended use of soft tissue approximation and the technology of using barbs instead of knots to hold the tissue in approximation. In addition, the Nonabsorbable Nylon Surgical Suture manufactured by Surgical Specialties Corp. and GramsMed, LLC have an intended use of soft tissue approximation.

Material Comparison:

The Quill® Nonabsorbable Nylon Barbed Suture is equivalent to the Nonabsorbable Nylon Surgical Suture manufactured by Surgical Specialties Corp. and GramsMed, LLC as identical materials (nylon fiber, needles & packaging materials) and sterilization method is utilized.

Section 5 - 510(k) Summary (continued)

The comparison of the predicate devices to the new device is summarized below:

	Quill® Nonabsorbable Nylon Barbed Sutures, 510(k) TBD	Surgical Specialties Corp. Nonabsorbable Nylon Surgical Sutures, 510(k) K930825	GramsMed, LLC, Grams Nylon Nonabsorbable Sutures, 510(k) K003000	Quill® Nonabsorbable Polypropylene Barbed Sutures, 510(k) K052373
Product Code	Identical – GAR	Identical – GAR	Identical – GAR	Different – GAW
Suture Characteristic	Identical – Nonabsorbable Monofilament	Identical – Nonabsorbable Monofilament	Identical – Nonabsorbable Monofilament	Identical – Nonabsorbable Monofilament
Intended Use	Identical – Soft tissue approximation with applicable Warnings	Identical – Soft tissue approximation	Identical – Soft tissue approximation	Identical – Soft tissue approximation with applicable Warnings
Technique of Deployment	Identical – Attached needles	Identical – Attached needles	Identical – Attached needles	Identical – Attached needles
Technological Characteristic	Identical – Bi- directional barbs along the long axis of the suture monofilament	Different – Suture monofilament that utilizes knots to secure the suture	Different – Suture monofilament that utilizes knots to secure the suture	Identical – Bi- directional barbs along the long axis of the suture monofilament
Material	Identical – Nylon (cleared per K930825 or K003000)	Identical – Nylon	Identical – Nylon	Different – Polypropylene
Sterilization	Identical – EO	Identical – EO	Identical – EO	Identical – EO
Packaging	Identical – Device wound onto cardboard inner support card and packaged in a Tyvek pouch)	Identical – Device wound onto cardboard inner support card and packaged in a Tyvek pouch)	Identical – Device wound onto cardboard inner support card and packaged in a Tyvek pouch)	Identical – Device wound onto cardboard inner support card and packaged in a Tyvek pouch)

7. Safety and Performance:

The difference between the Quill® Nonabsorbable Nylon Barbed Suture and the above mentioned predicate devices do not raise any questions regarding the safety and effectiveness of the barbed suture. The Quill® Nonabsorbable Nylon Barbed Suture employs the same technological characteristics to support the intended use of soft tissue approximation as the Quill® Nonabsorbable Polypropylene Barbed Suture. In addition, nylon (as used in the predicate device by Surgical Specialties Corp. and the predicate device by GramsMed, LLC) is commonly used in medical applications and has proven to be biocompatible. The device, as designed, is as safe and effective as its predicate devices.

Section 5 - 510(k) Summary (continued)**8. Conclusion**

Based on the design, material, function and intended use discussed herein, Quill Medical believes the Quill® Nonabsorbable Nylon Barbed Suture is substantially equivalent to predicate devices currently marketed under the Federal Food, Drug and Cosmetic Act.



NOV 28 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Lois Smart
Director, Quality Assurance and Regulatory Affairs
Quill Medical, Inc.
2505 Meridian Drive, Suite 150
Research Triangle Park, North Carolina 27713

Re: K052900
Trade/Device Name: Quill® Nonabsorbable Nylon Barbed Suture
Regulation Number: 21 CFR 878.5020
Regulation Name: Nonabsorbable polyamide surgical suture
Regulatory Class: II
Product Code: GAR
Dated: October 13, 2005
Received: October 14, 2005

Dear Ms. Smart:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Barbara Melkerson" with a stylized flourish at the end.

Mark N. Melkerson
Acting Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K052900

Section 4 - Indications for Use Statement

510k number if known: _____

Device Name: Quill® Nonabsorbable Nylon Barbed Suture

Indications for Use:

Quill® Nonabsorbable Nylon Barbed Sutures are indicated for soft tissue approximation excluding closure of the epidermis.

Prescription Use ✓
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

**Division of General, Restorative,
and Neurological Devices**

510(k) Number K052900