

DAVANTIS/ANTELYS Vertebral Body Replacement
510(k) Summary

K053292
Page 1 of 1

<u>Date</u>	November 21, 2005	
<u>Submitter</u>	Scient'x Batiment Calypso Parc Ariane 3 78284 Guyancourt FRANCE	JAN 19 2006
<u>Contact person</u>	J.D. Webb 1001 Oakwood Blvd Round Rock, TX 78681 512-388-0199	
<u>Trade Name</u>	DAVANTIS/ANTELYS Vertebral Body Replacement	
<u>Common name</u>	Vertebral Body Replacement Device	
<u>Classification name</u>	Class II per 21 CFR section 888.3060	
<u>Product Code</u>	MQP	
<u>Equivalent Device</u>	CO Vertebral Body Replacement (K050348)	

Device Description

The DAVANTIS/ANTELYS Vertebral Body Replacements are oval shaped. The interior of the spacers has three compartments to provide space for bone graft. The implant is always implanted in the vertical position. These components are available in various heights (9mm-75mm). The superior and inferior surfaces of these components have ridges to interface with the vertebral endplates to resist rotation and migration. The body has a multitude of holes to allow additional impaction of bone graft.

The DAVANTIS/ANTELYS components are fabricated from pure poly(ether ether ketone) (PEEK). This material closely matches the modulus of elasticity of cortical bone, improving the biomechanical interface and reducing the stress shielding effect. They are fully radio-translucent, which enables optimum follow-up with diagnostic imaging, as the interbody fusion progresses. Two metal wires at the opposite ends of the spacers allow radiological confirmation of the cage position post operatively.

Intended Use

DAVANTIS/ANTELYS Vertebral Body Replacements are vertebral body replacements for use in the lumbar and thoracic spine (T4-L5) to replace a damaged, collapsed or unstable vertebral body due to tumor or trauma (i.e. fracture). These are not stand-alone devices, ISOBAR Ø6.2 Hemispherical Screws with Offset Clamps and Ø5.5 Rods must be utilized to enhance the stability of the reconstruction in skeletally mature patients following full or partial corpectomy.

Summary Nonclinical Tests

The DAVANTIS/ANTELYS Vertebral Body Replacement is a modification to the Scient'x CO VBR (K050348). It is similar in design and strength, is manufactured from the same material and has the same indications.



JAN 19 2006

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Scient'x
c/o J. D. Webb
The OrthoMedix Group, Inc.
1001 Oakwood Boulevard
Round Rock, Texas 78681

Re: K053292
Trade/Device Name: Davantis/Antelys Vertebral Body Replacements
Regulation Number: 21 CFR 888.3060
Regulation Name: Spinal intervertebral body fixation orthosis
Regulatory Class: II
Product Code: MQP
Dated: December 22, 2005
Received: December 23, 2005

Dear Sir/Madam:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Rogers

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a small "for" written below the main signature.

Mark N. Melkerson,
Acting Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

cc: HFZ-401 DMC
HFZ-404 510(k) Staff
HFZ-410 Division
D.O.
f/t:SMC:rrr: 1/19/06

Indications for Use

510(k) Number (if known): K053292Device Name: Davantis/Antelys Vertebral Body Replacements

Indications for Use:

Davantis/Antelys Vertebral Body Replacements are vertebral body replacements for use in the lumbar and thoracic spine (T4-L5) to replace a damaged, collapsed or unstable vertebral body due to tumor or trauma (i.e. fracture). These are not stand-alone devices, ISOBAR Ø6.2 Hemispherical Screws with Offset Clamps and Ø5.5 Rods must be utilized to enhance the stability of the reconstruction in skeletally mature patients following full or partial corpectomy.

Prescription Use X

(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use _____

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE
OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

**Division of General, Restorative,
and Neurological Devices**

510(k) Number K053292