

510(k) Summary

510(k) Owner: Xodus Medical, Inc.
Westmoreland Business & Research Park
702 Prominence Drive
New Kensington, PA 15068
Phone: 724-337-5500
Fax: 724-337-0555
Contact: Craig Kaforey (President)

Establishment Registration Number: 2530138

Date Prepared: 12/05/05

Device Information

Trade/Device Name: Cautery Tip Cleaner
Common Name: Cautery Tip Cleaner
Classification Name: Electrosurgical, Cutting & Coagulation Device & Accessories
Regulation Number: 21 CFR 878.4400
Product Code: GEI
Regulatory Class: II

Predicate Device

Device Name: Electro-Cautery Tip Cleaner
Common Name: Cautery Tip Cleaner
510 (k) Number: K030339
510 (k) Owner: Aspen Surgical Products
Classification Name: Electrosurgical, Cutting & Coagulation Device & Accessories
Regulation Number: 21 CFR 878.4400
Product Code: GEI
Regulatory Class: II

510(k) Summary

Device Description

The disposable Cautery Tip Cleaner is a small 2" x 2" square polyurethane foam pad which features a textile abrasive layer with an adhesive backing. A radiopaque strip within the adhesive makes the device x-ray detectable. The adhesive backing allows for universal placement as well as allowing the device to remain in place while the tip of the electrosurgical cauterization device is "scratched" on the abrasive surface to remove eschar buildup.

Intended Use

The disposable Cautery Tip Cleaner is a single use sterile product. Its intended use is as an electrosurgical accessory to clean uncoated cautery blades that are part of an electrosurgical pencil. The cautery blade is "scratched" on the cautery tip cleaner to remove eschar build-up during surgical procedures to allow the cautery blade to function effectively throughout the procedure. The product is usually placed somewhere on the sterile field, typically on the mayo stand. This product does not come in contact with the patient.

Technological Characteristics Comparison

Xodus Medical Inc.'s cautery tip cleaner has the same physical characteristics, material and design as the predicate device. The Xodus Medical Cautery Tip Cleaner is a 2" X 2" x ¼" square polyurethane foam pad featuring a textile abrasive layer with an adhesive backing containing an x-ray detectable radiopaque strip. Aspen Surgicals' Cautery Tip Cleaner has the same design but is a little smaller measuring 1½" x 1½" x ¼". Both the predicate device and Xodus Medicals' device are used for the cleaning of electrosurgical cautery tips during surgical procedures.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 19 2006

Mr. Craig Kaforey
President
Xodus Medical, Inc.
702 Prominence Drive
New Kensington, Pennsylvania 15068

Re: K053433

Trade/Device Name: Cautery Tip Cleaner
Regulation Number: 21 CFR 878.4400
Regulation Name: Electrosurgical cutting and coagulation device and accessories
Regulatory Class: II
Product Code: GEI
Dated: December 5, 2005
Received: December 9, 2005

Dear Mr. Kaforey:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Acting Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K053433

Device Name: Cautery Tip Cleaner

Indications for Use:

The disposable Cautery Tip Cleaner is a single use sterile product. Its intended use is as an electrosurgical accessory to clean uncoated cautery blades that are part of an electrosurgical pencil. The cautery blade is "scratched" on the cautery tip cleaner to remove eschar build-up during surgical procedures to allow the cautery blade to function effectively throughout the procedure. The product is usually placed somewhere on the sterile field, typically on the mayo stand. This product does not come in contact with the patient.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

**Division of General, Restorative,
and Neurological Devices**

Page 2.1

510(k) Number K053433