

FEB 10 2006

K053602

1.4 510(k) Summary of Safety and Effectiveness

Submitted by: Herbert Crane
Director, Regulatory Affairs

Address: Nobel Biocare USA LLC
22715 Savi Ranch Parkway
Yorba Linda, CA 92887

Telephone: (714) 282-5074

Facsimile: (714) 282-9023

Date of Submission: December 22, 2005

Classification Name: Endosseous Implant (21 CFR 872.3640)

Trade or Proprietary
or Model Name: Procera® Software

Legally Marketed Device(s): Nobel Biocare – Teeth in an Hour (K030685)
Nobel Biocare – Guided Surgery (K050393)

Device Description:

The Procera Client Management System (PCMS) software is intended to be used by dental practitioners and dental labs. Through the use of the PCMS software, a customer can import scanned data from the Procera scanners and define the shapes of Procera products through the use of a 3D-CAD tool. The user can then place and manage orders of the resulting Procera dental products.

Indications for Use:

- Nobel Biocare's Procera Software imports patient specific data from scanners and defines the shapes of dental prosthetic devices such as dental abutments, copings, laminates, and bridges through the use of a 3D-CAD tool. The software also serves as a means of ordering and managing orders of Procera products.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB 10 2006

Nobel Biocare AB
C/O Mr. Herbert Crane
Regulatory Affairs
Nobel Biocare USA, LLC
22715 Savi Ranch Parkway
Yorba Linda, California 92887

Re: K053602

Trade/Device Name: Procera Software
Regulation Number: 21 CFR 872.3630
Regulation Name: Endosseous Dental Implant Abutment
Regulatory Class: II
Product Code: NHA
Dated: December 22, 2005
Received: December 23, 2005

Dear Mr. Crane:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

1.2

Indications for Use

510(k) Number (if known):

Device Name: Procera Software

Indications For Use:

Nobel Biocare's Procera Software imports patient specific data from scanners and defines the shapes of dental prosthetic devices such as dental abutments, copings, laminates, and bridges through the use of a 3D-CAD tool. The software also serves as a means of ordering and managing orders of Procera products.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division 300-01)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Clinics

510(k) Number: K053602

Page 1 of 1