



510(k) Summary

Preparation Date: February 2, 2006

MAR 30 2006

Applicant/Sponsor: Biomet Manufacturing Corp.

Contact Person: Susan Alexander

Proprietary Name: Vanguard™ PS Open Box Porous Femoral Components

Common Name: Porous coated knee replacement components

Classification Name: Cemented semi-constrained polymer/metal/polymer knee prosthesis (888.3560); Knee joint patellofemoral tibial metal/polymer porous-coated uncemented prosthesis (888.3565)

Legally Marketed Devices To Which Substantial Equivalence Is Claimed: Maxim® Accel Knee System (K023546); Non-Cemented Porous Coated Knees (K033489)

Device Description: The Vanguard™ PS Open Box Porous Femoral Components described in this submission have the same articulating surface as the predicate Maxim® Accel (Vanguard™) posterior stabilized (PS) Interlok® femoral components and the cruciate retaining (CR) porous femoral components cleared in K023546 and feature the exact same porous-coated inner surface as the cruciate retaining (CR) porous design cleared in K023546 and K033489.

Intended Use: The Vanguard™ PS Open Box Porous Femoral Components are indicated for cemented or non-cemented use in cases of:

1. Painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, traumatic arthritis where one or more compartments are involved.
2. Correction of varus, valgus, or posttraumatic deformity.
3. Correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure.

The device is a single use implant.

Summary of Technologies: The technological characteristics (material, design, sizing, indications) of the Vanguard™ PS Open Box Porous Femoral Components are similar to or identical to the predicate devices.

Non-Clinical Testing: To supplement previous testing, an engineering justification including a Finite Element Comparison was performed. The results indicated that the device was functional within its intended use.

Clinical Testing: None provided as a basis for substantial equivalence.

*All trademarks are property of Biomet with the exception of the following:
Insall/Burstein II® is a trademark of Zimmer, Inc.
Tyvek® is a trademark of E.I. dePont de Nemours and Company*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAR 30 2006

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Biomet Manufacturing Corp.
c/o Ms. Susan Alexander
Regulatory Specialist
P.O. Box 587
Warsaw, Indiana 46581-0587

Re: K060303

Trade/Device Name: Vanguard™ PS Open Box Porous Femoral Components

Regulation Number: 21 CFR 888.3565

Regulation Name: Knee joint patellofemorotibial metal/polymer porous-coated uncemented prosthesis

Regulatory Class: Class II

Product Code: MBH, JWH

Dated: February 3, 2006

Received: February 6, 2006

Dear Ms. Alexander:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

Page 2 – Ms. Susan Alexander

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,


for Mark N. Melkerson
Director

Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

K060303

Device Name: Vanguard™ PS Open Box Porous Femoral Components

Indications For Use:

The Vanguard™ PS Open Box Porous Femoral Components are indicated for cemented or non-cemented use in cases of:

1. Painful and disabled knee joint resulting from osteoarthritis, rheumatoid arthritis, traumatic arthritis where one or more compartments are involved.
2. Correction of varus, valgus, or posttraumatic deformity.
3. Correction or revision of unsuccessful osteotomy, arthrodesis, or failure of previous joint replacement procedure.

The device is a single use implant.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-the-Counter Use NO
(Part 21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
**Division of General, Restorative,
and Neurological Devices**

510(k) Number K060303