

K 060388

Airway Management Inc.
6116 North Central Expressway
Suite 605
Dallas, Texas 75206

Non-Confidential Summary of Safety and Effectiveness
May 7, 2006

Airway Management Inc.
6116 North Central Expressway
Suite 605
Dallas, Texas 75225

Tel- (972) 369-0978

Fax- (214) 691-3151

Official Contact	Darren Edward Henderson
Proprietary or Trade Name	TAP II
Common / Usual Name	Oral Appliance – anti snoring device.
Classification Name	Anti -Snoring device
Device:	TAP II
Predicate Device	Nellcor Puritan Bennett – TAP-K962516 Thornton Oral Appliance – TOA-K972061 OASYS Oral Airway System – K030440

Device Description:

The TAP II Anti-Snoring device is comprised of:

- ◆ Lower tray fitted over the lower teeth.
- ◆ Upper tray fitted over the upper teeth.
- ◆ Impression material
- ◆ Hook and Base mechanism to attach lower tray to upper tray.

Intended Use:

Indicated Use--	The TAP II is intended to reduce or alleviate night time snoring and mild to moderate obstructive sleep apnea, OSA.
Target population --	Adult patients 18 years and older
Environment of Use--	Home and sleep laboratories.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 17 2006

Mr. Darren Henderson
Quality Manager
Airway Management, Incorporated
6116 North Central Expressway, Suite 605
Dallas, Texas 75206

Re: K060388

Trade/Device Name: Tap II Anti-Snoring and Obstructive Sleep Apnea Oral Appli

Regulation Number: 872.5570

Regulation Name: Intraoral Devices for Snoring and Intraoral Devices for Snoring and
Obstructive Sleep Apnea

Regulatory Class: II

Product Code: LRK

Dated: May 10, 2006

Received: May 11, 2006

Dear Mr. Henderson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and
Radiological Health

Enclosure

SECTION 3

INDICATIONS FOR USE

510(k) Number: K060388 (To be assigned)

Device Name: TAP II

Intended Use: To reduce or alleviate night time snoring and mild to moderate obstructive sleep apnea (OSA).

Environment of use: Home and sleep laboratories

Disposable / Reusable: Single patient – multi – use

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Print Sign-Off)
Susan Runner, M.D., M.Sc.
Division of Anesthesiology, General Hospital,
Pain Control, Dental Devices

510(k) Number: K060388

Prescription Use ☒

or

Over-the-counter use ☐

(Per CFR 801.109)