

MAY 10 2006

## 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K060656.

<b>1. Submitted by:</b>	Sysmex America, Inc. One Nelson C. White Parkway Mundelein, IL 60060 Phone: (847) 996-4675; FAX: (847) 996-4655 Contact person: Nina Gamperling Date prepared: March 10, 2006
<b>2. Name of Device:</b>	<u>Trade or proprietary name:</u> Sysmex <sup>®</sup> XS, Automated Hematology Analyzer.  <u>Common name:</u> XS  <u>Classification name:</u> Automated Differential Cell Counter, Sysmex <sup>®</sup> XS-1000i and XS-800i (21 CFR 864.5220)
<b>3. Predicate Device:</b>	The Sysmex <sup>®</sup> XS, Automated Hematology Analyzer, is substantially equivalent to the Sysmex XE-2100 and Sysmex XT-Series, Automated Hematology Analyzers.
<b>4. Device Description:</b>	The XE-2000 and XT-Series are automated hematology analyzers previously cleared by the FDA.
<b>5. Intended Use:</b>	The Sysmex <sup>®</sup> XS is an automated hematology analyzer for <i>in vitro</i> diagnostic use in clinical laboratories.
<b>6. Substantial equivalence-similarities and differences</b>	The following table compares the XS with the predicate method.

## 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS (continued)

**Comparison Table to Predicate Method**

	<b>Sysmex XE-2100</b>	<b>Sysmex XS-Series</b>
	<b>Predicate</b>	<b>New Instrument</b>
<b>Intended Use</b>	The Sysmex™ XE-2100 is a multi-parameter hematology analyzer intended to classify formed elements in anti-coagulated blood. The XE-2100 can provide accurate and precise test results for up to 32 analysis parameters in whole blood.	The Sysmex® XS-Series is an automated hematology analyzer for <i>in vitro</i> diagnostic use in clinical laboratories.
<b>Methodology</b>	The XE-2100 performs hematology analyses using the following methods: RF/DC Detection Method, Sheath Flow DC Detection Method, Flow Cytometry Methods using a Semiconductor Laser and SLS-hemoglobin method.	The XS performs hematology analyses using the following methods: Sheath Flow DC Detection Method, Flow Cytometry Methods using a Semiconductor Laser and SLS-hemoglobin method.
<b>Type of Anticoagulant</b>	EDTA	EDTA
<b>Specimen Type</b>	Peripheral blood	Peripheral blood
<b>Accuracy</b>	Performance was established in the previous 510(k) submission.	Comparison to the XE-2100 and XT-Series analysers demonstrated excellent correlation.

### 7. Clinical Performance Data:

Studies were performed to evaluate the equivalency of XS to the predicate method.

### 8. Conclusions:

The performance data demonstrated substantial equivalence.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Nina M. Gamperling, MBA, MT (ASCP), RAC  
Manager, Regulatory Affairs  
Sysmex America, Inc.  
One Nelson C. White Parkway  
Mundelein, Illinois 60060

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

MAR 10 2006

Re: k060656  
Trade/Device Name: Sysmex® XS, Automated Hematology Analyzer  
Regulation Number: 21 CFR § 864.5220  
Regulation Name: Automated differential cell counter  
Regulatory Class: II  
Product Code: GKZ  
Dated: March 10, 2006  
Received: March 13, 2006

Dear Ms. Gamperling:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 –

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,

A handwritten signature in cursive script, reading "Robert L. Becker, Jr.", written in dark ink.

Robert L. Becker, Jr., MD, Ph.D

Director

Division of Immunology and Hematology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K060656

Device Name: Sysmex® XS™ Automated Hematology Analyzer

### Indications For Use:

The Sysmex® XS is an automated hematology analyzer for *in vitro* diagnostic use in clinical laboratories.

(PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CHRD, Office of Device Evaluation (ODE)

Prescription Use ☒

OR

Over-The-Counter Use ☐

Josephine Bantadi  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

Sysmex America, Inc.  
510(k) XS Submission

510(k) K060656