

K060705

**510(k) Summary**  
acc. to 807.92

**JUN 14 2006**

**Submitter's Name and Address:** Dräger Medical b.v.  
Kanaaldijk 29  
5683CR Best  
The Netherlands

**Contact Person:** Mr. Hans Venings  
Vice President Processes, Quality and Regulatory Affairs

Phone: +31 499 331331  
Fax: +31 499 331222

**Applicants US Contact Person:** Ms Monica Ferrante  
Director, Regulatory Affairs

Phone: (215) 721-5400  
Fax: (215) 721-5425

**Date submission was prepared:** 2006-03-10

**Device Name:**

Common Name:	Ventilator
Classification Name:	Ventilator, Continuous
Regulation Number:	21 CFR 868.5895
Class:	II

**Legally Marketed Device Identification:** Carina home

**Device Description:**

the Carina home is a pressure controlled and volume controlled (pressure regulated) ventilator. The device is intended for use with patients with a tidal volume of 100 mL upwards. It is used for patients with respiratory insufficiency, using room air (also with additional oxygen).

**Intended Use:**

The Intended use of the Carina home is round-the-clock pressure controlled and pressure regulated volume controlled ventilation of patients with respiratory insufficiency, using room air, also with additional oxygen. The Carina home is intended for use by people without medical qualifications (patients and caregivers) as well as qualified medical and technical personnel. The device is intended for use with patients with a tidal volume of 100 mL upwards. It can be used invasive and non-invasive, for both controlled and assisted ventilation. The device can be used in Home Care environments, in and out of hospital, in stationary position (lying or sitting) or in a vehicle (e.g. car).

type	release status	effective date	number	organization	page/of
TEMPLATE	RELEASED	30.09.2004	DMS PQ2160 A4	Dräger Medical	1/2

**Predicate Devices:**

510(k) Number	Device Name
K052554	iVent 201 Portable Ventilator
K040790	LTV 1000 Ventilator
K033008	Air Safety HEPA and non-HEPA Filters
K003068	Savina
K970839	Siemens Servo Ventilator 300A

**Substantial Equivalence:**

The Carina home is similar to the homecare features to the iVent 201 (K052554) and the LTV 1000 (K040790).

It incorporates a pressure regulated volume controlled ventilation mode like the PRCV mode of the Siemens Servo Ventilator 300A (K970839).

It uses an integrated HEPA filter similar to HEPA filters from Air Safety Limited (K033008).

The Carina home and the iVent 201 are both software controlled ventilators using the same operating principle and can be used for invasive as well as non invasive ventilation.

The performance data of both ventilators is comparable, with the exception of some parameters. However, for these parameters, performance data is comparable with the LTV 1000 ventilator (K040790)

**Summary of Performance Testing:**

Safety testing was conducted per IEC60601-1, IEC60601-1-2, ASTM F1246 and other applicable standards with respect to mechanical, electrical and biocompatibility.

The results of all verification and validation testing demonstrate that all system and design requirements for the Carina home device have been met.

Qualification included hazard analysis, system level qualification and verification / validation tests.

type	release status	effective date	number	organization	page/of
TEMPLATE	RELEASED	30.09.2004	DMS PQ2160 A4	Dräger Medical	2/2



JUN 14 2006

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Dräger Medical, B.V.  
C/O Ms. Monica Ferrante  
Director, Regulatory Affairs  
Dräger Medical, Incorporated  
3135 Quarry Road  
Telford, Pennsylvania 18969

Re: K060705  
Trade/Device Name: Carina™Home  
Regulation Number: 868.5895  
Regulation Name: Continuous Ventilator  
Regulatory Class: II  
Product Code: CBK  
Dated: March 10, 2006  
Received: March 16, 2006

Dear Ms. Ferrante:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", is positioned above the printed name.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Use of this document / contents is forbidden without express authority. All rights reserved. Refer to protection notice ISQ 16016. Dräger Medical.

INDICATIONS FOR USE STATEMENT


Indications for Use

510(k) Number (if known): \_\_\_\_\_

Device Name: Carina<sup>TM</sup>home

Indications for Use:

the Carina *home* is a pressure controlled and volume controlled (pressure regulated) ventilator. The device is intended for use with patients with a tidal volume of 100 mL upwards. It is used for patients with respiratory insufficiency, using room air (also with additional oxygen). Federal law restricts this device to sale by or on the order of a physician.

  
Ann Suleman  
(Non Sign-Off)  
Department of Anesthesiology, General Hospital,  
Pain Control, Dental Devices  
(510) Number: K060705

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF  
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

template / rev. - / -	CVi HC USA Indicated Use Statement Carina home USA	reference STD321672	revision 2	page / of 3 / 3
--------------------------	---	------------------------	---------------	--------------------