

K060746

510 (k) Summary

(As required by 21 CFR 807.92 and 21 CFR 807.93)

MAY 25 2006

NAME OF SPONSOR: DePuy Orthopaedics, Inc.
700 Orthopaedic Drive
Warsaw, Indiana 46581-0988
Establishment Registration Number: 1818910

510(K) CONTACT: Rhonda Myer
Regulatory Affairs
Telephone: (574) 371-4927
Facsimile: (574) 371-4987
Electronic Mail: Rmyer7@dpyus.jnj.com

DATE PREPARED: February 23, 2006

PROPRIETARY NAME: Memory™ Staple

COMMON NAME: Bone Fixation Staple

CLASSIFICATION: Class II device per 21 CFR 888.3040: Smooth or threaded metallic bone fixation fastener

DEVICE PRODUCT CODE: 87 JDR

SUBSTANTIALLY EQUIVALENT DEVICES: DePuy Memory Staple (Size 12), K964226, cleared on August 5, 1997

DEVICE DESCRIPTION:

The Memory Staple is a single-use, bone fixation appliance intended to be implanted in the foot. It is a bicortical compression staple manufactured from a Nickel-Titanium alloy.

INTENDED USE AND INDICATIONS:

Intended Use:

The Memory Staple is intended to be implanted for bone fixation in the foot.

Indications for Use:

Memory 12 staples are indicated for osteotomies of the first phalanx of the foot.
Memory 20 staples are indicated for arthrodesis of the first metatarsal phalangeal joint.

BASIS OF SUBSTANTIAL EQUIVALENCE:

The substantial equivalence of the DePuy Memory Staple is substantiated by its similarity in intended use, indications for use, materials and design to the existing DePuy Memory Staple, cleared in K964226, on August 5, 1997.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 25 2006

DePuy Orthopaedics, Inc
% Ms. Rhonda Myer
Regulatory Affairs Associate
700 Orthopaedic Drive
Warsaw, Indiana 46581-0988

Re: K060746

Trade/Device Name: Memory™ Staple
Regulation Number: CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: JDR
Dated: March 17, 2006
Received: March 20, 2006

Dear Ms. Myer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Rhonda Myer

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



DePuy Orthopaedics, Inc.

PO Box 988
700 Orthopaedic Drive
Warsaw, Indiana 46581-0988
USA

Tel: +1 (574) 267 8143

Indications for Use Statement

510 (k) Number (if known): K060746

Device Name: _____

Indications for Use:

Memory 12 and Memory 20 staples are implants intended for fixation of the foot.

Memory 12 staples are indicated for osteotomies of the first phalanx of the foot.


Memory 20 staples are indicated for arthrodesis of the first metatarsal phalangeal joint.

The patient's anatomy and skeleton must be capable of receiving the selected implant.

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use: _____ x _____ OR Over-The-Counter-Use: _____

(Please do not write below this line – continue on another page if necessary)


(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

510(k) Number K060746