

K060979

Special 510(k) Summary of Safety and Effectiveness:

Line Extension to the Xia[®] and Xia[®] 4.5 Spinal Systems

Submission Information

Name and Address of the Sponsor of the 510(k) Submission: Stryker Spine
2 Pearl Court
Allendale, NJ 07401

Contact Person: Simona Voic
Regulatory Affairs Project Manager
Stryker Spine
2 Pearl Court, Allendale, NJ 07401
Tel: (201) 760 - 8145

Date of Summary Preparation: April 6, 2006

Device Identification

Proprietary Name: Xia[®] Titanium Spinal System and Xia[®] 4.5 Spinal System

Common Name: Spinal Fixation Appliances

Classification Name and Reference: Spinal Interlaminar Fixation Orthosis,
21 CFR §888.3050
Spinal Intervertebral Body Fixation Orthosis
21 CFR §888.3060
Pedicle Screw Spinal System
21 CFR §888.3070(b)(1) and (b)(2)

Device Product Code: NKB, KWP, KWQ, MNH, and MNI

01

K060979

Predicate Device Information:

K060361 – Stryker Spine Xia® Spinal System and Stryker Spine Xia® 4.5 Spinal System

K020709 – Medtronic Sofamor Danek CD Legacy 4.5 Spinal System

K043488 – Medtronic Sofamor Danek CD Horizon® Spinal System

Predicate Device Identification

The Stryker Spine Xia® Spinal System consists of Ø6mm rods, Monoaxial and Polyaxial screws, Hooks, Blockers, and Monoaxial and Polyaxial Cross Connectors. The Xia® 4.5 Spinal System is comprised of Ø 4.5 mm rods, Polyaxial and Monoaxial bone screws, Blockers, Hooks, Dual Staples, and Connectors. The components of both systems are available in a variety of diameters and lengths in order to accommodate patient anatomy and are fabricated from titanium alloy. The implants are provided non-sterile and are used for either posterior or anterior/anterolateral non-cervical spinal fixation.

Description of Device Modification

This submission is intended to address a line extension to both the Xia® Titanium and Xia® 4.5 Spinal Systems (K060361). The line extension for the Xia® Titanium Spinal System includes:

- A new Ø6mm spinal rod fabricated from Cobalt-Chromium-Molybdenum Alloy (Vitallium®).

The line extension for the Xia® 4.5 Spinal System includes:

- New Ø4mm monoaxial titanium alloy screw in lengths of 20mm to 40mm in five millimeter increments,
- New Ø4mm polyaxial titanium alloy screw in lengths of 20mm to 40mm in five millimeter increments, and
- New Ø 4.5mm spinal rod fabricated from Cobalt-Chromium-Molybdenum Alloy (Vitallium®).

Intended Use:

The Xia® Spinal System and Xia® 4.5 Spinal System are intended for anterior/anterolateral and posterior, noncervical pedicle and non-pedicle fixation for the following indications: degenerative disc disease (DDD) (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion.

The 6 mm diameter rods from the DIAPASON™ Spinal System and OPUS™ Spinal System are intended to be used with the other components of the Xia® Titanium Spinal System. The Titanium Multi-Axial Cross Connectors are intended to be used with the other components of the Xia® Titanium Spinal System.

Statement of Technological Comparison:

The subject components share the same intended use and basic design concepts as that of the predicate devices: Xia® Spinal System (K060361) and the Xia® 4.5 Spinal System (K060361). Mechanical testing also demonstrated comparable mechanical properties to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 19 2006

Stryker Spine
% Ms. Simona Voic
Regulatory Affairs Project Manager
2 Pearl Court
Allendale, New Jersey 07401

Re: K060979

Trade/Device Name: Stryker Spine Xia® 4.5 Spinal System and Xia® 4.5 Spinal System

Regulation Number: 21 CFR 888.3070

Regulation Name: Pedicle screw spinal system

Regulatory Class: Class III

Product Code: NKB, KWP, KWQ, MNH, MNI

Dated: May 9, 2006

Received: May 10, 2006

Dear Ms. Voic:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Ms. Simona Voic

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,


Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: _____

Indications For Use:

Prescription Use _____ AND/OR Over-The-Counter Use _____
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Indications for Use

510(k) Number (if known): K 060979

Device Name: Stryker Spine Xia® Spinal System and Xia® 4.5 Spinal System

Indications For Use:

The Xia® Spinal System and Xia® 4.5 Spinal System are intended for anterior/anterolateral and posterior, noncervical pedicle and non-pedicle fixation for the following indications: degenerative disc disease (DDD) (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion.

The 6 mm diameter rods from the DIAPASON™ Spinal System and OPUS™ Spinal System are intended to be used with the other components of the Xia® Titanium Spinal System. The Titanium Multi-Axial Cross Connectors are intended to be used with the other components of the Xia® Titanium Spinal System.

Prescription Use X
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1


(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

510(k) Number K060979