

K061202

### III. 510(K) Summary

**SUBMITTED BY:**

Globus Medical Inc.  
303 Schell Lane  
Phoenixville, PA 19460  
(610) 415-9000 x218  
Contact: Kelly J. Baker

JUL 20 2006

**DEVICE NAME:**

REVERE™ Stabilization System

**CLASSIFICATION:**

Per 21 CFR as follows:

- §888.3050 Spinal Interlaminar Fixation Orthosis
- §888.3060 Spinal Intervertebral Body Fixation Orthosis
- §888.3070 Pedicle Screw Spinal System
- §888.3070 Spondylolisthesis Spinal Fixation Device System

Product Codes MNH, MNI, KWQ, KWP, NKB.

Regulatory Class II and III.

Panel code 87.

**PREDICATE DEVICES:**

Globus Medical PROTEX Stabilization System:

PROTEX K040442, SE date May 20, 2004

PROTEX 6.0mm rods, K052069, SE date August 17, 2005

Stryker XIA Stainless Steel System

K012870, SE date September 24, 2001

Line extensions K031090, SE date April 24, 2003

Line extensions K053115, SE date December 6, 2005

Product Codes MNH, MNI, KWQ, KWP, NKB.

Regulatory Class II and III.

**DEVICE DESCRIPTION:**

The REVERE™ Stabilization Systems consist of a variety of shapes and sizes of rods, hooks, monoaxial screws, polyaxial screws, locking caps, t-connectors, staples, and associated manual surgical instruments. Implant components can be rigidly locked into a variety of configurations for the individual patient and surgical condition. Polyaxial screws, hooks, and t-connectors are intended for posterior use only. Staples are intended for anterior use only. Rods and monoaxial screws may be used anteriorly or posteriorly. Locking caps are used to connect screws or hooks to the rod.

The rods are composed of titanium alloy, commercially pure titanium, or stainless steel, as specified in ASTM F136, F1295, F67, and F138. All other implants are composed of titanium alloy or stainless steel, as specified in ASTM F136, F1295, and F138.

**INTENDED USE:**

The REVERE™ Stabilization System, when used as posterior pedicle screw systems, are intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine: degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), degenerative spondylolisthesis with objective evidence of neurologic impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, pseudoarthrosis and failed previous fusion.

In addition, the REVERE™ Stabilization Systems are intended for treatment of severe spondylolisthesis (Grades 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft, having implants attached to the lumbosacral spine and/or ilium with removal of the implants after attainment of a solid fusion. Levels of pedicle screw fixation for these patients are L3-sacrum/ilium.

When used as posterior non-pedicle screw fixation systems, the REVERE™ Stabilization Systems are intended for the treatment of degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), spinal stenosis, spondylolisthesis, spinal deformities (i.e. scoliosis, kyphosis, and/or lordosis, Scheuermann's disease), fracture, pseudarthrosis, tumor resection, and/or failed previous fusion. Overall levels of fixation are T1-sacrum/ilium.

When used as anterolateral thoracolumbar systems, the REVERE™ Stabilization Systems are intended for anterolateral screw (with or without staple) fixation for the following indications: degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), spinal stenosis, spondylolisthesis, spinal deformities (i.e. scoliosis, kyphosis, and/or lordosis), fracture or dislocation of the thoracolumbar spine, pseudoarthrosis, tumor resection, and/or failed previous fusion. Levels of screw fixation are T8-L5.

**PERFORMANCE DATA:**

Mechanical testing in accordance with the "Guidance for Industry and FDA Staff, Guidance for Spinal System 510(k)s", May 3, 2004 is presented.

**BASIS OF SUBSTANTIAL EQUIVALENCE:**

The REVERE™ Stabilization System implants are similar to the predicate PROTEX™ (K040442, K052069) and Stryker Xia (K012870, K031090, K053115) thoracolumbar stabilization system with respect to technical characteristics, performance, and intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUL 20 2006

Globus Medical, Inc.  
% Ms. Kelly Baker  
Director, Regulatory and Clinical Affairs  
303 Schell Lane  
Phoenixville, Pennsylvania 19460

Re: K061202

Trade/Device Name: REVERE™ Stabilization System  
Regulation Number: 21 CFR 888.3070, 21 CFR 888.3060, 21 CFR 888.3050  
Regulation Name: Spinal Interlaminar Fixation Orthosis  
Regulatory Class: Class III  
Product Code: NKB, MNI, MNH, KWQ, KWP  
Dated: June 20, 2006  
Received: June 21, 2006

Dear Ms. Baker:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2- Ms. Kelly Baker

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a small "for" written below it.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**II. Indications for Use Statement**

510(k) Number: \_\_\_\_\_

Device Name: REVERE™ Stabilization System**Indications:**

The REVERE™ Stabilization System, when used as a posterior pedicle screw system, is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine: degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), degenerative spondylolisthesis with objective evidence of neurologic impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, pseudoarthrosis and failed previous fusion.

In addition, the REVERE™ Stabilization System is intended for treatment of severe spondylolisthesis (Grades 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft, having implants attached to the lumbosacral spine and/or ilium with removal of the implants after attainment of a solid fusion. Levels of pedicle screw fixation for these patients are L3-sacrum/ilium.

When used as a posterior non-pedicle screw fixation system, the REVERE™ Stabilization System is intended for the treatment of degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), spinal stenosis, spondylolisthesis, spinal deformities (i.e. scoliosis, kyphosis, and/or lordosis, Scheuermann's disease), fracture, pseudarthrosis, tumor resection, and/or failed previous fusion. Overall levels of fixation are T1-sacrum/ilium.

When used as an anterolateral thoracolumbar system, the REVERE™ Stabilization System is intended for anterolateral screw (with or without staple) fixation for the following indications: degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), spinal stenosis, spondylolisthesis, spinal deformities (i.e. scoliosis, kyphosis, and/or lordosis), fracture or dislocation of the thoracolumbar spine, pseudoarthrosis, tumor resection, and/or failed previous fusion. Levels of screw fixation are T8-L5.

Prescription Use   X   OR Over-The-Counter Use       
(Per 21 CFR §801.109)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Barbara Buckley, MD  
(Division Sign-Off) \_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K061202