

COOK®

K061670

William Cook Europe ApS  
Sandet 6, 4632 Bjaeverskov  
DENMARK  
Phone: +45 56 86 86 86  
Fax: +45 56 86 86 96  
CVR No. 83 74 23 13  
www.cookgroup.com

JAN 19 2007

**510(k) Summary**

**Submitted By:**

Tina H. Andersen  
WILLIAM COOK EUROPE APS  
Sandet 6,  
DK-4632 Bjaeverskov,  
DENMARK  
+45 56 86 87 60

June 09, 2006

**Device:**

Trade Name: Lunderquist Wire Guide/Lunderquist DC Wire Guide  
Proposed Classification: 870.1330 DQX  
Class II, Cardiovascular

**Predicate Devices:**

The *Lunderquist Wire Guide* and *Lunderquist DC Wire Guide* are similar in terms of intended use, materials of construction and technological characteristics to predicate devices designed for diagnostic and interventional procedures.

**Device Description:**

The *Lunderquist Wire Guide* is a PTFE-coated stainless steel wire guide with an outer diameter of .035 inches and comes in 90 to 300 cm lengths. For lengths 260-300 cm, a radiopaque gold coil is provided on the flexible distal tip. The distal tip has 4 or 7 cm of flexibility comes straight or as a J-curve with a 3 or 7.5 mm radius.

The *Lunderquist DC Wire Guide* is a PTFE-coated stainless steel wire guide with an outer diameter of .035 inches and comes in 260-300 cm lengths. A radiopaque gold coil is provided on the flexible distal tip. The *Lunderquist DC Wire Guide* is a PTFE-coated stainless steel wire guide with a double curved tip design. The double curved distal tip has 4 cm of tip flexibility.

**Substantial Equivalence:**

These devices will be manufactured according to specified process controls and a Quality Assurance Program. These devices will undergo packaging similar to the devices currently marketed and distributed by William Cook Europe ApS. These devices will undergo sterilization similar to the devices currently marketed and distributed. Being similar with respect to indications for use, materials and physical construction to predicate devices, these devices meets the requirements for section 510(k) substantial equivalence.

**Test Data:**

The *Lunderquist Wire Guide* and *Lunderquist DC Wire Guide* were subjected to the following tests to assure reliable design and performance under the specified testing parameters. These tests were comprised of:

1. Tensile Testing
2. Flexible Testing
3. Functional Testing
4. Corrosion Testing
5. Radiodetectability

The results of these tests provide reasonable assurance that the devices have been designed and tested to assure conformance to the requirements for their use as guide wires.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

William Cook Europe ApS  
c/o Ms. Tina H. Andersen  
Regulatory Affairs Coordinator  
Sandet 6, 4632 Bjaeverskov  
DENMARK

JAN 19 2007

Re: K061670

Trade/Device Name: Lunderquist Wire Guide and Lunderquist DC Wire Guide  
Regulation Number: 21 CFR 870.1330  
Regulation Name: Catheter guide wire  
Regulatory Class: Class II  
Product Code: DQX  
Dated: January 8, 2007  
Received: January 10, 2007

Dear Ms. Andersen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Tina H. Andersen

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Bram D. Zuckerman, M.D.

Director

Division of Cardiovascular Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

### Instructions for Use

510(k) Number (if known): K061670

Device Name: Lunderquist Wire Guide and Lunderquist DC Wire Guide

**Indications for Use:**

The Lunderquist Wire Guide and the Lunderquist DC Wire Guide are intended for complex diagnostic and interventional procedures where increased body, flexibility, and low surface friction of the wire guide are needed.

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-the-Counter Use \_\_\_\_\_

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED:

Concurrence of CDRH, Office of Device Evaluation (ODE)

Diana R. Volney

Division Sign-Off)

Division of Cardiovascular Devices

510(k) Number K061670