

K061946

510(K) SUMMARY

JUL 21 2006

Siemens Acuson X300 Diagnostic Ultrasound system

This summary of safety and effectiveness is provided as part of this Premarket Notification in compliance with the Safe Medical Device Act of 1990 revisions to 21 CFR, Part 807.92, Content and Format of a 510(k) Summary.

1. **Submitted By:**

Siemens Medical Solutions USA, Inc., Ultrasound Division
1230 Shorebird Way
Mountain View, California 94043

Contact Person:

Sheila W. Pickering
Regulatory Affairs
Phone (650) 943 7187
Fax (650) 943 7053

Date Prepared:

June 16, 2006

2. **Proprietary Names:**

Siemens Acuson X300™ Diagnostic Ultrasound System
Sonoline Premier Plus Diagnostic Ultrasound System
Sonavista X 300 Diagnostic Ultrasound System

Common/ Usual Name:

Diagnostic Ultrasound System with Accessories

Classification Name:

21 CFR 892.1550

Ultrasonic Pulsed Doppler Imaging System	FR # 892.1550	Product Code 90-IYN
Ultrasonic Pulsed Echo Imaging System	FR # 892.1560	Product Code 90-IYO
Diagnostic Ultrasound Transducer	FR # 892.1570	Product Code 90-ITX

3. **Predicate Device:**

- K050034, Antares Diagnostic Ultrasound System
- K052894, Sonoline G60 S Diagnostic Ultrasound System
- K043016, Sonoline Orchid Diagnostic Ultrasound System
- K042044, Acuson CV70 Cardiovascular System

4. **Device Description:**

The Siemens Acuson X300 is a general purpose, mobile, software-controlled, diagnostic ultrasound system with an on-screen display for thermal and mechanical indices related to potential bioeffect mechanisms. Its function is to acquire primary or secondary harmonic ultrasound echo data and display it in B-Mode, M-Mode, Pulsed (PW) Doppler Mode, Continuous (CW) Doppler Mode, Color Doppler Mode, Amplitude Doppler Mode, a combination of modes, 3D Imaging, or Harmonic Imaging and 4D imaging on a FPD(Flat Panel Display) display.

The Siemens Acuson X300, has been designed to meet the following product safety standards:

- UL 60601-1, Safety Requirements for Medical Equipment
- CSA C22.2 No. 601-1, Safety Requirements for Medical Equipment
- AIUM/NEMA, 1998, Standard for Real Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment

- AIUM/NEMA UD-2, 1998 Acoustic Output Measurement Standard for Diagnostic Ultrasound
- 93/42/EEC Medical Devices Directive
- Safety and EMC Requirements for Medical Equipment
 - EN 60601-1
 - EN 60601-1-1
 - EN 60601-1-2
- IEC 61157 Declaration of Acoustic Power
- ISO 10993 Biocompatibility

5. Intended Uses:

The Siemens Acuson X300 ultrasound imaging system is intended for the following applications: General Radiology, Fetal, Abdominal, Intraoperative, Pediatric, Small Parts, Neonatal/Adult Cephalic, Cardiac, Transesophageal, Pelvic, Transcranial, OB/GYN, Urology, Vascular, Musculoskeletal, Superficial Musculoskeletal, and Peripheral Vascular applications.

The system also provides for the measurement of anatomical structures and for analysis packages that provide information that is used for clinical diagnosis purposes.

6. Technological Comparison to Predicate Device:

The Siemens Acuson X300 is substantially equivalent to the Sonoline Antares (K050034), Sonoline G 60 S (K052894), Sonoline Orchid (K043016) and the Acuson CV70 (K042044). All systems transmit ultrasonic energy into patients, then perform post processing of received echoes to generate on-screen display of anatomic structures and fluid flow within the body. All systems allow for specialized measurements of structures and flow, and calculations.

End of 510(k) Summary



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 21 2006

Siemens Medical Solutions USA, Inc. Ultrasound Group
% Mr. Mark Job
Responsible Third Party Official
Regulatory Technology Services LLC
1394 25th Street NW
BUFFALO MN 55313

Re: K061946

Trade Name: ACUSON X300 Diagnostic Ultrasound System
Regulation Number: 21 CFR 892.1550
Regulation Name: Ultrasonic pulsed doppler imaging system
Regulation Number: 21 CFR 892.1560
Regulation Name: Ultrasonic pulsed echo imaging system
Regulation Number: 21 CFR 892.1570
Regulation Name: Diagnostic ultrasonic transducer
Regulatory Class: II
Product Code: IYN, IYO, and ITX
Dated: July 6, 2006
Received: July 10, 2006

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the ACUSON C300 Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Number

P4-2
CH5-2
VF10-5
L9-5
EC9-4
EV9-4
VF13-5
P8-4

BE9-4
CW2
CW5
Acu Nav 8F
Acu Nav 10F
V5Ms TEE
4V1c
VF13-5SP

C8-5
8L3
10V4
C7F2
EV9F4
L13F5 3D/4D

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:


Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Boulevard
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

If you have any questions regarding the content of this letter, please contact Andrew Kang, M.D. at (301) 594-1212.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure(s)

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Abdominal		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Intraoperative Neurological		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Pediatric		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Cardiac		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Transesophageal		N	N	N	N	N	N		BMDC	Note 2,3,7,8,9
Transrectal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Transvaginal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Transurethral										
Intravascular		N	N	N	N	N	N		BMDC	Note 2,3,7,8,9,10
Peripheral vessel		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Other (specify)		N	N	N	N	N	N		BMDC	Note 2,3,7,8,9,10

N = new indication; P = previously cleared by FDA; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
- Note 2 Ensemble tissue harmonic imaging
- Note 3 3D imaging
- Note 4 B&W SieScape panoramic imaging
- Note 5 Power SieScape panoramic imaging
- Note 6 For example: abdominal, vascular
- Note 7 Contrast agent imaging
- Note 8 SieClear multi-view spatial compounding
- Note 9 Tissue Equalization Technology
- Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Diagnostic Ultrasound Indications for Use Form

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

510(k) Number (if known):

K061946

Device Name:

P4-2 Phased Sector Array Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P	P	P	P		BMDC	Note 2,3,5,6,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)										
Neonatal Cephalic		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Cardiac		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new Indication; P = previously cleared under K043016; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

CH5-2 Convex Array Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
- Note 2 Ensemble tissue harmonic imaging
- Note 3 3D imaging
- Note 4 B&W SieScape panoramic imaging
- Note 5 Power SieScape panoramic imaging
- Note 6 For example: abdominal, vascular
- Note 7 Contrast agent imaging
- Note 8 SieClear multi-view spatial compounding
- Note 9 Tissue Equalization Technology
- Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

VF10-5 Linear Array Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D Imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

L9-5 Linear Array Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
- Note 2 Ensemble tissue harmonic imaging
- Note 3 3D imaging
- Note 4 B&W SieScape panoramic imaging
- Note 5 Power SieScape panoramic imaging
- Note 6 For example: abdominal, vascular
- Note 7 Contrast agent imaging
- Note 8 SieClear multi-view spatial compounding
- Note 9 Tissue Equalization Technology
- Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**EC9-4 Convex Array Endocavity Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transvaginal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transurethral										
Intravascular										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**EV9-4 Convex Array Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transvaginal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transurethral										
Intravascular										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D Imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
 (Division Sign-Off)
 Division: Reproductive, Abdominal,
 and Radiological Devices

510(k)

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K 061946

Device Name:

VF13-5 Linear Array Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Abdominal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D Imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K 061946

Device Name:

P8-4 Phase Array Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)										
Neonatal Cephalic		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K043016(P7-4); E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D Imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K 061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**BE9-4 Convex Array Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transvaginal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transurethral										
Intravascular										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K043016; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**CW2 Continuous Wave Doppler Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal					P					
Abdominal					P					
Intraoperative (Note 6)					P					
Intraoperative Neurological										
Pediatric					P					
Small Organ (Note 1)					P					
Neonatal Cephalic					P					
Adult Cephalic					P					
Cardiac					P					
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel					P					
Laparoscopic										
Musculo-skeletal Conventional					P					
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K001400; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
 (Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Radiological Devices

510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**CW5 Continuous Wave Doppler Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal					P					
Abdominal					P					
Intraoperative (Note 6)					P					
Intraoperative Neurological										
Pediatric					P					
Small Organ (Note 1)					P					
Neonatal Cephalic					P					
Adult Cephalic					P					
Cardiac					P					
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel					P					
Laparoscopic										
Musculo-skeletal Conventional					P					
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K023720; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

Acu Nav 8F Intracardiac Transducer for use with :

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9,10
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Intra-Cardiac)		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9,10

N = new indication; P = previously cleared under K042593; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**Acu Nav 10F Intracardiac Transducer for use with :
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9,10
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Intra-Cardiac)		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9,10

N = new indication; P = previously cleared under K033650; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number *K061946*

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

V5Ms TEE Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal		P	P	P	P	P	P		BMDC	Note 2,3,7,8,9
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K022567; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

4V1c Phased Array Transducer for use with :
ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P	P	P	P		BMDC	Note 2,3,5,6,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)										
Neonatal Cephalic		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Cardiac		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K022567; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D Imaging
 Note 4 B&W SieScape panoramic Imaging
 Note 5 Power SieScape panoramic Imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
 (Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Radiological Devices
 510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**VF13-5SP Linear array Transducer for use with :
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intraoperative (Note 6)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative Neurological		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Other (specify)										

N = new indication; P = previously cleared under K023720; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D Imaging
 Note 4 B&W SieScape panoramic Imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
 (Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Radiological Devices

510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**C8-5 Tight Curved Array Transducer for use with:
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,5,6,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P			
Neonatal Cephalic		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Cardiac		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P			
Musculo-skeletal Superficial		N	N	N		N	N			
Other (specify)										

N = new indication; P = previously cleared under K020353; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Shirley C. Kaplan
 (Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number *K061946*

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

8L3 Linear "Regel" Transducer for use with :
ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Other (specify)										

N = new indication; P = previously cleared under K020353(5.0L45); E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

10V4 Phased Array Neonatal High Frequency Transducer for use with:

ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)										
Neonatal Cephalic		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac		P	P	P	P	P	P		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		N	N	N	N	N	N		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K022567; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D Imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K 061946

Device Name:

C7F2 Curved array mechanical 3D/4D Transducer for use with :
ACUSON X300 Diagnostic Ultrasound Systems

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K023720; E = added under Appendix E

Note 1 For example: breast, testes, thyroid, penis, prostate, etc.

Note 2 Ensemble tissue harmonic imaging

Note 3 3D imaging

Note 4 B&W SieScape panoramic imaging

Note 5 Power SieScape panoramic imaging

Note 6 For example: abdominal, vascular

Note 7 Contrast agent imaging

Note 8 SieClear multi-view spatial compounding

Note 9 Tissue Equalization Technology

Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, and Radiological Devices

510(k) Number

K061946

Diagnostic Ultrasound Indications for Use Form

510(k) Number (if known):

K061946

Device Name:

**EV9F4 Curved array mechanical 3D/4D Transducer for use with :
ACUSON X300 Diagnostic Ultrasound Systems**

Intended Use:

Diagnostic imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal										
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric										
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal		N	N	N		N	N		BMDC	Note 2,3,4,5,7,8,9
Transvaginal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Transurethral										
Intravascular										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (specify)										

N = new indication; P = previously cleared under K050034; E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
 Note 2 Ensemble tissue harmonic imaging
 Note 3 3D imaging
 Note 4 B&W SieScape panoramic imaging
 Note 5 Power SieScape panoramic imaging
 Note 6 For example: abdominal, vascular
 Note 7 Contrast agent imaging
 Note 8 SieClear multi-view spatial compounding
 Note 9 Tissue Equalization Technology
 Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

[Signature]
 (Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Radiological Devices
 510(k) Number *K061946*

Diagnostic Ultrasound Indications for Use Form

K061946

Device Name: **L13F5 3D/4D mechanical wobbler linear transducer for use with:**
ACUSON X300 Diagnostic Ultrasound Systems
 Intended Use: **Diagnostic imaging or fluid flow analysis of the human body as follows:**

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (Specify)	Other (Specify)
Ophthalmic										
Fetal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Abdominal		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Intraoperative (Note 6)										
Intraoperative Neurological										
Pediatric		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Small Organ (Note 1)		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Neonatal Cephalic		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral vessel		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Laparoscopic										
Musculo-skeletal Conventional		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Musculo-skeletal Superficial		P	P	P		P	P		BMDC	Note 2,3,4,5,7,8,9
Other (specify)										

N = new indication; P = previously cleared under K023720(C7F2 Wobbler) and K043016(VF13-5 Linear);
 E = added under Appendix E

- Note 1 For example: breast, testes, thyroid, penis, prostate, etc.
- Note 2 Ensemble tissue harmonic imaging
- Note 3 3D imaging
- Note 4 B&W SieScape panoramic imaging
- Note 5 Power SieScape panoramic imaging
- Note 6 For example: abdominal, vascular
- Note 7 Contrast agent imaging
- Note 8 SieClear multi-view spatial compounding
- Note 9 Tissue Equalization Technology
- Note 10 Intracardiac imaging

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Radiological Devices

510(k) Number

K061946