

510(k) Summary  
as required by 807.92

K062053

AUG - 7 2006

1. Company Identification

EIZO NANA CORPORATION

153 Shimokashiwano-cho, Hakusan-shi, Ishikawa-ken, 924-8566, Japan

Tel: +81-76-274-2468

Fax: +81-76-274-2484

2. Official Correspondent

Hiroaki Hashimoto (Mr.)

Manager of Engineering Management Section

3. Date of Submission

June 28, 2006

4. Device Trade name

Monochrome LCD Monitor, RadiForce GS320

5. Common/Usual Name

Image display system, medical image workstation, image monitor/display, and others

6. Classification Number

Medical displays classified in Class II per 21 CFR 892.2050.

7. Predicate Device

Manufacturer : EIZO NANA CORPORATION

Device Name : Monochrome LCD Monitor

Model Name : RadiForce G310

510(k) No. : K060845

8. Description of Device

RadiForce GS320 is a 54cm (21.3") Monochrome LCD display for medical image viewing. GS310 displays high-definition medical imaging.

9. Intended Use

RadiForce GS320 is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device is not specified for digital mammography system.

10. Technological Characteristics

GS320 employs smaller grayscale tones than that of GS310. The only modification is that the panel size became big with 54cm(21.3") from 53cm(20.8"). Comparison table of the principal characteristics of 2 devices in Attachment 1 shows that new and predicate devices are substantially equivalent in the areas of technical characteristics, general functions.

## Appendix 1: Comparison Table with Predicate Device

Items	GS310	GS320
510(k) Number	K060845	Not known
Panel Size and Type	53 cm (20.8") TFT Monochrome LCD panel	54 cm (21.3") TFT Monochrome LCD panel
Cabinet Color	Black	It is same as the following.
Pixel Pitch	0.207 x 0.207mm	0.2115 x 0.2115mm
Ratio of Sub-pixel Opening Areas	4: 4: 4	It is same as the following.
Grayscale Tones	1,024 from a pallet of 8,161	It is same as the following.
Viewing Angles	H: 170°, V: 170°	It is same as the following.
Scanning Frequency (H, V)	31-100kHz, 48-71.5Hz (VGA Text: 69-71Hz) Frame synchronous mode: 59-61Hz	It is same as the following.
Native Resolutions	1536 x 2048	It is same as the following.
Brightness	700 cd/m <sup>2</sup> (Typical)	It is same as the following.
Contrast Ratio	900: 1 (typical)	850: 1 (typical)
DOT Clock	165 MHz	It is same as the following.
Response Time	50 ms (typical)	It is same as the following.
Active Display Size (H x V)	318 x 424 mm	324 x 433 mm
Viewable Image Size	529 mm (20.8") (diagonal)	541mm (21.3") (diagonal)
Luminance Calibration	Built-in swing calibration sensor provided.	It is same as the following.
Input Signals	DVI Standard 1.0	It is same as the following.
Input Terminals	DVU-D 24 pin	It is same as the following.
USB Ports / Standard	1 upstream, 2 downstream	1 upstream, 2 downstream/Standard Rev.2.0
Power	AC100-120V, 200-240V, 50/60Hz	It is same as the following.
Power Management	DVI-DMPM	It is same as the following.
Dimensions (W x H x D)	With Stand: 368 x 515.5 mm – 597.5 x 209 mm Without Stand: 368 x 486 x 90 mm	With Stand: 376 x 522.5 mm – 604.5x 208.5 mm Without Stand: 376 x 500 x 92 mm
Certifications & Standards	TUV/GM, CE Medical Device Directive, CB (EN60601-1), cTUVus (UL2601-1, CSA C22.2 No. 601-1), VCCI-B, FCC-B, Canadian ICES-003-A, CCC	TUV/GM, CE Medical Device Directive, CB (EN60601-1), cTUVus (UL2601-1, CSA C22.2 No. 601-1), VCCI-B, FCC-B, Canadian ICES-003-A, CCC



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
9200 Corporate Blvd.  
Rockville MD 20850

AUG - 7 2006

Mr. Hiraoki Hashimoto  
Manager of Engineering Management Section  
EIZO NANA O Corporation  
153 Shimokashiwano, Hakusan,  
Ishikawa 924-8566  
JAPAN

Re: K062053

Trade/Device Name: Monochrome LCD Monitor, RadiForce GS320  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and communications system  
Regulatory Class: II  
Product Code: LLZ  
Dated: July 3, 2006  
Received: July 20, 2006

Dear Mr. Hashimoto:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.



*Protecting and Promoting Public Health*

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

21 CFR 876.xxx	(Gastroenterology/Renal/Urology)	240-276-0115
21 CFR 884.xxx	(Obstetrics/Gynecology)	240-276-0115
21 CFR 894.xxx	(Radiology)	240-276-0120
Other		240-276-0100

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): ~~Not known~~

K062053

Device Name : Monochrome LCD Monitor, RadiForce GS320

## Indications For Use:

RadiForce GS320 is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners. The device is not specified for digital mammography system.

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

David M. Chapman

(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K062053