

SEP - 1 2006

K062307

1/1

Section 5 – 510(k) Summary

Applicant: Anulex Technologies, Inc.
5600 Rowland Road, Suite 280
Minnetonka, MN 55343

Contact Person: Tim Miller
Vice President, Regulatory and Clinical Affairs
Telephone: 952.224.4003
Fax: 952.224.4040
e-mail: tmiller@anulex.com

Date Prepared: August 7, 2006

Trade Name: Xclose™ Tissue Repair System

Product Classification and Code: 21 CFR §878.5000
Class: II
Product Code: GAT

Predicate Device: K061386 - Anulex Anchor Band Suturing System

Device Description: The Xclose™ Tissue Repair System consists of two (2) non-absorbable braided surgical 3-0 suture and T-anchor assemblies, connected together with a loop of green 2-0 suture. The 2-0 suture loop is used to facilitate tightening, drawing the 3-0 suture assemblies together, thereby re-approximating the tissue. The device construct is composed of polyethylene terephthalate (PET) and conforms to USP requirements. The construct is provided sterile and preloaded on a disposable delivery instrument.

Intended Use: The Xclose™ Tissue Repair System is indicated for use in soft tissue approximation for procedures such as general and orthopedic surgery.

Summary of Technological Characteristics: The modifications to the Anchor Band Suturing System were conducted in accordance with the Anulex Design Control System. Accordingly, the risk analysis identified necessary design verification and validation activities. As a result of this analysis, tensile testing was performed to confirm compliance to USP suture requirements.

Conclusion: The Xclose™ Tissue Repair System is substantially equivalent to the Anchor Band Suturing System in regards to the indications for use, the basic operating principle, materials, sterilization, packaging and shelf-life.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

SEP - 1 2006

Anulex Technologies, Inc.
% Mr. Tim Miller
VP Regulatory and Clinical Affairs
5600 Rowland Road, Suite 280
Minnetonka, Minnesota 55343

Re: K062307

Trade/Device Name: Xclose™ Tissue Repair System
Regulation Number: 21 CFR 878.5000
Regulation Name: Nonabsorbable poly(ethylene terephthalate) surgical suture
Regulatory Class: Class II
Product Code: GAT
Dated: August 7, 2006
Received: August 8, 2006

Dear Mr. Miller:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA

Page 2 - Mr. Tim Miller

finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Section 4 – Indications For Use Statement

510(k) Number (if known): K062307

Device Name: Xclose™ Tissue Repair System

Indications for Use:

The Xclose™ Tissue Repair System is indicated for use in soft tissue approximation for procedures such as general and orthopedic surgery.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

**Division of General, Restorative
and Neurological Devices**

510(k) Number K062307