

K062408 (p 1 of 1)

510(K) Summary  
Smith & Nephew Modular Femoral (Hemi) Heads

SEP 12 2006

SUBMITTER'S NAME: Smith & Nephew, Inc., Orthopaedic Division  
SUBMITTER'S ADDRESS: 1450 East Brooks Road, Memphis, TN 38116  
SUBMITTER'S TELEPHONE NUMBER: 901-399-6707  
CONTACT PERSON: Gino J. Rouss  
DATE SUMMARY PREPARED: August 15, 2006  
TRADE OR PROPRIETARY DEVICE NAME: Smith & Nephew Modular Femoral (Hemi) Heads  
COMMON OR USUAL NAME: Artificial Hip Replacement Prosthesis  
CLASSIFICATION NAME: Hip joint femoral (hemi-hip) metallic cemented or uncemented prosthesis, 21 CFR 888.3360  
DEVICE CLASS: Class II  
PANEL CODE: KWL – prosthesis, hip, hemi-, femoral, metal  
Orthopedics Panel/87

**A. INTENDED USE:**

The Smith & Nephew Modular Femoral Heads are indicated for the following:

- Noninflammatory degenerative joint disease including osteoarthritis, osteonecrosis, avascular necrosis and post traumatic arthritis;
- rheumatoid arthritis;
- arthritis secondary to a variety of diseases and anomalies and correction of functional deformity such as congenital hip dysplasia or ankylosing spondylitis;
- revision procedures where other treatments have failed; and
- treatment of proximal femoral non-union, femoral neck fracture, and trochanteric fractures of the proximal femur with head involvement.

The modular femoral heads are for single use only and are intended to be used as part of a hemi-hip replacement system when articulating against the natural acetabulum.

**B. DEVICE DESCRIPTION:**

New cobalt chrome (CoCr) modular femoral heads have been designed and developed by Smith & Nephew Orthopaedics. The subject devices are offered in sizes ranging from 38-58mm and feature a female Morse-type taper that has been modified to accept taper sleeves that create a variety of neck length offsets. The taper sleeves feature a 12/14 female taper for attachment to the trunnion of a commercially-available Smith & Nephew femoral stem. The overall design of the component is based upon the existing Smith & Nephew Modular Heads and Uni-polar implants subject of K061243 and K896580, respectively.

**C. SUBSTANTIAL EQUIVALENCE INFORMATION:**

The Smith & Nephew Modular Femoral (Hemi) Heads are similar to the following commercially available devices regarding design features, overall indications, and materials:

- Smith & Nephew Modular Femoral Heads (K061243)
- Smith & Nephew Uni-Polar Head (K896580)



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

SEP 12 2006

Smith & Nephew Orthopaedics  
c/o Mr. Gino J. Rouss  
Project Manager, Regulatory Affairs  
1450 East Brooks Road  
Memphis, Tennessee 38116

Re: K062408

Trade/Device Name: Smith & Nephew Modular Femoral (Hemi) Heads  
Regulation Number: 21 CFR 888.3360  
Regulation Name: Hip joint femoral (hemi-hip) metallic cemented or uncemented prosthesis  
Regulatory Class: Class II  
Product Code: KWL  
Dated: August 15, 2006  
Received: August 17, 2006

Dear Mr. Rouss:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Gino J. Rouss

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or 240-276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a small "for" written below the name.

Mark N. Melkerson

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K062408

Device Name: Smith & Nephew Modular Femoral Heads

### Indications for Use:

The Smith & Nephew Modular Femoral Heads are indicated for the following:

- Noninflammatory degenerative joint disease including osteoarthritis, osteonecrosis, avascular necrosis and post traumatic arthritis;
- rheumatoid arthritis;
- arthritis secondary to a variety of diseases and anomalies and correction of functional deformity such as congenital hip dysplasia or ankylosing spondylitis;
- revision procedures where other treatments have failed; and
- treatment of proximal femoral non-union, femoral neck fracture, and trochanteric fractures of the proximal femur with head involvement.

The modular femoral heads are for single use only and are intended to be used as part of a hemi-hip replacement system when articulating against the natural acetabulum.


Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use             
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -- CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K062408

Page 1 of