

APPENDIX A. 510(k) SUMMARY

OCT 17 2006

Sponsor/Submitter: Abbott Laboratories
Abbott Vascular Inc.
400 Saginaw Drive
Redwood City, CA 94063

Contact Person: Kathryn Marchel
Regulatory Affairs Coordinator
Phone:(650) 474-6341
Fax:(650) 474-3041

Date of Submission: September 21, 2006

Device Trade Name: Fox SV PTA Catheter

Device Classification: Class II

Regulation Number: 21 CFR 870.1250

Classification Name: Percutaneous Transluminal Angioplasty Catheter

Product Code: LIT

Predicate Device: Fox SV PTA Catheter (K051519)

Intended Use: The Fox SV PTA Catheter is intended for dilatation of lesions in the femoral, renal, iliac, popliteal, peroneal, and profunda arteries and native or synthetic arteriovenous dialysis fistulae. This catheter is not intended for the delivery or expansion of stents.

Device Description: The Fox SV PTA Catheter is a standard over-the-wire PTA catheter. The double lumen catheter has a balloon located near the distal tip. One lumen is used for inflation of the balloon, while the second lumen allows access to the distal tip of the catheter for guidewire insertion (max 0.018"). The balloon material expands to a known diameter at specific pressure.

Summary of Substantial Equivalence: The Fox SV PTA Catheter is substantially equivalent to the predicate device. Substantial equivalence was confirmed through non-clinical testing.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

OCT 17 2006

Kathryn Marchel
Regulatory Affairs Coordinator
500 Saginaw Drive
Redwood City, CA 9406
Ph: 650-474-6341
Fax: 650-474-3041

Re: K062843
Trade/Device Name: Fox SV PTA Transluminal Angioplasty Catheter
Regulation Number: 21 CFR 870.1250
Regulation Name: Angioplasty Catheter
Regulatory Class: II (two)
Product Code: LIT
Dated: September 21, 2006
Received: September 22, 2006

Dear Ms. Marchel;

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4586. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

APPENDIX B. INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K062843

Device Name: Fox SV PTA Catheter

Indications For Use: The Fox SV PTA Catheter is intended for dilatation of lesions in the femoral, renal, iliac, popliteal, peroneal, and profunda arteries and native or synthetic arteriovenous dialysis fistulae.

This catheter is not intended for the expansion or delivery of stents.

Prescription Use ☒
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐
(Part 21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Dana R. Vachner
(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K062843