

K063006

Special 510(k) Premarket Notification
GE LOGIQ Works Ultrasound Workstation - Medical Image Management Device
September 28, 2006

Attachment B:
Summary of Safety and Effectiveness
Prepared in accordance with 21 CFR Part 807.92(c).

NOV - 2 2006



GE Healthcare

General Electric Company
P.O. Box 414, Milwaukee, WI 53201

Section a):

1. Submitter: GE Medical Systems, Ultrasound and Primary Care Diagnostics
PO Box 414
Milwaukee, WI 53201

Contact Person: Allen Schuh,
Manager, Safety and Regulatory Engineering
Telephone: 414-721-3992; Fax: 414-721-3899

Date Prepared: September 28, 2006
2. Device Name: GE LOGIQworks™ Ultrasound Image Workstation
System, Image Processing, Radiological, 21 CFR 892.2050, 90-LLZ
3. Marketed Device: RadWorks Review K963699, Platinum PACS Workstation K981217,
Centricity PACS Plus K023557 and Centricity RA 600 K042525; all are 90-LLZ
4. Device Description: The GE LOGIQ Works Image Analysis and Review Workstation is optimized for ultrasound images that are acquired primarily via the GE LOGIQ family of diagnostic ultrasound systems. It is sold with or without the designated computer hardware.
5. Indications for Use: To retrieve and display for reading, diagnostic review, analysis and measurement of digital medical images acquired from DICOM compliant acquisition systems such as CT, MR, PET, X-ray and Ultrasound. Also, for images acquired from GE LOGIQ ultrasound, to process, analyze and manipulate images with the same post acquisition capability as that provided on the GE LOGIQ ultrasound scanner.
6. Comparison with Predicate Device: The GE LOGIQworks™ is of a comparable type and substantially equivalent to the above listed marketed devices with the addition of post acquisition image manipulation and processing capability of the GE LOGIQ 9 Ultrasound. It has equivalent technological characteristics, key safety and effectiveness features, physical design, construction, and materials, and has essentially the same intended uses as the predicate device.

Section b):

1. Non-clinical Tests: The device has been evaluated for conformance to its design specifications and applicable industry standards for software development. It is further verified for system compatibility with the devices with which it communicates. Computer hardware is certified to applicable safety standards.
2. Clinical Tests: None required to confirm safety and effectiveness. However, evaluation in a clinical setting is performed to help assure reliability and compatibility within the intended network environment.
3. Conclusion: Intended uses and other key features of the device are consistent with traditional clinical practice, FDA guidelines and established methods of handling patient examination images and data. The design and development process of the manufacturer conforms with 21 CFR 820, ISO 9001:2000 and ISO 13485:2000 quality systems. The device conforms to applicable medical device safety standards and compliance is verified through internal and independent quality system audit. PACS devices and medical information management systems in general have accumulated a long history of safe and effective performance. Therefore, it is the opinion of GE Healthcare that the GE LOGIQworks™ is substantially equivalent with respect to safety and effectiveness to devices currently cleared for market.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration
9200 Corporate Blvd.
Rockville MD 20850

Mr. Allen Schuh
Manager, GE Ultrasound Safety
and Regulatory Engineering
General Electric Company
GE Medical Systems, Ultrasound and
Primary Care Diagnostics, LLC
9900 Innovation Drive
WAUWATOSA WI 53226

NOV - 2 2006

Re: K063006
Trade/Device Name: GE LOGIQworks™
Regulation Number: 21 CFR §892.2050
Regulation Name: Picture archiving and communications systems
Regulatory Class: II
Product Code: LLZ
Dated: September 28, 2006
Received: October 21, 2006

Dear Mr. Schuh:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.



Protecting and Promoting Public Health

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

| | | |
|----------------|----------------------------------|--------------|
| 21 CFR 876.xxx | (Gastroenterology/Renal/Urology) | 240-276-0115 |
| 21 CFR 884.xxx | (Obstetrics/Gynecology) | 240-276-0115 |
| 21 CFR 894.xxx | (Radiology) | 240-276-0120 |
| Other | | 240-276-0100 |

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Attachment E

Indications for Use

510(k) Number (if known): K063006

Device Name: GE LOGIQworks™

Indications For Use: To retrieve and display for reading, diagnostic review, analysis and measurement of digital medical images acquired from DICOM compliant acquisition systems such as CT, MR, PET, X-ray and Ultrasound.

Also, for images acquired from GE LOGIQ ultrasound, to process, analyze and manipulate images with the same post acquisition capability as that provided on the GE LOGIQ ultrasound scanner.

Prescription Use XXX
(Part 21 CFR 801 Subpart D)

AND/OR

~~Over-The-Counter Use~~ _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

David R. Beggs
(Division Sign-Off)

Division of Reproductive, Abdominal, and
Radiological Devices

510(k) Number K063006