

K063130

510(k) Summary of Safety and Effectiveness (21 CFR 807.92)

Submitter's Name: Toshiba America Medical Systems, Inc.
Address: P.O. Box 2068,
2441 Michelle Drive
Tustin, CA 92781-2068
Contact: Paul Biggins, Senior Regulatory Affairs Manager
Telephone No.: (714) 730-5000

NOV - 2 2006

Device Proprietary Name: SSA-790A, Aplio™ XG Version 1.00
Common Name: Ultrasound Imaging System

Classification:
Regulatory Class: II
Review Category: Tier II

- Ultrasonic Pulsed Doppler Imaging System – Product Code: 90-IYN
(21 CFR 892.1550)
- Ultrasonic Pulsed Echo Imaging System – Product Code: 90-IYO
(21 CFR 892.1560)
- Diagnostic Ultrasonic Transducer – Product Code: 90-ITX
(21 CFR No. 892 1570)

Identification of Predicate Devices:

Toshiba America Medical System believes that this device is substantially equivalent to the following device:

Toshiba SSA-770 Aplio Version 5.5 Diagnostic Ultrasound, K032281

Device Description:

The Aplio XG Ultrasound System is a mobile system. It is a Track 3 device that employs a wide array of probes that include flat linear array, convex linear array, and sector array with a frequency range of 2 MHz to 12 MHz

Intended Use:

The Aplio XG System is intended to be used for the following types of studies: feta, abdominal, intraoperative, pediatric, small organs, neonatal cephalic, adult cephalic, cardiac, transrectal, transvaginal, transesophageal, peripheral vascular and musculoskeletal (both conventional and superficial).

Safety Considerations:

This device is designed and manufactured in conjunction with the Quality System Regulations IEC-60601 (applicable portions), and IEC 60601-1-2 (applicable portions, IEC60601-2-37 (applicable portions), and the AIUM-NEMA UD2 Output Measurement Standard as applied to Track Ultrasound Systems and the AIUM-NEMA UD3 Output Display Standard.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Toshiba Medical Systems Corporation
c/o Mr. Mark Job
Responsible Third Party Official
Regulatory Technology Services LLC
1394 25th Street, N.W.
BUFFALO MN 55313

NOV - 2 2006

Re: K063130

Trade Name: SSA-790A, Aplio™ XG Version 1.00
Regulation Number: 21 CFR §892.1550
Regulation Name: Ultrasonic pulsed doppler imaging system
Regulatory Class: II
Product Code: IYN, IYO and ITX
Dated: October 10, 2006
Received: October 13, 2006

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the SSA-790A, Aplio™ XG Version 1.00 System, as described in your premarket notification:

Transducer Model Numbers

PVT-375BT PVT-661VT PLT-1202S PC-20M PET-510MB PST-25BT PLT-604AT

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

If you have any questions regarding the content of this letter, please contact Ralph Schuping, Sc.D., at (240) 276- 3666.

Sincerely yours,

A handwritten signature in cursive script that reads "Nancy C. Brogdon".

Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosures

Diagnostic Ultrasound Indications For Use Form

System X Transducer _____

Model SSA-790A

510(k) Number(s) _____

K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal	N	N	N	N	N	N		N				N
Abdominal	N	N	N	N	N	N		N	N			N
Intraoperative (Specify)	N	N	N	N	N			N				N
Intraoperative Neurological												
Pediatric	N	N	N	N	N	N	N	N	N			N
Small Organ (Specify)*	N	N	N	N	N			N				N
Neonatal Cephalic	N	N	N	N	N	N		N	N			N
Adult Cephalic	N	N	N	N	N	N		N	N			N
Cardiac	N	N	N	N	N	N	N	N	N	N		N
Transesophageal	N	N	N	N			N	N	N			N
Transrectal	N	N	N	N	N	N		N				N
Transvaginal	N	N	N	N	N	N		N				N
Transurethral												
Intravascular												
Peripheral Vascular	N	N	N	N	N			N				N
Laparoscopic												
Musculo-skeletal Superficial	N	N	N	N	N			N				N
Musculo-skeletal Conventional	N	N	N	N	N			N				N

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

* : For example: thyroid, parathyroid, breast, scrotum and penis

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Ancyc Brogdon

(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices
510(k) Number K063130

Diagnostic Ultrasound Indications For Use Form

System _____ Transducer X

Model PVT-375BT

510(k) Number(s) K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal	P	P	P	P	P	P		P				P
Abdominal	P	P	P	P	P	P		P				P
Intraoperative (Specify)												
Intraoperative Neurological												
Pediatric	P	P	P	P	P	P		P				P
Small Organ (Specify)*												
Neonatal Cephalic												
Adult Cephalic												
Cardiac												
Transesophageal												
Transrectal												
Transvaginal												
Transurethral												
Intravascular												
Peripheral Vascular												
Laparoscopic												
Musculo-skeletal Superficial												
Musculo-skeletal Conventional												

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

510(k) for this transducer k041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
Concurrence of CDRII, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Nancy C. Brogdon
(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices
510(k) Number K063130

Diagnostic Ultrasound Indications For Use Form

System ____ Transducer X

Model PVT-661VT

510(k) Number(s)

K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal												
Abdominal												
Intraoperative (Specify)												
Intraoperative Neurological												
Pediatric												
Small Organ (Specify)*												
Neonatal Cephalic												
Adult Cephalic												
Cardiac												
Transesophageal												
Transrectal	P	P	P	P	P	P		P				P
Transvaginal	P	P	P	P	P	P		P				P
Transurethral												
Intravascular												
Peripheral Vascular												
Laparoscopic												
Musculo-skeletal Superficial												
Musculo-skeletal Conventional												

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

510(k) for this transducer k041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Nancy C Brogdon
(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices
510(k) Number K063130

Diagnostic Ultrasound Indications For Use Form

System _____ Transducer X

Model PLT-1202S

510(k) Number(s) _____

K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal												
Abdominal												
Intraoperative (Specify)	E	E	E	E	E			E				E
Intraoperative Neurological												
Pediatric												
Small Organ (Specify)*	E	E	E	E	E			E				E
Neonatal Cephalic												
Adult Cephalic												
Cardiac												
Transesophageal												
Transrectal												
Transvaginal												
Transurethral												
Intravascular												
Peripheral Vascular	E	E	E	E	E			E				E
Laparoscopic												
Musculo-skeletal Superficial	E	E	E	E	E			E				E
Musculo-skeletal Conventional	E	E	E	E	E			E				E

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

510(k) for this transducer k041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Nancy C. Brogdon
(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number _____

K063130

Diagnostic Ultrasound Indications For Use Form

System _____ Transducer X

Model PC-20M

510(k) Number(s) _____

K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal												
Abdominal												
Intraoperative (Specify)												
Intraoperative Neurological												
Pediatric									P			
Small Organ (Specify)*												
Neonatal Cephalic												
Adult Cephalic												
Cardiac									P			
Transesophageal												
Transrectal												
Transvaginal												
Transurethral												
Intravascular												
Peripheral Vascular									P			
Laparoscopic												
Musculo-skeletal Superficial												
Musculo-skeletal Conventional												

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

510(k) for this transducer k041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number _____

K063130

Diagnostic Ultrasound Indications For Use Form

System Transducer X

Model PET-510MB

510(k) Number(s) K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal												
Abdominal												
Intraoperative (Specify)												
Intraoperative Neurological												
Pediatric												
Small Organ (Specify)*												
Neonatal Cephalic												
Adult Cephalic												
Cardiac												
Transesophageal	P	P	P	P			P	P	P			P
Transrectal												
Transvaginal												
Transurethral												
Intravascular												
Peripheral Vascular												
Laparoscopic												
Musculo-skeletal Superficial												
Musculo-skeletal Conventional												

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

510(k) for this transducer k041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
 Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Nancy Brozdon
 (Division Sign-Off)
 Division of Reproductive, Abdominal,
 and Urological Devices
 510(k) Number K063130

Diagnostic Ultrasound Indications For Use Form

System _____ Transducer X

Model PST-25BT

510(k) Number(s) _____

K06 3130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal												
Abdominal	E	E	E	E	E	E	E	E	E			E
Intraoperative (Specify)												
Intraoperative Neurological												
Pediatric	E	E	E	E	E	E	E	E	E			E
Small Organ (Specify)*												
Neonatal Cephalic	E	E	E	E	E	E	E	E	E			E
Adult Cephalic	E	E	E	E	E	E	E	E	E			E
Cardiac	E	E	E	E	E	E	E	E	E	E		E
Transesophageal												
Transrectal												
Transvaginal												
Transurethral												
Intravascular												
Peripheral Vascular												
Laparoscopic												
Musculo-skeletal												
Superficial												
Musculo-skeletal Conventional												

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;

BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;


CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

Letter to file against K041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


 (Division Sign-Off)
 Division of Reproductive, Abdominal,
 & Neurological Devices
 510(k) Number K063130

Diagnostic Ultrasound Indications For Use Form

System _____ Transducer X

Model PLT-604AT

510(k) Number(s) _____

K063130

Clinical Application	Mode of Operation											
	B	THI	M	Color Doppler	Power	Dynamic Flow	TDI	PW	CW	CHI 2D	CHI Dynamic Flow	Combined (Specify)
Ophthalmic												
Fetal												
Abdominal												
Intraoperative (Specify)												
Intraoperative Neurological												
Pediatric												
Small Organ (Specify)*	P	P	P	P	P	P		P				P
Neonatal Cephalic												
Adult Cephalic												
Cardiac												
Transesophageal												
Transrectal												
Transvaginal												
Transurethral												
Intravascular												
Peripheral Vascular	P	P	P	P	P	P		P				P
Laparoscopic												
Musculo-skeletal Superficial	P	P	P	P	P	P		P				P
Musculo-skeletal Conventional	P	P	P	P	P	P		P				P

N= new indication; P = Previously Cleared by FDA; E = Added under Appendix E (LTF)

Additional Comments: _____ Combined Modes: B/M; B/PWD;
BDF/PWD; BDF/MDF; BDF/MDF/PWD; B-TDI; M-TDI; 2D/CWD; BDF/CWD;
CHI/2D; FEI/2D; CHI/BDF; FEI/BDF

510(k) for this transducer k041499

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON OTHER PAGES IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Nancy C Brogdon
(Division Sign-Off)
Division of Reproductive, Abdominal,
Neurological Devices
File Number K063130