NOV 1 5 2006

510(k) Summary As Required by 21 section 807.92 (c)

1-Submitter Name:

SHANGHAI DENTAL INSTRUMENT FACTORY CO., LTD

("SDIF")

2-Address:

820 LINGSHI ROAD, SHANGHAI CHINA 200072

3-Phone:

4-Fax:

Tel +86 21 66255988. Fax +86 21 56956197

5-Contact Person:

Ms Chen Yimei, Vice General Manager

6- Consultant: Jay Mansour, Mansour Consulting LLC, 845 Aronson Lake Court, Roswell, GA 30075 USA. Tel 678-908-8180. Fax 678-623-3765

7-Date summary prepared: October 26th, 2006

8-Device Trade or Proprietary Name: SDIF Dental chair with Operative Unit and

accessories

9-Device Common or usual name: Dental chair with operative unit 10-Device Classification Name: Chair, Dental, with Operative unit

11-Substantial Equivalency is claimed against the following device: K052470

12-Description of the Device:

SDIF dental chairs with operative units provide patient comfort and dentists' air and water supplies for dental instruments and procedures. Water pressure: 200~400Kpa. Air pressure: 600~800Kpa. The units can be powered from ~110/220/230V. The chair is motorized. The system is compatible with instruments from various manufacturers with standard fittings.

13-Intended use of the device:

SDIF dental chairs with Operative Unit is intended for use to supply power to and serve as a base for dental devices and accessories. This product includes a dental chair. The unit is intended for use in the dental clinic/office environment and used by trained dentists and/or dental technicians and assistants.

14-Safety and Effectiveness of the device:

This device is safe and effective as the other predicate device cited above.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Shanghai Fimet Medical Instrument Company, Limited C/O Mr. Mark Job Responsible Third Party Official Regulatory Technology Services, LLC 1394 25th Street, NW Buffalo, Minnesota 55313

NOV 15 2006

Re: K063302

Trade/Device Name: SDIF Dental Chairs with Operative Unit

Regulation Number: 872.6640

Regulation Name: Dental Operative Unit and Accessories

Regulatory Class: I Product Code: EIA

Dated: October 30, 2006 Received: November 1, 2006

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices Office of Device Evaluation

quite y. Michian Om.s.

Center for Devices and Radiological Health

Enclosure

Indications for Use

Device Name: SDIF DENTAL CHAIRS WITH OPERATIVE UNIT

510(k) Number (if known):

Indications For Use:
SDIF Dental chairs with operative unit is indicated for use to supply power to and serve as a base for dental devices and accessories. This product includes a dental chair. The unit is indicated for use in the dental clinic/office environment used by trained dentists and/or dentatechnicians and assistants
Prescription Use AND/OR Over-The-Counter Use (21 CFR 801 Subpart C)
(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CDRH, Office of Device Evaluation (ODE)
Sysn Rumer
Cign-Off) Ten of Anesthesiology, General Hospital, Control, Dental Devices Page 1 of 1
1) THE OF MANAGEMENT AND ADDRESS OF THE PARTY OF THE PART