



K063326

Page 1 of 3

11311 Concept Boulevard • Largo, FL 33773-4908 • 727-392-6464 • www.linvatec.com

January 4, 2007

JAN 31 2007

**510(k) SUMMARY**  
***ConMed Linvatec 24k Irrigation System***

In accordance with the requirements of the Safe Medical Device Act of 1990 and 21 CFR 807.92, ConMed Linvatec is hereby submitting the 510(k) Summary of Safety and Effectiveness for the 510(k) Number \_\_\_\_\_.

**A. Submitter**

ConMed Linvatec  
11311 Concept Boulevard  
Largo, Florida 33773-4908  
Registration Number: 1017294

**B. Company Contact**

Sue F. Dauterman  
Regulatory Affairs Specialist  
(727) 399-5321 Telephone  
(727) 399-5264 FAX

**C. Device Name**

Trade Name: *ConMed Linvatec 24k Irrigation System*

Common Name: Irrigation System

Classification Name: Arthroscope, 888.1100

Proposed Class/Device: Class II

Product Code: HRX

K063326

Page 2 of 3



11311 Concept Boulevard • Largo, FL 33773-4908 • 727-392-6464 • www.linvatec.com

510(k) Summary  
24k Irrigation System  
January 4, 2007

**D. Predicate/Legally Marketed Devices**

- 10k Irrigation System 510(k) # K033573  
Linvatec Corporation
- 3M Fluid Irrigation System 510(k) #K941388  
Linvatec Corporation (acquired from 3M Healthcare)
- Nextra Arthroscopic Pump and Shaver System 510(k) #K954465  
Future Medical Systems, Inc.
- Arthroscopy Pump, Model A107 510(k) #K030402  
W.O.M. World of Medicine, Ag  
Marketed by Stryker Endoscopy as FloControl
- Hall Irrigation System 510(k) # K852143  
Linvatec Corporation (acquired from Hall Surgical)
- HydroFlex LI Laparoscopic Irrigation System 510(k) # K961224  
Davol, Inc.

**E. Device Description**

The *ConMed Linvatec 24k Irrigation System* is a peristaltic pump system designed to automatically provide and control distention and irrigation of the operative site during arthroscopic procedures and fluid irrigation of the operative site during laparoscopic procedures using sterile fluids. The pump is used in conjunction with specific tubing sets designed for either arthroscopic or laparoscopic procedures.

**F. Intended Use**

The *ConMed Linvatec 24k Irrigation System* is a dual-head pump system intended to provide fluid distention, irrigation, and suction during diagnostic and operative knee, shoulder, hip, and small bone arthroscopic procedures and fluid irrigation during laparoscopic procedures.



1063326

Page 3 of 3

11311 Concept Boulevard • Largo, FL 33773-4908 • 727-392-6464 • www.linvatec.com

510(k) Summary  
24k Irrigation System  
January 4, 2007

### **G. Substantial Equivalence**

The *ConMed Linvatec 24k Irrigation System* is substantially equivalent in intended use, design and technological characteristics to the 10k Irrigation System (Linvatec Corporation), 3M Fluid Irrigation System (Linvatec Corporation), Nextra Arthroscopic Pump and Shaver System (Future Medical Systems, Inc.), Arthroscopy Pump, Model A107 (W.O.M. World of Medicine, Ag), HydroFlex LI Laparoscopic Irrigation System (Davol, Inc.) and Hall Irrigation System (Linvatec Corporation). Testing conducted prior to product release will assure the new device does not raise any new issues of safety and efficacy.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

ConMed Linvatec  
% Ms. Sue F. Dauterman  
Regulatory Affairs Specialist  
11311 Concept Boulevard  
Largo, Florida 33773-4908

JAN 8 1 2007

Re: K063326

Trade/Device Name: *ConMed Linvatec 24k Irrigation System*  
Regulation Number: 21 CFR 888.1100  
Regulation Name: Arthroscope  
Regulatory Class: II  
Product Code: HRX  
Dated: January 4, 2007  
Received: January 9, 2007

Dear Ms. Dauterman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Sue F. Dauterman

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a large, sweeping flourish extending to the right.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

16063326



11311 Concept Boulevard • Largo, FL 33773-4908 • 727-392-6464 • www.linvatec.com

October 16, 2006

510(k) Number (if known):

Device Name: *ConMed Linvatec 24k Irrigation System*

Indications for Use:

The *ConMed Linvatec 24k Irrigation System* is a dual-head pump system intended to provide fluid distention, irrigation, and suction during diagnostic and operative knee, shoulder, hip, and small bone arthroscopic procedures and fluid irrigation during laparoscopic procedures.

Prescription Use  AND/OR  
(Part 21 CFR 801 Subpart D)

Over-the-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

**(Division Sign-Off)**  
**Division of General, Restorative,**  
**and Neurological Devices**

510(k) Number 16063326