

## Summary of Safety and Effectiveness

K063414

DEC 12 2006

Date: November 9, 2006

Manufacturer:

Encore Medical, L.P.  
9800 Metric Blvd  
Austin, TX 78758

Trade Name: FMP Metal/Metal Acetabular  
Insert

Common Name: Prosthesis, hip, semi-  
constrained

Contact Person:

Teffany Hutto  
Regulatory Affairs Specialist  
Phone: (512) 834-6255  
Fax: (512) 834-6313  
Email: Teffany\_Hutto@encoremed.com

Classification Name: Hip joint metal/metal  
semi-constrained, with an uncemented  
acetabular component, prosthesis, 21CFR  
888.3330

Description: The modification to the system consists of a change in the method or porous coating of the acetabular shells from a two layer process to a three layer process utilizing a smaller bead size and smaller pore size.

Intended Use: The FMP Metal/Metal Acetabular Insert used in total hip is intended for conditions of degenerative joint disease including osteoarthritis and avascular necrosis; rheumatoid arthritis; correction of functional deformity; revision procedures where other treatments or devices have failed; and treatment of nonunion, femoral neck and trochanteric fractures of the proximal femur with head involvement, which is unmanageable using other techniques.

Comparable Features to Predicate Device(s): Features comparable to predicate devices include the same materials, design, indications, packaging, labeling, and sterilization,



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 12 2006

Encore Medical, L.P.  
% Ms. Teffany Hutto  
Regulatory Affairs Specialist  
9800 Metric Boulevard  
Austin, Texas 78758

Re: K063414  
Trade/Device Name: FMP Metal/Metal Acetabular Insert  
Regulation Number: 21 CFR 888.3330  
Regulation Name: Hip joint metal/metal semi-constrained, with an uncemented acetabular component, prosthesis  
Regulatory Class: Class III  
Product Code: KWA  
Dated: November 9, 2006  
Received: November 21, 2006

Dear Ms. Hutto:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

Page 2 – Ms. Teffany Hutto

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", written over a horizontal line.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): \_\_\_\_\_

Device Name: FMP Metal/Metal Acetabular Insert

Indications for Use:

**FMP Metal/Metal Acetabular Insert  
Indications for Use**

Indications for use in total hip replacement include: degenerative joint disease including osteoarthritis and avascular necrosis; rheumatoid arthritis; correction of functional deformity; revision procedures where other treatments or devices have failed; and treatment of nonunion, femoral neck and trochanteric fractures of the proximal femur with head involvement, which is unmanageable using other techniques.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Barbara Bruchman Foxman  
(Division Sign-Off)  
Division of General, Restorative,  
and Neurological Devices

510(k) Number K063414