


## SECTION 4

K063453

AUG 30 2007

### SECTION 4 – 510(k) SUMMARY

[As required by 21CFR807.92]



ELITE SURGICAL SUPPLIES (PTY) LTD

P.O. Box 26115  
Arcadia  
Pretoria  
0007  
South Africa

Contact Person: John O'Toole  
Contact Title: Managing Director  
Tel: +27 12 386 0012  
Fax: +27 12 386 2745  
Email: john@elitesurgical.com

## SECTION 4

### 4.1 Date Prepared [21CFR807.92(a)(1)]

October 4, 2006

### 4.2 Submitter's Information [21CFR807.92(a)(1)]

<b>Company Name:</b>	Elite Surgical US
<b>Street Address:</b>	503 Commerce Park Drive, Suite I
<b>City:</b>	Marietta
<b>State/Province:</b>	Georgia
<b>Country:</b>	USA
<b>Telephone:</b>	404-408-2396
<b>Facsimile:</b>	770-590-5152
<b>Contact Person:</b>	Charise de Barros
<b>Contact Title:</b>	Agent
<b>Contact Email:</b>	Charise@EliteSurgicalUSA.com

### 4.3 Trade Name, Common Name, Classification [21CFR807.92(a)(2)]

<b>Trade Name:</b>	Vertefix® Pedicle Screw Spinal System
<b>Common Name:</b>	Pedicle Screw Spinal Fixation System
<b>Classification Name:</b>	Spinal Pedicle Screw Fixation Orthosis(MNI) per 21 CFR § 888.3070
	Spondylolisthesis Spinal Fixation Orthosis (MNH) per 21 CFR § 888.3070
	Spinal Interlaminar Fixation Orthosis (KWP) Per 21 CFR § 888.3050
<b>Device Class:</b>	Class II
<b>Product Code:</b>	MNI, MNH, KWP

## SECTION 4

### 4.4 Identification of Predicate Device(s) [21CFR807.92(a)(3)]

PREDICATE DEVICES
U and I Corporation, OPTIMA™ Spinal System MNI, MNH, KWQ (K051971)
K2M Inc, Denali™ Spine System KWP, MNH, MNI (K042635)
Howmedica Osteonics Corp, XIA™ Spine System Hooks MNH (K001319)

There are no significant differences between the Vertefix® Pedicle Screw Spinal System and the OPTIMA™ Spinal System, the Denali™ Spine System, the XIA™ Spine System Hooks, or other spinal fixation systems currently being marketed, which would adversely affect the use of the product. It is substantially equivalent to these other devices in design, function, materials, operational principles and intended use.

### 4.5 Description of the Device [21CFR807.92(a)(4)]

The Vertefix® Pedicle Screw Spinal System is a multiple component, posterior spinal fixation system which consists of pedicle screws, rods, set screws, connectors, hooks, and a transverse (cross) linking mechanism.

The Vertefix® Pedicle Screw Spinal System will allow surgeons to build a spinal implant construct to stabilize and promote spinal fusion. The Vertefix® Pedicle Screw Spinal System components are supplied non-sterile, and are single use.

### 4.6 Intended Use [21CFR807.92(a)(5)]

The Vertefix® Pedicle Screw Spinal System is intended for posterior, non cervical pedicle fixation in skeletally mature patients receiving fusion by autogenous bone graft in the thoracolumbar spine for the following indications: Spondylolisthesis, Trauma (e.g. fracture or dislocation), Spinal stenosis, deformities or curvatures (scoliosis, kyphosis or lordosis), tumor, pseudoarthrosis, and failed previous fusion.

Removal of the implants after the attainment of a solid fusion is optional.

## **SECTION 4**

### **4.7 Technological Characteristics [21CFR807.92(a)(6)]**

The Vertefix® Pedicle Screw Spinal System components were biomechanically tested in accordance with ASTM F1717. This testing demonstrated substantial equivalence to the predicate systems and other marketed systems. Further comparison to the predicates also demonstrated substantial equivalence in terms of intended use, operating principle, materials, shelf life and design. There are no significant differences between the Vertefix® spinal system, and other systems currently marketed which would adversely affect the use of the product. The Vertefix® Spinal System components are fabricated from titanium alloy (Ti-6Al-4v ELI) that conforms to ASTM F 136. Various sizes of these implants are available for the application.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Elite Surgical Supplies (Pty) Ltd.  
% Ms. Charise De Barros  
Agent  
503 Commerce Park Drive SE, Suite i  
Marietta, GA 30060

AUG 30 2007

Re: K063453  
Trade/Device Name: VERTEFIX® Pedicle Screw Spinal System  
Regulation Number: 21 CFR 888.3070  
Regulation Name: Pedicle screw spinal system  
Regulatory Class: II  
Product Code: MNI, MNH, KWP  
Dated: August 6, 2007  
Received: August 9, 2007

Dear Ms. De Barros:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson". The signature is fluid and cursive, with the first name "Mark" being the most prominent.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## SECTION 3

### SECTION 3 – INDICATIONS FOR USE

Applicant Elite Surgical Supplies (Pty) Ltd  
510(k) Number (if known) K063453  
Device Name Vertefix® Pedicle Screw Spinal System

**Indications for Use:** The Vertefix® Pedicle Screw System is intended for posterior, non-cervical pedicle fixation in skeletally mature patients receiving fusion by autogenous bone graft in the thoracolumbar spine, for the following indications:

- Spondylolisthesis;
- Trauma (e.g. fracture or dislocation);
- Spinal stenosis;
- Deformities or curvatures (scoliosis, kyphosis and lordosis);
- Tumor;
- Pseudoarthrosis; and
- Failed previous fusion.

Prescription Use ☒ OR Over-the-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109)  
(Optional Format 1-2-96)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER  
PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of General, Restorative,  
and Neurological Devices

510(k) Number K 063453  
Elite Surgical Supplies  
Vertefix® Pedicle Screw Spinal System