SHANGHAI MOTEX HEATHCARE CO., LTD.

No.369, Jiasong Zhong Road, Huaxin, Qingpu, Shanghai, 201708, China Telephone: 86-21-5979 9888 Fax: 86-21-5979 9728 E-mail: motex@public1.sta.net.cn

"___510(k) SUMMARY "

The assigned 510(k) number is : <u>K063757</u>

Submitter's Name: SHANGHAI MOTEX HEALTHCARE CO., LTD.

No.369, Jiasong Zhong Road, Huaxin, Qingpu, Shanghai, 201708, China

Date summary prepared:

November 26, 2006

Device Name:

Classification name:

Surgeon's Gloves

• Classification number:

KGO, Class I

• Regulation Number:

878.4460

• Proprietary name:

Motex Powder-Free Surgical Gloves & Powdered Latex

Surgical Gloves

• Predicate Device:

K050071. MEDISPO powdered surgeon's gloves and

MEDISPO-PF powder-free surgeon's gloves

• Official Correspondent:

Dr. Jen, Ke-Min

E-mail: ceirs.jen@msa.hint.net (Tel) 886-3-5208829; (Fax) 886-3-5209783

Address: No.58, Fu Chiun Street, Hsin Chu City, 30067, Taiwan, ROC

Intended Use:

Motex Powder-Free Surgical Gloves are sterile disposable devices made of natural rubber latex (that may bear a trace amount of glove powder) and is intended to be worn on the hands, usually in surgical settings, to provide a barrier against potentially infectious materials and other contaminants.

SHANGHAI MOTEX HEATHCARE CO., LTD.

No.369, Jiasong Zhong Road, Huaxin, Qingpu, Shanghai, 201708, China Telephone: 86-21-5979 9888 Fax: 86-21-5979 9728 E-mail: motex@public1.sta.net.cn

Motex Powdered Latex Surgical Gloves are sterile disposable devices made of natural rubber latex that bears powder to facilitate donning, and it is intended to be worn on the hands, usually in surgical settings, to provide a barrier against potentially infectious materials and other contaminants..

• Technological Characteristics:

The Motex Powder-Free Surgical Gloves & Powdered Latex Surgical Gloves characteristics are summarized below compared to ASTM and ISO standards to the predicate device:

<u>Characteristic</u> <u>Standard</u>

Dimensions meets ASTM D 3577-06,
Physical Properties meets ASTM D 3577-06,
Freedom from Holes meets ASTM D 3577-06,

Biocompatibility meets ISO10993-5/-10

Sterilization Validation meets ISO11137

Clinical Data:

Not applicable.

Conclusions:

The Motex Powder-Free Surgical Gloves & Powdered Latex Surgical Gloves and predicate devices also meet the technological characteristics of ASTM D 3577-06. ISO10993-5/-10, and ISO11137 standards. Besides, for Powder-Free SURGICAL GLOVES contain no more than 2mg powder and no more than 50ug/dm2 extractable protein, and for Powdered Latex SURGICAL GLOVES contain no more than 120mg powder and no more than 120ug/dm2 extractable protein claim.

Thus the new devices are substantially equivalent to the predicate devices.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

JUL 2 7 2007

Shanghai Motex Healthcare Company Limited C/O Dr. Jen Ke-Min Official Correspondent ROC Chinese-European Industrial No. 58, Fu Chiun Street Hsin Chu City, Taiwan 30067 CHINA

Re: K063757

Trade/Device Name: Motex Powder-Free Surgical Gloves and Powdered

Latex Surgical Gloves

Regulation Number: 878.4460 Regulation Name: Surgeon's Glove

Regulatory Class: I Product Code: KGO Dated: June 15, 2007 Received: June 29, 2007

Dear Dr. Jen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

Office of Device Evaluation

outey Michael MA

Center for Devices and Radiological Health

Indications for Use

510 (K) Number (If Known	():K0637	(5)
Device Name: <u>Motex Powo</u>	der-Free Surg	ical Gloves and Powdered Late
Surgical Glov	<u>′es</u>	
Indications for Use:		
latex (that may bear a trace am	nount of glove po tings, to provide	lisposable devices made of natural rubbowder) and is intended to be worn on t a barrier against potentially infectio
Motex Powdered Latex Surgice		erile disposable devices made of natur
	tings, to provide	ing, and it is intended to be worn on t a barrier against potentially infectio
hands, usually in surgical sett materials and other contaminant	tings, to provide ts.	a barrier against potentially infectio
hands, usually in surgical sett	tings, to provide	
hands, usually in surgical sett materials and other contaminant. Prescription Use (Part 21 CFR 801 Subpart D)	tings, to provide ts. AND/OR	a barrier against potentially infection \mathbf{v} Over-The-Counter Use \mathbf{v}
hands, usually in surgical sett materials and other contaminant. Prescription Use	tings, to provide ts. AND/OR OW THIS LINE	a barrier against potentially infection Over-The-Counter Use √ (21 CFR 807 Subpart C)