

K 070063

510(K) Summary

Submitter:	Cynosure, Inc. 5 Carlisle Road Westford, MA 01886	JAN 30 2007
Contact:	George Cho Senior Vice President of Medical Technology	
Date Summary Prepared:	January 4, 2007	
Device Trade Name:	Cynosure YAG MIR II Laser	
Common Name:	Medical Laser System	
Classification Name:	Instrument, surgical, powered, laser 79-GEX 21 CFR 878.4810	
Equivalent Device:	Cynosure YAG Family laser, Cynosure YAG MIR Family laser	
Device Description:	The Cynosure YAG MIR II laser is a Nd:YAG laser, having a ND:YAG crystal rod as a lasing medium. Laser activation is by footswitch or finger switch. Overall weight of the laser is 285lbs, and the size is 41"x18"x32" (HxWxD). Electrical requirement is 220 VAC, 30A, 50-60 Hz, single phase.	
Intended Use:	The Cynosure YAG MIR II Laser is indicated for the treatment of vascular and pigmented lesions and wrinkles.	
Comparison:	The Cynosure YAG MIR II Laser has the same indications for use, the same principle of operation, and same wavelengths and pulse energy range as the predicate device(s).	
Nonclinical Performance Data:	none	
Clinical Performance Data:	none	
Conclusion:	The Cynosure YAG MIR II Laser is a safe and effective device for the indications specified.	
Additional Information:	none	



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Cynosure Incorporated
% Mr. George Cho
Senior Vice President
5 Carlisle Road
Westford, Massachusetts 01886

JAN 30 2007

Re: K070063

Trade/Device Name: Cynosure YAG MIR II Laser

Regulation Number: 21 CFR 878.4810,

Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology.

Regulatory Class: Class II

Product Code: GEX

Dated: January 4, 2007

Received: January 8, 2007

Dear Mr. Cho:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. George Cho

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Mark N. Melkerson".

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

5PB

510(k) Number (if known): K 070063

Device Name: Cynosure YAG MIR II Laser

Indications For Use:

1064nm: The Cynosure YAG MIR II laser is intended for the coagulation and hemostasis of benign vascular lesions, such as , but not limited to, port wine stains, hemangiomas, warts, telangiectasia, rosacea, venus lake, leg veins, spider veins and poikiloderma of civatte; and treatment of benign cutaneous lesions such as warts, scars, striae and psoriasis. The laser is also intended for the treatment of benign pigmented lesions such as, but not limited to, lentigos (age spots), solar lentigos (sun spots), Café au lait macules, seborrheic keratoses, nevi, chloasma, verrucea, skin tags, keratoses, tattoos (significant reduction in the intensity of black and/or blue/black tattoos) and plaques.

The laser is also indicated for the treatment of wrinkles such as, but not limited to, periocular and perioral wrinkles.

1320nm: The Cynosure YAG MIR II laser is indicated for use in general surgery and dermatology for the incision, excision, ablation, vaporization, coagulation and hemostasis of soft tissue. It is also indicated for the treatment of periorbital and perioral wrinkles. It is also indicated for the treatment of fine lines and wrinkles.

1440nm: The Cynosure YAG MIR II laser is indicated for use in general surgery and dermatology for the incision, excision, ablation, vaporization, coagulation and hemostasis of soft tissue. It is also indicated for the treatment of periorbital and perioral wrinkles and pigmented lesions.

Prescriptive Use X OR Over-The-Counter Use _____
(Part 21 CFR 801 Subpart D) (Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

**Division of General, Restorative,
and Neurological Devices**

510(k) Number

16040063