

k070181

APR 18 2007

**DSC/ALF Spinal System  
510(k) Summary  
December 2006**

- I. Company:** Sintea Biotech, Inc.  
407 Lincoln Rd. Suite 10L  
Miami Beach, FL 33139  
(305) 673-6226
- II. Proprietary Trade Name:** DSC/ALF Spinal System
- Regulation Number:** 888.3060
- Regulation Name:** Spinal Intervertebral Body Fixation Orthosis
- Product Code:** MQP, KWQ

**III. Product Description**

The Sintea Biotech DSC/ALF system is a dorsolumbar somatic vertebral body replacement device with its own supplemental fixation. The DSC/ALF Spinal System provides two basic components: A VBR (DSC) and an anterior lateral fixation plate (ALF). The DSC has a hollow internal module that can slide in relation to an external tightening module, which enables the former to be locked in the most suitable position by means of two screws for each individual case. A safety screw is also provided. Both modules can be completed with extremity covers that optimise the contact with the vertebral bodies between which the device is inserted, thus ensuring better primary stability. The ALF is composed of two plates that connect to the anterior side of a vertebra by two screws. The two plates are connected by a rod that ranges from 32mm to 92mm in length. All components of the DSC/ALF Spinal System are made of medical grade titanium alloy (Ti-6Al-4V) as described by ASTM Standard F136.

**IV. Indications**

The DSC/ALF Spinal System is a vertebral body replacement system intended for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma (i.e., fracture). The DSC/ALF Spinal System is to be used with supplemental fixation. Specifically, the DSC/ALF system is to be used with the anterior lateral plate that is included in the system, and may also be used with the Sintea Biotech PLS Spinal System. Additionally, the use of bone grafting material with the DSC/ALF System is optional.

K070181

## **V. Performance Data**

Mechanical testing in accordance with the "Guidance for Industry and FDA Staff, Guidance for Spinal Systems 510(k)s", May 2004 is presented in the body of the 510(k).

## **VI. Substantial Equivalence**

Sintea Biotech, Inc. believes that the DSC/ALF Spinal System is substantially equivalent to the VBR™ (K003155) Made by Ulrich GmbH & Co, with respect to functional design, indications for use, principles of operation, and performance. The material used in the Sintea Biotech DSC/ALF Spinal System is the same as the predicate device.

Page 7 of 7



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Sintea Biotech  
% Mr. Gustavo A. Rios  
Regulatory Affairs  
407 Lincoln Rd. Suite 10L  
Miami Beach, Florida 33139

APR 18 2007

Re: K070181  
Trade/Device Name: DSC/ALF Spinal System  
Regulation Number: 21 CFR 888.3060  
Regulation Name: Spinal intervertebral body fixation orthosis  
Regulatory Class: Class II  
Product Code: MQP, KWQ  
Dated: January 8, 2007  
Received: January 19, 2007

Dear Mr. Rios:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at 240-276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K070181

Device Name: DSC/ALF Spinal System

### Indications for Use:

The DSC/ALF Spinal System is a vertebral body replacement system intended for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma (i.e., fracture). The DSC/ALF Spinal System is to be used with supplemental fixation. Specifically, the DSC/ALF system is to be used with the anterior lateral plate that is included in the system, and may also be used with the Sintea Biotech PLS Spinal System. Additionally, the DSC/ALF system is intended to be used with bone graft.

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)



(Division Sign-Off)

**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K070181

page 1 of 1