K07115 pag 1/2

510(k) SUMMARY

Bilok® ST Screw

MAY 17 2007

Applicant

Biocomposites Ltd

Keele Science Park

Keele

Staffordshire England ST5 5NL

Contact Person

Mr Simon Fitzer

Tel: +44 (0) 1782 338580 Fax +44 (0) 1782 338599 Email: sf@biocomposites.com

Classification Name:

Screw, fixation, bone

Common/Usual Name:

Bone screw

Trade/Proprietary Name

Bilok® ST Screw

Product Code

HWC

CFR Section

21CFR888.3040

Device Description

The Bilok® ST Screw is a cannulated, sterile, single use bone screw manufactured from a composite mix of calcium phosphate and poly L-lactic acid (PLLA).

Intended Use / Indications

The Bilok® ST Screw is indicated for use in anterior cruciate ligament (ACL) reconstruction procedures.

The Bilok® ST Screw is used to provide suspensary fixation during femoral fixation in ACL reconstruction using double looped (semitendinosis/gracilis) or quadruple (semitendinosis) graft.

P88242

Summary of Technology

The Bilok® ST Screw has the same technological characteristics as the legally marketed predicate device and any differences do not raise any concerns regarding safety and effectiveness.

Non Clinical Testing

Documentation provided demonstrates that the Bilok® ST Screw is substantially equivalent to the legally marketed predicate device and any differences do not raise any concerns regarding safety and effectiveness.

Substantial Equivalence

Documentation provided demonstrates that the Bilok® ST Screw is substantially equivalent to the legally marketed predicate device and any differences do not raise any concerns regarding safety and effectiveness.

Safety and Performance

Documentation provided demonstrates that the Bilok® ST Screw is substantially equivalent to the legally marketed predicate device and any differences do not raise any concerns regarding safety and effectiveness.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Biocomposites Ltd % Mr. Simon Fitzer Quality and Regulatory Affairs Manager Keele Science Park Keele, Staffordshire England, ST5 5NL

MAY 17 2007

Re: K071115

Trade/Device Name: Bilok® ST Screw Regulation Number: 21 CFR 888.3040

Regulation Name: Smooth or threaded metallic bone fixation fastener

Regulatory Class: Class II Product Code: HWC Dated: April 12, 2007 Received: April 20, 2007

Dear Mr. Fitzer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

Page 2 – Mr. Simon Fitzer

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Mark N. Melkerson

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if knov	vn): K071115	
Device Name:	Bilok® ST Screw	
Indications For Use:		
reconstruction procedu The Bilok®ST Screw is	res. sused to provide suspruction using double le	anterior cruciate ligament (ACL) ensary fixation during femoral coped (semitendinosis/gracilis) or
Prescription Use	OR part D)	Over-The-Counter use
PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED		
Concurre	nce of CDRH, Office of De	vice Evaluation (ODE)
(Division Signature) Division of (General, Restorati	Dana d at d
and Neurological Devices		
510(k) Num	ber_ <u>K071115</u>	