

---

**510 (k) SUMMARY- CryoSeal® FS Applicator System**

K071126

Submitter Name: ThermoGenesis Corp.  
Submitter Address: 2711 Citrus Road  
Rancho Cordova, CA 95742

AUG 16 2007

Contact Person: John Chapman PhD  
VP, Scientific Affairs  
Phone Number: 916.858.5132  
Fax Number: 916.858.5199

Date Prepared: April 19, 2007

Device Trade Name: CryoSeal® FS Applicator System

Classification Number: 21 CFR 880.5860  
Classification Name: Syringe, Piston  
Product Code: FMF

Predicate Devices: Baxter Healthcare Corp., DuploReach, K014088  
Baxter Healthcare Corp., Duploject Easy-Prep System, K020666

Statement of Intended Use: The CryoSeal® FS Applicator System is intended for the simultaneous application of the two components of CryoSeal FS (fibrin sealant) onto the incised liver surface in patients undergoing liver resection.

Device Description: The CryoSeal® FS Applicator System is composed of a thermoplastic polystyrene applicator handle, and two types of dispensing tips with a Y-body manifolds for application of the fibrin sealant: line/drop tips (5 and 10 cm lengths) and spray tips (3 configurations).

Comparison to the Predicate Devices: This device, with respect to material composition, device characteristics, principles of operation, and intended use, is substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

John Chapman, Ph.D.  
Vice President, Scientific Affairs  
ThermoGenesis Corporation  
2711 Citrus Road  
Rancho Cordova, California 95742

AUG 16 2007

Re: K071126  
Trade/Device Name: CryoSeal® FS Applicator System  
Regulation Number: 21 CFR 880.5860  
Regulation Name: Piston Syringe  
Regulatory Class: II  
Product Code: FMF  
Dated: July 27, 2007  
Received: July 30, 2007

Dear Dr. Chapman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Chiu Lin', with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known):

K071126

Device Name:

CryoSeal® FS Applicator System

**Indications for Use:**

The CryoSeal® FS Applicator System is intended for the simultaneous application of the two components of CryoSeal FS (fibrin sealant) onto the incised liver surface in patients undergoing liver resection.

Prescription Use   x    
(Part 21 CFR 801 Subpart D)

AND/OR


Over-The-Counter Use             
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*(Posted November 13, 2003)*

  
\_\_\_\_\_  
Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K071126