

510(k) Summary

K071181

As required by 807.92

MAY 30 2007

1. Company Identification

Konica Minolta Medical & Graphic, Inc.
2970 Ishikawa-machi, Hachioji-shi, Tokyo 192-8505, Japan
Phone: 81-42-660-9607
Fax: 81-42-660-9588

2. Official Correspondent

Koji Matsushima (Mr.)
Manager
Standards & Regulations Department, Quality Assurance Center

3. Date of Submission

April 26, 2007

4. Establishment Registration No.

3003769120

5. Device Trade Name

Direct Digitizer, REGIUS Model 110

6. Common Name

Computed Radiography Image Reader

7. Classification

Class II, 90 MQB, 21 CFR 892.1630.

8. Predicate Device

Direct Digitizer, REGIUS Model 190, 510(k) number: K052095

9. Description of Device

The Direct Digitizer, REGIUS Model 110 is an X-ray image reader which uses a stimuable phosphor plate (Plate) as X-ray detector installed in a separate cassette, and reads the image recorded on the Plate by inserting a cassette in the entrance slot of this device. By means of laser scan and photoelectric method, this device reads the X-ray image data created in form of a latent image on the Plate exposed by an external X-ray generating device, and converts the read data into digital. REGIUS Model 110 also reads the image

data of long areas of anatomy and to verify the position for a radiotherapy location.

The Regius 110 is really just a minor modification of our 510k cleared Direct Digitizer, REGIUS Model 190, 510(k) number: K052095. The features in the Regius 110 that are different than the Regius 190 are minor and are features that are present in the 510k cleared Model DD341 / REGIUS150, 510(k) number: K0990359 which we have included on the Comparison Table.

The modification was made to the Regius 110 to achieve a more compact design. To achieve the compact design, the Maximum Image Matrix Size and Sampling Pitch are decreased from the Size and Pitch in the REGIUS Model 190 and the transport is slightly different. However, these differences are the same technology as found in the 510k cleared Model DD341 / REGIUS150.

Please refer to section 2, Substantial Equivalence Comparison Table for details. The basic operations of REGIUS Model 110 such as a starting, a shut down, a registration-of-patient, a setting of a various condition are operated with the optional Image Processing Work Station, REGIUS CONSOLE CS-2000/3000 (510(k) cleared, K051523, July 20, 2005) and so on. Then the image data transfer to an externally connected device such as a host computer, an order input device, an image display device, a printer, an image data filing device, and other image reproduction devices.

This device is not intended for use with digital mammography system.

10. Intended Use

The Direct Digitizer, REGIUS Model 110 is an X-ray image reader which uses a stimuable phosphor plate (Plate) as X-ray detector installed in a separate cassette. It reads the image recorded on the Plate and transfers the image data to an externally connected device such as a host computer, an order input device, an image display device, a printer, an image data filing device, and other image reproduction devices. REGIUS Model 110 also reads the image data of long areas of anatomy and to verify the position for a radiotherapy location. It is designed and intended to be used by trained medical personnel in a clinic, a radiology department in a hospital and in other medical facilities.

This device is not intended for use with digital mammography system.

11. Substantial Equivalence to Predicate Device

The Direct Digitizer, REGIUS Model 110 is substantially equivalent to our Direct Digitizer REGIUS Model 190, 510(k) number: K052095. Comparison of the principal characteristics is shown in the Section 2.

12. Compliance Standard

Safety standard : IEC60601-1 Ed.2(1988)+ A1(1991)+A2(1995)
Electromagnetic Compatibility : IEC60601-1-2 Ed.2(2001)+A1(2004)
Radiation safety : 21 CFR 1040.10, IEC60825-1(1993)+A1(1997)+A2:2001

13. Conclusion

The Direct Digitizer, REGIUS Model 110 has the same intended use and basically the same technological characteristic as the predicate device which is approved 510(k) number: K052095. This Special 510(k) has demonstrated substantial equivalence as the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room – WO66-G609
Silver Spring, MD 20993-0002

Konica Minolta Medical & Graphic, Inc.
% Mr. Russell Munves
Official Correspondent
Storch, Amini & Munves, P.C.
140 East 45th Street, 25th Floor Two Grand Central Tower
NEW YORK NY 10017

AUG 23 2013

Re: K071181
Trade/Device Name: Direct Digitizer, REGIUS Model 110
Regulation Number: 21 CFR 892.1680
Regulation Name: Stationary x-ray system
Regulatory Class: II
Product Code: MQB
Dated: April 26, 2007
Received: April 30, 2007

Dear Mr. Munves:

This letter corrects our substantially equivalent letter of May 30, 2007.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of

medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely Yours,

A handwritten signature in black ink, appearing to read "Janine M. Morris", is written over the typed name and title.

Janine M. Morris
Acting Director
Division of Radiological Devices
Office of In Vitro Diagnostic Device
Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K071181

Device Name

:

Direct Digitizer, REGIUS Model 110

Indications for Use:

The Direct Digitizer, REGIUS Model 110 is an X-ray image reader which uses a stimuable phosphor plate (Plate) as X-ray detector installed in a separate cassette. It reads the image recorded on the Plate and transfers the image data to an externally connected device such as a host computer, an order input device, an image display device, a printer, an image data filing device, and other image reproduction devices. REGIUS Model 110 also reads the image data of long areas of anatomy and to verify the position for a radiotherapy location. It is designed and intended to be used by trained medical personnel in a clinic, a radiology department in a hospital and in other medical facilities.

This device is not intended for use with digital mammography system

Prescription Use ☒
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Reproductive, Abdominal, and
Radiological Devices

510(k) Number

K071181

Page 1 of _____