Diad Colour, Chiefficted Lowder Libertaine Examination Closes

\$70(k) SUMMARY \$67/740 BLUE COLOUR, CHLORINATED POWDER FREE NITRILE EXAMINATION GLOVES

OCT 4 2007		
Submitter's Name	MEDTEXX MANUFACTURING SDN. BHD.	
Submitter's Address	PT 4004, Kamunting Industrial Estate 34600 Taiping, Perak, Malaysia	
Submitter's Phone Number	605-891 1111 / 605-891 5555	
Submitter's Fax Number	605-891 1088	
Name of Contact Person	Ooi Loon Seng	
Date of Preparation		
Name of Device		
Trade Name :	BLUE COLOUR, CHLORINATED POWDER FREE NITRILE EXAMINATION GLOVES	
Common Name :	Nitrile Examination Gloves	
Classification Name :	Patient Examination Gloves	
Legally Marketed Device to which Equivalency is Being Claimed	Blue Colour, Chlorinated Powder Free Nitrile Examination Gloves as described in this 510 K Notification is substantially equivalent to the current Class 1 Patient Examination glove bearing the product code 80LZA (21 CFR 880.6250). It meets all the current specifications listed under the ASTM Specification D 6319 – 00a Standard Specification for Nitrile Examination Gloves for Medical Application.	
Description of the Device	The Blue Colour, Chlorinated Powder Free Nitrile Examination Gloves is substantially equivalent to the Class 1 patient examination glove bearing the product code 80LZA (21 CFR 880.6250). It	

	meets all the current specifications listed under the ASTM Specification D 6319 – 00a Standard Specification for Nitrile Examination Gloves for Medical Application ¹ . They are made from nitrile compound(dispersion of butadiene acrylonitrile copolymer. They are blue in color and are powder free.	
Intended Use of the Device	The Blue Colour, Chlorinated Powder Free Nitrile Examination Gloves are intended for single use for medical purposes that is worn on the hand of health care and similar personnel to prevent contamination between the health care personnel and the patients.	
Summary of Technological Characteristic Compared to the Predicate Device	There is no different technological characteristic. Gloves are made from nitrile compound (dispersion of butadiene acrylonitrile copolymer) and the initial products are powder free nitrile examination gloves.	
Brief Description of Non-Clinical Tests	Testing performed per ASTM D 6319 – 00a Standard Specification for Nitrile Examination Gloves for Medical Application and 21 CFR 800.20. Gloves meet all the current ASTM D 6319-00a. Primary skin irritation testing in the rabbit and delayed dermal contact sensitization study in the guinea pigs indicate no irritation or sensitization.	
Brief description of Clinical Tests	No new clinical tests were conducted under this 510(k).	
Conclusions Drawn from the Non-Clinical and Clinical Tests	Non-Clinical laboratory and animal based test data indicate that the powder free product meets all performance and biocompatibility requirements.	
Other Information Deemed Necessary by FDA	Not Applicable	







Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mdm. Ooi Loon Seng Regulatory Affairs Manager Medtexx Manufaxturing Sdn. Bhd. PT 4004 Kamunting Industrial Estate 34600 Taiping, Perak MALAYSIA

OCT 4 2007

Re: K071740

Trade/Device Name: Blue Colour, Chlorinated Powder Free Nitrile Examination

Gloves

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient Examination Glove

Regulatory Class: I Product Code: LZA

Dated: September 19, 2007 Received: September 24, 2007

Dear Mdm. Seng:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4618. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices Office of Device Evaluation Center for Devices and

Radiological Health

INDICATIONS FOR USE

An	n	lica	nt
, 10	\sim	IIYU	111

MEDTEXX MANUFACTURING SDN. BHD.

PT 4004, Kamunting Industrial Estate,

34600 Taiping, Perak,

Malaysia.

510(k) Number

(if known)

Device Name

BLUE COLOUR, CHLORINATED

POWDER FREE NITRILE EXAMINATION

GLOVES

Indications For Use

Blue Colour, Chlorinated Powder Free Nitrile Examination Glove is a single use device intended for medical purposes that is worn on the hand of healthcare and similar personnel to prevent contamination between the healthcare personnel and the patient.

Concurrence of CDRH Office of Device Evaluation (ODE)

Prescription Use Per 21 CFR 801.109

(Division Sign-Off)

OR Over-The-Counter

Division of Anesthesiology, General Hospital GCJ for 88 m

Infection Control, Dental Devices

510(k) Number: <u>K071740</u>