

ATTACHMENT III – 510(k) SUMMARY, REVISED**Mondeal® Hand Contour System
510(k) Summary of Safety and Effectiveness**

SEP 17 2007

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is K071797. This Summary was prepared on August 27, 2007

GENERAL INFORMATION**Manufacturer and Applicant Information:**

Mondeal Medical Systems GmbH
Moltkestr. 39
Tuttlingen, GERMANY 78532
Contact: Jay Evans
Telephone: 858-901-4123
Fax: 858-225-0311

Trade Name:

Device Name: Mondeal® Hand Contour System
Common Name: Screw, Fixation, Bone and Plate, Fixation, Bone
Classification Name: Plate, Fixation, Bone, Class II
and Regulation: 21 CFR 888.3030, 87HRS

Substantial Equivalence:

Mondeal Medical Systems GmbH, claims substantial equivalence to the Howmedica Profyle ® Hand and Small Fragment System, K961497.

Indications for Use:

The Mondeal® HAND CONTOUR System is intended to be used for the internal fixation of small bones including the hand and the foot.

Description of the Device

The Mondeal® Hand Contour System consists of titanium plates with shapes and sizes designed for internal fixation of small bones including the hand and foot, and screws of varying lengths from 4 to 23 mm and 1.2, 1.7, and 2.3 mm in diameter, supplied nonsterile packaged together in either tempered plastic or stainless steel trays suitable for recommended steam sterilisation, and also individually for single implantable use. The plates include straight, "T" shaped and condylar configurations. Manual reusable surgical instruments are supplied to facilitate implantation

ATTACHMENT III – 510(k) SUMMARY, REVISED

Characteristic	Mondeal® Hand Contour System (subject of this 510(k))	Howmedica Profyle Hand and Small Fragment System
Names	Mondeal® Hand Contour System	Howmedica Profyle Hand and Small Fragment System
510(k) number	K071797	K961497
Indications for use	The Mondeal® Hand Contour System is intended to be used for the internal fixation of small bones of the hand and foot.	Howmedica Profyle Hand and Small Fragment System is intended for use in internal fixation of small bones including the hand, foot, and cranialmaxillofacial skeleton.
Design and Device Characteristics		
Technology	Titanium plates and screws	Titanium plates and screws
Application	Internal Fixation of small bones of hand and foot	Internal Fixation of small bones of hand and foot
Design / Components	CP Titanium Grade 2 plates and Ti-6Al-4V ELI screws	CP Titanium Grade 2 plates and Ti-6Al-4V ELI screws
Performance Specifications		
Corrosion resistance	Identical	Identical
Mechanical properties	Similar hardness, yield and tensile strength, elongation, reduction in area, chemical content	Similar hardness, yield and tensile strength, elongation, reduction in area, chemical content
Sterilization Method	Steam Autoclave	Steam Autoclave
Packaging	Tempered plastic and or stainless steel trays suitable for steam sterilization including plate, screw, and hand tools (Class I) assortment, plates and screws also packaged individually, all non-sterils, intended for sterilization by purchaser	Tempered plastic and or stainless steel trays suitable for steam sterilization including plate, screw, and hand tools (Class I) assortment, plates and screws also packaged individually, all non-sterils, intended for sterilization by purchaser

Conclusions

Mondeal Medical Systems GmbH considers the Mondeal® HAND CONTOUR System to be substantially equivalent to the aforementioned predicate devices with regard to intended use, materials, biocompatibility, and overall performance characteristics in accordance with the above comparison summary.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

SEP 17 2007

Mondeal Medical Systems, GmbH
% Mr. Jay Evans
Mondeal North America, Inc.
13566 Freeport Road
San Diego, California 92133

Re: K071797
Trade/Device Name: Mondeal® Hand Contour System
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/multiple component metallic bone fixation appliances and accessories
Regulatory Class: Class II
Product Code: HRS
Dated: August 28, 2007
Received: September 4, 2007

Dear Mr. Evans:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number: K071797

Device Name: Mondeal® Hand Contour System

Indications For Use:

The Mondeal® HAND CONTOUR System is intended to be used for the internal fixation of small bones including the hand and the foot.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of General, Restorative,
and Neurological Devices

510(k) Number K071797

Page 1 of 1