

SEP 14 2007

510(k) SUMMARY

"This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92."

"The assigned 510(k) number is: K072246."

Submitter: Maine Standards Company
Address: 765 Roosevelt Trail
Windham, ME 04062
Telephone: 207-892-1300
Fax: 207-892-2266
Contact: Holly A. Cressman, Mgr. QA/RA

Summary prepared on: August 01, 2007

Device classification name: Multi-Analyte Controls, All Kinds (Assayed and Unassayed)
Device description: Quality control material (assayed and unassayed)
Proprietary Name: VALIDATE® SP1 Calibration Verification Test Set
Regulation Number: 21 CFR 862.1660
Product Code: JJY
Regulatory Class: Class I

Predicate Device:

Cliniqa LiniCAL® Calibration Verifier Protein 1 (K031575) and Cliniqa LiniCAL® Calibration Verifier Protein 2 (K031577) Cliniqa, San Marcos, California.

Device description: VALIDATE® SP1 Calibration Verification Test Sets are human serum calibration verification / linearity materials containing multiple levels used to establish the relationship between theoretical operation and actual performance of the included analytes. Each test set consists of one bottle each of five (5) levels. Each bottle of Levels 1 through 5 contains 1.0 milliliter. There exists a linear relationship between Levels 1 through 5.

Intended use: The VALIDATE® SP1 Calibration Verification Test Set solutions are intended for *in vitro* diagnostic use in the quantitative determination of linearity, calibration verification and verification of reportable range in automated, semi-automated and manual chemistry systems for the following analytes: Immunoglobulin A (IgA), Immunoglobulin G (IgG), Immunoglobulin M (IgM), Complement C3 (C3), Complement C4 (C4), α 1-Antitrypsin (AAT), and Transferrin (TRF). VALIDATE® SP1 Calibration Verification Test Set solutions are not intended for use as routine quality control materials or as calibration materials.

Summary:

The information provided in this pre-market notification demonstrates that the performance of VALIDATE® SP1 Calibration Verification Test Sets is substantially equivalent in form and function to the predicate devices for its stated intended use.



SEP 14 2007

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Maine Standards Co.
c/o Ms. Holly Cressman
Manager, QA/RA
765 Roosevelt Trail
Windham, ME 04062-5365

Re: k072246

Trade/Device Name: VALIDATE® SP1 Calibration Verification Test Set
Regulation Number: 21 CFR 862.1660
Regulation Name: Quality control material (assayed and unassayed)
Regulatory Class: Class I
Product Code: JJY
Dated: August 08, 2007
Received: August 13, 2007

Dear Ms. Cressman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The

Page 2 –

FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Robert L. Becker, Jr.", is positioned above the typed name.

Robert L. Becker, Jr., M.D., Ph.D.
Director

Division of Immunology and Hematology Devices
Office of In Vitro Diagnostic Device Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number: K072246

Device Name: VALIDATE® SP1 Calibration Verification Test Set

Indications For Use:

VALIDATE® SP1 Calibration Verification Test Set solutions are intended for *in vitro* diagnostic use in the quantitative determination of linearity, calibration verification and verification of reportable range in automated, semi-automated and manual chemistry systems for the following analytes: Immunoglobulin A (IgA), Immunoglobulin G (IgG), Immunoglobulin M (IgM), Complement C3 (C3), Complement C4 (C4), α 1-Antitrypsin (AAT), and Transferrin (TRF).

Prescription Use ☒
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Maria M Chan
Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety
510(k) K072246

Page 1 of _____