WUNIEULLUV 9 4140 K 072337

NO. 225, YUAN-PIER ST., HSIN CHU CITY, CHINA (TAIWAN) TEL: 886-3-5382105 FAX: 886-3-5382191

Homepage: www.wustech.com.tw Email: wustis@ms45.hinet.net

SEP 2 1 2007

"__510(k) SUMMARY"

Submitter's Name: WU'S TECH CO., LTD.

NO. 225, YUAN-PIER ST., HSIN CHU CITY, 30093, CHINA (TAIWAN)

Date summary prepared:

August 6, 2007

Device Name:

Proprietary Name:

WU'S 4-WHEELED NEO SCOOTER, WT-M4A

Common or Usual Name:

SCOOTER

Classification Name:

MOTORIZED 3-WHEELED VEHICLE, Class II,

21 CFR 890.3800

Indications for Use:

The device is intended for medical purposes to provide mobility to persons restricted to a seated position.

Description of the device:

The WU'S SCOOTER WT-M4A is an indoor / outdoor electric scooter that is battery operated. It has a base with four-wheeled with a seat, armrests, and a front basket. The movement of the scooter is controlled by the rider who uses hand controls located at the top of the steering column. The device can be disassembled for transport and is provided with an onboard battery charger.

Performance Testing:

EMC Report ANSI / RESNA WC/Vol.2-1998, CISPR 11: 1990, EN61000-3-2: 1995, IEC61000-3-3: 1995 (Electrically powered wheelchairs, scooters, and their chargers – requirements and test methods)

Legally marketed device for substantial equivalence comparison:

WU'S SCOOTER WT-M4 (K023166)

THUR LUCLE VUO LLU O

NO. 225, YUAN-PIER ST., HSIN CHU CITY, CHINA (TAIWAN) TEL: 886-3-5382105 FAX: 886-3-5382191

Homepage: www.wustech.com.tw Email: wustis@ms45.hinet.net

Summary for substantial equivalence comparison:

The intended use between the two devices is the same. Mainframes of two devices are fixed. Mainframes materials of the two devices all meet the strength and fatigue tests and they are similar for the material aspects. The overall dimensions are similar. The two devices use the same suspension of cross brace, same size of wheels, and same type of armrest, footplate, wheel lock, and same warranty. Besides, the back upholstery materials are also the same fabric and passed the resistance ignition test.

Especially the electronic systems between two devices are the same suppliers, and all passed by the UL certificated, for instance the electronic controller, batteries, motors, and recharge. Thus the same safety level for the two devices is assured.

Owing to the new device WT-M4A is huger than the predicate device WT-M4 that means the weight capability of the new device is 205kgs and the predicate device is 130kgs. Thus the major differences exiting the two devices are the cruising range per charge and incline angle. This means the cruising range of the new device is 10~15 miles and the incline angle is 8 degrees, and the predicate device is 19 miles and 12 degrees for the incline angle. Certainly the real range depends on the practical environments, i.e., weight, surface, incline, and temperature. For the real life use, the two devices are substantially equivalent.







SEP 2 1 2007

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

WU'S TECH Co., Ltd. % ROC Chinese-European Industrial Research Society Dr. Jen, Ke-Min No. 58, Fu-Chiun St. Hsin-Chu City, 30067, Taiwan, ROC

Re: K072337

Trade/Device Name: Neo Scooter WT-M4A Regulation Number: 21 CFR 890.3800

Regulation Name: Motorized three-wheeled vehicle

Regulatory Class: Class II

Product Code: INI Dated: August 6, 2007 Received: August 20, 2007

Dear Dr. Jen, Ke-Min:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Dr. Jen, Ke-Min

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120 Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Mark N. Melkerson

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Indications for Use

| 510 (K) Number (If Known): <u>K</u> | |
|--|---|
| Device Name: <u>WU'S NEO SCOO</u> | TER, WT-M4A |
| | |
| Indications for Use: | |
| The device is intended for medical purpose | s to provide mobility to persons restricted to |
| a sitting position. | |
| | |
| | |
| | |
| | |
| | |
| | |
| Prescription Use AND/OI | R Over-The-Counter Use |
| (Part 21 CFR 801 Subpart D) | (21 CFR 807 Subpart C) |
| | |
| | |
| (PLEASE DO NOT WRITE BELOW THIS LIN | E-CONTINUE ON ANOTHER PAGE IF NEEDED) |
| | |
| | |
| | ffice of Device Evaluation (ODE) |
| (Division'\S | a T aran Sanara and a |
| · | General, Restorative |
| and Neuro | logical Devices Page 1 of 1 |
| 518(k) Nur | $(-c \circ) ? \mathcal{L} \mathcal{L}$ |