

**510(k) Summary - TxControl**K072615  
page 1 of 2

Date Prepared 10<sup>th</sup> September 2007

Official Contact Dr Lionel King  
V.P., Quality Assurance & Regulatory Affairs  
ResMed Ltd  
1 Elizabeth Macarthur Drive  
Bella Vista, NSW 2153  
Australia  
Tel: +61 (2) 8884 1000  
Fax: +61 (2) 8884 2021

Classification Reference 21 CFR 868.5905

Product Code 73 BZD

Common/Usual Name Noncontinuous ventilator (IPPB).

Proprietary Name TxControl

Predicate Device(s) ResControl II (K040944)

Reason for submission New Device

JAN - 4 2008

## Intended Use

The TxControl™ is a software application intended to be used by clinicians with ResMed flow generators that incorporate ResMed's proprietary communication protocol. TxControl provides real-time data display directly from the flow generator or via ResControl II™.

TxControl can also provide flow generator setting changes, remotely within a clinical environment.

## Device Description

ResMed's TxControl™ is a PC-based software application that enables clinicians to monitor real-time patient and flow generator information and adjust flow generator settings as required from the control room within the sleep lab clinical setting.

TxControl is also designed to transfer data from a CPAP or bilevel flow generator to a polysomnograph (PSG) display system such as Somnologica (K971813) when used in conjunction with a ResControl II device. The performance and functional characteristics of the TxControl includes similar user friendly features of the predicate device, ResControl II (K040944).

## Substantial Equivalence

The new device has the following similarities to the previously cleared predicate device.

- Similar intended use
- Similar operating principle
- Similar technologies
- Simplified manufacturing process

Design and Verification activities were performed on the TxControl as a result of the risk analysis and product requirements. All tests confirmed the product met the predetermined acceptance criteria. ResMed has determined that the new device is Substantially Equivalent to the predicate device. The transfer of ResControl II hardware features to the TxControl software application has not altered the safety and effectiveness when used primarily in the management of patients with Obstructive Sleep Apnea (OSA). The new device complies with the applicable standards and requirements referenced in the FDA guidance documents:

- FDA Draft Reviewer Guidance for Ventilators (July 1995)
  - FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices (May 11, 2005)
  - FDA Off-the-Shelf Software Use in Medical Devices (September 9, 1999)
  - FDA Guidance for Industry - Cybersecurity for Networked Medical Devices Containing Off-the-Shelf (OTS) Software (January 14, 2005)
-



JAN - 4 2008

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

ResMed, Limited  
C/O Mr. David D'Cruz  
Vice President Clinical & Regulatory Affairs  
ResMed Corporation  
14040 Danielson Street  
Poway, California 92064-6857

Re: K072615  
Trade/Device Name: TxControl™  
Regulation Number: 21 CFR 868.5905  
Regulation Name: Noncontinuous Ventilator (IPPB)  
Regulatory Class: II  
Product Code: BZD  
Dated: December 18, 2007  
Received: December 26, 2007

Dear Mr. D'Cruz:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin, Ph.D.", is positioned above the printed name.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

### Indication for Use

510(k) Number (if known):

Device Name: TxControl™

Indication for Use

The TxControl™ is a software application intended to be used by clinicians with ResMed flow generators that incorporate ResMed's proprietary communication protocol. TxControl provides real-time data display directly from the flow generator or via ResControl II™.

TxControl can also provide flow generator setting changes, remotely within a clinical environment.

Prescription Use   X  

AND/OR

Over-The-Counter Use       


(Part 21 CFR 801 Subpart D)

(Part 21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of Device Evaluation (ODE)

Page 1 of   1  

  
(Division Sign-Off) *acting B.S.P.*  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number:   K072615  

*for*