

Premarket Notification
510(k) Summary
(As Required by 21 CFR 807.93)

K072751

This 510(k) Summary of safety and effectiveness for the New Star Model CoolTouch NS160 CoolLipo Nd:YAG Surgical Laser system is submitted in accordance with the requirements of the SMDA 1990 and following guidance concerning the organization and content of a 510(k) summary.

Submitter: New Star Lasers, Inc. d.b.a. CoolTouch, Inc.

Address: 9085 Foothills Boulevard
Roseville, CA 95747

Contact Person: Natalie Vollrath
Quality Assurance Manager

JAN 02 2007

Telephone/Fax/Email: (916) 677-1912 – Phone
(916) 677-1901 – Fax
nvollrath@newstarlasers.com - Email

Date prepared: October 26, 2007

Device Trade Name: CoolTouch NS160 CoolLipo Nd:YAG Surgical Laser

Common Name: Nd: YAG Surgical Laser

Classification Name: Instrument, surgical, powered, laser
79-GEX
21 CFR §878.4810

Legally Marketed Predicate Devices: Cynosure SmartLipo Nd:YAG laser system
K062321

Device Description: The CoolTouch NS160 CoolLipo Nd:YAG Surgical Laser System is an Nd:YAG laser producing laser emission at 1320 nm. The laser consists of a cabinet which houses the power supply, the cooling system, microcontroller, laser, foot switch, and the fiber optic for delivery of the laser energy with microcannula setup.

Intended Use: The Cooltouch NS160 Cool Lipo is indicated for laser-assisted lipolysis.

Comparison:	The Cooltouch NS160 CoolLipo has the same indication for use, the same principle of operation, and essentially the same wavelength and pulse energy rate as the predicate device.
Nonclinical Performance Data	Nonclinical performance data produced results that indicate the Cooltouch NS160 CoolLipo laser system is effective for laser-assisted lipolysis.
Clinical Performance Data:	None
Conclusion:	The CoolTouch NC160 CoolLipo Nd:YAG Surgical Laser System is substantially equivalent to the predicate device and is indicated for laser-assisted lipolysis.
Additional Information:	None requested at this time.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JAN 02 2007

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

New Star Lasers, Inc.
% Ms. Natalie R. Vollrath
Quality Assurance Manager
9085 Foothills Boulevard
Roseville, California 95747

Re: K072751

Trade/Device Name: CoolTouch NS160 CoolLipo Nd:YAG Laser System
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general and plastic surgery and
in dermatology
Regulatory Class: II
Product Code: GEX
Dated: September 26, 2007
Received: September 27, 2007

Dear Ms. Vollrath:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,




Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number Pending K072751

Device Name CoolTouch NS160 CoolLipo Nd:YAG Laser System

Indications for Use In addition to previously cleared indications for use, the CoolTouch Model NS160 CoolLipo Nd:YAG Surgical Laser is indicated for laser-assisted lipolysis.

 FOR M. MELKERSEN
(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

510(k) Number K072751

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-the-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)