

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR Part 807.93

| | |
|------------------------------|--|
| Submitter | EBI Inc. #289, Wuisong-ri, Apyang-myeon Kyeongsan-city, Kyeongsangbook-do 712-825, Korea Telephone 82-53-817-7767 Fax: 82-53-817-7768 |
| Contact | Myung-bae Jegal, Managing Director |
| Date Prepared | February 21, 2008 |
| Device Name | EBI Internal Implant |
| Classification Name | Implant, endosseous, root-form DZE and Abutment implant, dental, endosseous NHA |
| Device Classification | Class II Dental Devices Panel 21 CFR 872.3640 and 872.3630 |
| Predicate Devices | <ul style="list-style-type: none">• AVANA Dental Implant System K051576• ITI Dental implant system K003552 |
| Performance | Performance standards have not been established by the FDA under Section 514 of the Federal Food, Drug, and Cosmetic Act. |
| Device Description | The EBI Internal Implant system is a threaded, root form endosseous implant comprised of an internal octagon, 8° morse tapered device manufactured from pure Titanium ASTM F-67 and the Abutment from Titanium 6A-4V alloy E 11 –ASTM 136. The implant's external threaded surfaces are roughened to facilitate tissue and bone integration. The implant's self-tapping feature may be used with or without pre-tapping the bone. System components for restorative purposes include: screw retained and cement retained abutments, cover screws, healing abutments, and angled abutments. |
| Indications for Use | EBI dental implants are intended for immediate, delayed, or conventional placement in the maxillary and/or Mandibular arches to support crowns, bridges or overdentures in edentulous or partially edentulous patients. |

Technical Characteristics

| Characteristics | Predicate Devices | | |
|-----------------------------|--|---|---|
| | EBI Internal Implant System | AVANA Dental Implant K051576 OSSTEM Co., Ltd. | ITI Dental Implant K003552 Institut Straumann AG |
| 510(k) | - | | |
| Material composition | Implant: Pure Titanium Abutment: Ti Alloy | Implant: Pure Titanium Abutment: Ti Alloy | Implant: Pure Titanium Abutment: Ti Alloy |
| Design | Internal Morse Taper, Internal Octagon | Internal Connection Morse Taper, Internal Octagon | Internal Connection Morse Taper, Internal Octagon |
| Implant Diameter | 3.3-4 mm | 4.1-4.8 mm | 3.3 – 4 mm |
| Implant Length | 6-15 mm | 8.4-15 mm | 6 – 14mm |
| Sterilization Method | Gamma | Gamma | Gamma |
| Indications for use | Intended for immediate delayed or conventional placement in the maxillary and/or Mandibular arches to support crowns, bridges, or overdentures in edentulous or partially edentulous patients. | Intended for use in partially or fully edentulous mandibles and maxillae, in support of single or multiple-unit restorations including cement retained, screw retained, or overdenture restorations, and terminal or intermediate abutment support of fixed bridgework. | Intended for immediate, delayed, or conventional placement in the maxillary and/or Mandibular arches to support crowns, bridges or overdentures in edentulous or partially edentulous patients. |
| Surface treatment | Resorbable Blasting media (RBM) | Resorbable Blasting media (RBM) | Resorbable Blasting media (RBM) |

Conclusion: In all respects, the EBI Internal Implants components are the equivalent of currently marked devices. They are made of the same materials and have similar dimensions and characteristics. This device is substantially equivalent in design, material, intended use and function to the products listed in the table above.



FEB 28 2008

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

EBI, Incorporated
C/O Mr. Daniel Kamm
Principal Consultant
Kamm & Associates
P.O. Box 7007
Deerfield, Illinois 60015

Re: K073116
Trade/Device Name: EBI Internal Implant System
Regulation Number: 872.3640
Regulation Name: Endosseous Dental Implant Abutment
Regulatory Class: II
Product Code: DZE, NHA
Dated: February 21, 2008
Received: February 25, 2008

Dear Mr. Kamm:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin, Ph.D.", with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K073116

10/1

Indications for Use

510(k) Number (if known): K07

Device Name: EBI Internal Implant System

Indications For Use: Intended for immediate, delayed, or conventional placement in the maxillary and/or mandibular arches to support crowns, bridges or overdentures in edentulous or partially edentulous patients.

Prescription Use ☒
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Runner

(Division Sign-Off)

Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

Page 1 of 1

510(k) Number: K073116