

K073286 1/2

Summary of Safety and Effectiveness

MAR - 7 2008

Submitter: Zimmer, Inc.
P.O. Box 708
Warsaw, IN 46581-0708

Contact Person: Brandon Hipsher, RAC
Associate Manager, Corporate Regulatory Affairs
Telephone: (574) 371-8083
Fax: (574) 372-4605

Date: March 3, 2008

Trade Name: *Gender Solutions*[™] *Natural-Knee*[®] Flex System

Common Name: Total Knee Prosthesis

Classification Name and Reference: Knee joint patellofemorotibial
polymer/metal/polymer semiconstrained cemented
prosthesis
21 CFR § 888.3560

Knee joint patellofemorotibial metal/polymer
porous-coated uncemented prosthesis
21 CFR § 888.3565

Knee joint patellofemoral polymer/metal
semiconstrained cemented prosthesis
21 CFR § 888.3540

Predicate Devices: *Gender Solutions Natural-Knee* Flex System,
manufactured by Zimmer, Inc., K070214, cleared
March 16, 2007.
Zimmer Patellofemoral Joint Prosthesis,
manufactured by Zimmer, Inc., K070695, cleared
June 7, 2007.

Device Description: The *Gender Solutions Natural-Knee* Flex (N-K
Flex) components are semiconstrained, nonlinked
condylar knee prosthesis that are designed to have a
maximum active flexion of 155 degrees. The N-K
Flex femoral provides increased flexion capability

for patients who have both the flexibility and desire to increase their flexion range.

Intended Use:

N-K Flex System:

- Components with *CSTi* porous coating are indicated for uncemented or cemented use in skeletally mature individuals with intact medial and lateral collateral ligaments undergoing primary surgery for rehabilitating knees damaged as a result of Noninflammatory Degenerative Joint Disease (NIDJD) or Inflammatory Joint Disease (IJD).
- Components without *CSTi* porous coating are indicated for cemented use in skeletally mature individuals with intact medial and lateral collateral ligaments with conditions of Noninflammatory Degenerative Joint Disease (NIDJD) or Inflammatory Joint Disease (IJD), correctable varus-valgus deformity and moderate flexion contracture, or failed previous surgery where pain, deformity or dysfunction persists.

Comparison to Predicate Device:

Except for dimensional and material modifications, the subject components are identical to the predicate device.

Performance Data (Nonclinical and/or Clinical):

Non-Clinical Performance and Conclusions:

Mechanical testing of the subject devices demonstrates that they are substantially equivalent to the predicate devices.

Clinical Performance and Conclusions:

Clinical data and conclusions were not needed for these devices.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Zimmer, Inc.
% Mr. Brandon Hipsher
P.O. Box 708
Warsaw, Indiana 46581-0708

MAR - 7 2008

Re: K073286

Trade/Device Name: *Gender Solutions™ Natural-Knee®* Flex System

Regulation Number: 21 CFR 888.3560

Regulation Name: Knee joint patellofemorotibial polymer/metal/polymer semi-constrained cemented prosthesis

Regulatory Class: II

Product Code: JWH, MBH, KRR

Dated: March 3, 2008

Received: March 5, 2008

Dear Mr. Hipsher:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Brandon Hipsher

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K073286

Device Name:

Gender Solutions[™] *Natural-Knee*[®] Flex System

Indications for Use:

Gender Solutions[™] *Natural-Knee*[®] Flex System:

- Components with *CSTi* porous coating are indicated for uncemented or cemented use in skeletally mature individuals with intact medial and lateral collateral ligaments undergoing primary surgery for rehabilitating knees damaged as a result of Noninflammatory Degenerative Joint Disease (NIDJD) or Inflammatory Joint Disease (IJD).
- Components without *CSTi* are indicated for cemented use in skeletally mature individuals with intact medial and lateral collateral ligaments with conditions of Noninflammatory Degenerative Joint Disease (NIDJD) or Inflammatory Joint Disease (IJD), correctable varus-valgus deformity and moderate flexion contracture, or failed previous surgery where pain, deformity or dysfunction persists.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(Please do not write below this line – Continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neil R. Olsen
(Division Sign-Off)

Page 1 of 1

**Division of General, Restorative,
and Neurological Devices**

510(k) Number K073286