

FEB 22 2008

K073322

510(k) SUMMARY – Fox Q-810, Q-980 and Q-1064 Laser

Applicant Name: A.R.C. Laser GmbH
Bessemerstr. 14, D-90411 Nurnberg, Germany

Contact Person: Reinhardt Thyzel, President

Date Prepared: February 8, 2008

Device Trade Name: Fox Q-810, Q-980 and Q-1064 Lasers

Device Common Name: Diode Laser

Classification Name: Laser Surgical Instrument

Predicate Devices: Fox Q-810 (K062619), Ceralas 980 (K072106) and Laserscope 1064 (K990903)

Device Description: Fox Q-810, Q-980 and Q-1064 are standard diode medical lasers with 810, 980nm and 1064nm wavelength, respectively

Intended Use: Fox Q-810:
indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology , general surgery, neurosurgery, otolaryngology, ophthalmology, and pulmonology .

Fox Q-980:
indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, neurosurgery, otolaryngology, orthopedics,

ophthalmology, pulmonology, and thoracic surgery.

Fox Q-1064:

indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, general surgery, genitourinary, neurosurgery, otolaryngology, orthopedics, ophthalmology, pulmonology, and thoracic surgery.

Device Technological
Characteristics and
Comparison to Predicate
Device(s):

The Fox Q-980 uses diodes to generate energy in the 980 nm range. (Commercially available) Fibers deliver energy to the tissue. The Fox Q-1064 Laser is the same system, but generates energy in the 1064 nm range. The new Fox Q-810 is also the same, producing energy in the 810 nm range.

Predicates: The Fox Q-810 Laser (cleared under K062619) is the same system, but generates energy in the 810 nm range (at a lower output than the new 810 nm device). The Ceralas 980 is also a diode laser producing energy in the 980 nm range, and the Laserscope system produces energy in the 1064 nm range.

Performance Standards:

The Fox Q-810, Q-980 and Q-1064 Lasers comply with the performance requirements of 21CFR 1040.10 and 1040.11, with permissible deviations defined in Laser Notice 50, dated July 26, 2001. The diode lasers also comply with IEC 60601-1:1998 including amendment 1, IEC 60601-2-22:1995, and IEC 60825-1:1993 including amendments 1 and 2.

Conclusion:

The Fox 810, 980 and 1064 Lasers are substantially equivalent to the predicate devices. They have similar intended uses and comply with the same safety and performance standards.



FEB 22 2008

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

A.R.C. Laser GmbH
% PPD Medical Device
Ms. Kirsten H. Paulson
3202 Tower Oaks Blvd.
Rockville, MD 20852

Re: K073322
Trade/Device Name: Fox 1-980, Fox Q-1064, Fox Q-810
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general and
plastic surgery and in dermatology
Regulatory Class: II
Product Code: GEX
Dated: November 21, 2007
Received: November 26, 2007

Dear Ms. Paulson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K073322

Device Name: Fox Q-810, Q-980 and Q-1064 Diode Laser

Indications for Use:

Fox Q-810:

indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, neurosurgery, otolaryngology, ophthalmology, and pulmonology.

Fox Q-980:

indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, gastroenterology, general surgery, genitourinary, gynecology, neurosurgery, otolaryngology, orthopedics, ophthalmology, pulmonology, and thoracic surgery.


Fox Q-1064:

indicated for surgical applications requiring the ablation, vaporization, excision, incision, hemostasis, or coagulation of soft tissues in medical specialties including dermatology, dentistry, general surgery, genitourinary, neurosurgery, otolaryngology, orthopedics, ophthalmology, pulmonology, and thoracic surgery.

Prescription Use X AND/OR Over-The-Counter Use
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices
510(k) Number K073322