

510(K) Summary
PhotoNova Family of Pulsed Light Systems

K073477

This 510(K) Summary of safety and effectiveness for the Photonova Family of Pulsed Light Systems is submitted in accordance with the requirements of the SMDA 1990 and following guidance concerning the organization and content of a 510(K) summary.

Applicant: Photonova of Sweden AB

Address: Hulda Lindgrens Gata 6C
42131 Vastra Frolunda
SWEDEN

Contact Person: Anders Bonde

MAY - 9 2008

Telephone / Fax / Email 46-31-742-1281- Phone
abonde@photonova.com

Preparation Date: October 31, 2007

Device Trade Name: Photonova Family of Pulsed Light Systems

Common Name: Intense Pulsed Light

Classification Name: Instrument, Surgical, Powered, laser
79-GEX, 21 CFR 878-48

Legally Marketed Predicate Device: OmniLight K032191
Lux V K040081
StarLux K03349

Description of the Photonova Family of Pulsed Light Systems

The Photonova Pulsed Light Systems delivers pulsed light at wavelengths starting at 515 nanometers. The device consists of three interconnected sections: The cabinet which houses the power supply, the cooling system and the microcontroller, the umbilical to the handpiece, and the handpiece, which houses the waveguide

Intended use of the Photonova Family of Pulsed Light Systems

The Photonova Pulsed Light System is indicated for the following:

- The removal of unwanted hair from all skin types and to effect stable long-term or permanent hair reduction
- The treatment of benign pigmented lesions, including lentigines, nevi, melasma and café-au-lait
- The treatment of vascular lesions, including port wine stains, hemangiomas, angiomas, telangiectasias, rosacea, facial and leg veins
- Treatment of inflammatory acne (acne vulgaris)

Performance Data:

None

510(K) Summary

PhotoNova Family of Pulsed Light Systems

Conclusion:

The PhotoNova Family of Pulsed Light Systems is substantially equivalent to other existing pulsed light systems in commercial distribution in Dermatology and Plastic Surgery.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY - 9 2008

Photonova of Sweden AB
% Ms. Connie Hoy
908 Stetson Street
Woodland, California 95776

Re: K073477

Trade/Device Name: PhotoNova Family of Pulsed Light Systems

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser surgical instrument for use in general and plastic surgery and
in dermatology

Regulatory Class: II

Product Code: GEX

Dated: April 30, 2008

Received: May 5, 2008

Dear Ms. Hoy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Connie Hoy

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Mark N. Melkerson
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Attachment 4
Indications For Use Statement

510(k) Number (if Known): K073477

Device Name: **PhotoNova Family of Pulsed Light Systems**

Indications for Use:

The PhotoNova Family of Pulsed Light Systems is designed for use in dermatology applications including:

- The removal of unwanted hair from all skin types and to effect stable long-term or permanent hair reduction
- The treatment of benign pigmented lesions, including lentigines, nevi, melasma and café-au-lait
- The treatment of vascular lesions, including port wine stains, hemangiomas, angiomas, telangiectasias, rosacea, facial and leg veins
- Treatment of inflammatory acne (acne vulgaris)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED)


Concurrence of Neil H. Dyl for for me
(Division Sign-Off) Office of Device Evaluation (ODE)

**Division of General, Restorative,
and Neurological Devices**

510(k) Number K073477

Prescription Use

X

OR

Over-The-Counter Use