

**SECTION IV****510(k) SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION**

as required by the Safe Medical Devices Act of 1990 and codified in 21 CFR 807.92 upon which the substantial equivalence is based.

**Smith & Nephew TWINFIX PK FP Suture Anchor**

Date Prepared: December 12, 2007

MAR - 4 2008

**A. Submitter's Name:**

Smith & Nephew, Inc., Endoscopy Division  
150 Minuteman Road  
Andover MA, 01810

**B. Company Contact**

Deana Boushell  
Principle Regulatory Affairs Specialist  
Phone: (508) 337-4036  
FAX: (508) 261-3620

**C. Device Name**

Trade Name: TWINFIX PK FP Suture Anchor  
Common Name: Fastener, fixation, non-degradable, soft tissue  
Classification Name: Fastener, fixation, non-degradable, soft tissue

**D. Predicate Devices**

The Smith & Nephew TWINFIX PK FP Anchor is substantially equivalent in Intended Use and Fundamental Scientific Technology to the following legally marketed device in commercial distribution: Arthrex® PushLock™

**E. Description of Device**

The TWINFIX PK FP Anchor is a suture anchor manufactured from PEEK polymer. The pound in anchor is not preloaded with suture, incorporates a two piece anchor body and facilitates arthroscopic repair of soft tissue to bone. The

design allows the surgeon to adjust the tension on the tissue intraoperatively and secure the repair.

#### F. Intended Use

The Smith & Nephew TWINFIX PK FP Suture Anchor is intended for use for the reattachment of soft tissue to bone for the following indications:

**Shoulder:** Rotator cuff repair, Bankart repair, Slap lesion repair, Biceps tenodesis, Acromio-Clavicular separation, Deltoid repair, and Capsular shift or Capsulolabral reconstruction.

**Foot/Ankle:** Lateral stabilization, Medial stabilization, Achilles tendon repair, Hallux valgus reconstruction, Mid-foot reconstruction, Metatarsal ligament repair.

**Knee:** Medial collateral ligament repair, Lateral collateral ligament repair, Patellar tendon repair, Posterior oblique ligament repair, Illoiotibial band tenodesis.

**Hand/Wrist:** Scapholunate ligament reconstruction, Ulnar collateral ligament reconstruction, Radial collateral ligament reconstruction.

**Elbow:** Biceps Tendon reattachment, Ulnar or radial collateral ligament reconstruction.

#### G. Comparison of Technological Characteristics

The Smith & Nephew TWINFIX PK FP Suture Anchor is substantially equivalent in design, materials, function and intended use to the Arthrex® PushLock™ suture anchor, cleared in K051219. The proposed and the predicate devices both have the same intended use, indications for use, anchor material.

#### H. Summary Performance Data

The performance testing conducted demonstrates substantial equivalence to the Arthrex® PushLock™ suture anchor, cleared in K051219. The testing also demonstrates that the differences in the new device and the predicate device do not raise any new issues of safety and efficacy.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Smith & Nephew, Inc., Endoscopy Division  
% Ms. Deana Boushell  
Regulatory Affairs Specialist  
150 Minuteman Road  
Andover, MA 01810

MAR - 4 2008

Re: K073509  
Trade/Device Name: Twinfix PK FP Anchor  
Regulation Number: 21 CFR 888.3040  
Regulation Name: Smooth or threaded metallic bone fixation fastener  
Regulatory Class: Class II  
Product Code: HWC, MBI, JDR  
Dated: December 12, 2007  
Received: December 13, 2007

Dear Ms. Boushell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Deana Boushell

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (240) 276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at (240) 276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at toll-free number (800) 638-2041 or (240) 276-3150 or the Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson", with a stylized flourish at the end.

Mark N. Melkerson  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Page 141

§1.D

**Indications for Use**

510(k) Number (if known): K073509

Device Name: Smith & Nephew TWINFIX PK FP Suture Anchor

**Indications For Use:**

The Smith & Nephew TWINFIX PK FP Suture Anchor is intended for use for the reattachment of soft tissue to bone for the following indications:

**Shoulder:** Rotator cuff repair, Bankart repair, Slap lesion repair, Biceps tenodesis, Acromio-Clavicular separation, Deltoid repair, and Capsular shift or Capsulolabral reconstruction.

**Foot/Ankle:** Lateral stabilization, Medial stabilization, Achilles tendon repair, Hallux valgus reconstruction, Mid-foot reconstruction, Metatarsal ligament repair.

**Knee:** Medial collateral ligament repair, Lateral collateral ligament repair, Patellar tendon repair, Posterior oblique ligament repair, Illoiotibial band tenodesis.

**Hand/Wrist:** Scapholunate ligament reconstruction, Ulnar collateral ligament reconstruction, Radial collateral ligament reconstruction.

**Elbow:** Biceps Tendon reattachment, Ulnar or radial collateral ligament reconstruction.

Prescription Use   x  

AND/OR

Over-The-Counter Use   No  

(Per 21 CFR 801 Subpart D)

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neil R. P. [Signature]  
(Division Sign-Off)

**Division of General, Restorative,  
and Neurological Devices**

510(k) Number K073509