



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

AUG 15 2008

Ms. Linda Nelson
Regulatory Affairs Manager
SteriTec Products Manufacture Company, Incorporated
599 Topeka Way Suite 400
Castle Rock, Colorado 80109

Re: K080136

Trade/Device Name: Instant Readout Chemical Integrator Challenge Pack,
Model No.: LT 105
Regulation Number: 21 CFR 880.2800
Regulation Name: Chemical Indicator
Regulatory Class: II
Product Code: JOJ
Dated: August 05, 2008
Received: August 12, 2008

Dear Ms. Nelson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K080136

Device Name: SteriTec Instant Readout Chemical Integrator Challenge Pack, Model No.: LT 105

Indications For Use:

The SteriTec Instant Readout Chemical Integrator Challenge Pack is a steam sterilization process challenge device that meets the FDA Chemical Indicator Guidelines performance requirements for chemical integrator test packs. The Integrator Challenge Pack monitors pre-vacuum steam sterilizers operating at 132°C (270°F) or 135°C (275°F). When exposed to steam sterilization conditions the Instant Readout Chemical Integrator card gives an integrated response by changing from purple to green. The integrator card inside the test pack also meets emulator performance standards according to ISO 11140-1:2005.

Critical Parameters in hospital steam sterilizers

Temperature/ Sterilizer Type	Indicator "A"	Indicator "B"	Indicator "C"
132°C (270°F)/ Pre-Vacuum	4 minutes	10 minutes	20 minutes
135°C (275°F)/ Pre-Vacuum	3 minutes	8 minutes	16 minutes


Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use X
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

Page 1 of 1

510(k) Number: K 080136